


*Luis M. de la Maza*  
*Marie T. Pezzlo*  
*Ellen Jo Baron*

# Color Atlas of Diagnostic Microbiology



# Color Atlas of Diagnostic Microbiology

**Luis M. de la Maza, M.D., Ph.D.**

Professor Department of Pathology  
Director Division of Medical Microbiology  
University of California, Irvine Medical Center  
Orange, California

**Marie T. Puzlo, M.A., F(AAM)**

Administrative Director  
Department of Pathology  
University of California, Irvine Medical Center  
Orange, California

**Ellen Jo Barren, Ph.D., F(AAM)**

Adjunct Associate Professor of Medicine  
University of California, Los Angeles  
Clinical Associate Professor of Molecular Microbiology and Immunology  
University of Southern California  
Los Angeles, California

 **Mosby**

St. Louis Denver Boston Chicago Dallas Houston New York Philadelphia Toronto  
London Madrid Mexico City Singapore Taipei Seoul Tokyo Washington



978 0 7163 0621 0  
M 0 000000

The Author/Publisher Relationship  
Editor/Book Editor  
Book Review Editor/Book Editor  
Project Manager Linda McKinley  
Production Editor Lynn Becken  
Designer Elizabeth Fox  
Manufacturing Executive Theresa Fackel

Copyright © 1997 by Mosby-Year Book, Inc.

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without written permission of the publisher.

Permission to photocopy or reproduce solely for internal or personal use, or limited for internal or personal use, registered with the Copyright Clearance Center, provided that the base fee of \$4.00 per chapter plus \$1.00 per page is paid directly to the Copyright Clearance Center, 27 Congress Street, Salem, MA 01970. This consent does not extend to other kinds of copying, such as copying for general distribution, for advertising or promotional purposes, for creating new collective works, or for resale.

Printed in the United States of America  
Composition and typography by/under film by JPS Graphics  
Printing/Binding by Wakefield Press, Inc.

Mosby-Year Book, Inc.  
1180 Phoenix Industrial Drive  
St. Louis, Missouri 63144

Accession Number: Mosby Book Number  
Q-0191 0621 0

## Preface

**I**nfectious diseases caused more than 17 million deaths in 1995, almost one third of all deaths that occurred worldwide that year. The World Health Organization has listed 18 "new" infectious diseases identified during the last 12 years that are contributing to these statistics. This ominous rise in infectious diseases is accompanied by an almost equally alarming decrease in the numbers of programs available for training clinical laboratory scientists and microbiologists.

Diagnostic microbiology is a visual science; ability to recognize characteristics both macroscopically and microscopically is the hallmark of competence in the field. A consistent and inclusive source of microbiological color images is virtually nonexistent. We have sought to remedy that situation with this volume.

The book contains 735 full-color images, encompassing laboratory safety practices, specimen collection and transport devices, specimen handling, colony morphologies, microscopic images, biochemical and other tests, and susceptibility testing. The microbes that cause human disease are found among viruses, fungi, parasites, and bacteria, and this book includes representative images from each group. Emphasis has been placed on the pictures and text is minimal, with just enough information to place the photographs *in context*. The book is meant to be a collection of some of the most important visual images in human diagnostic microbiology.

All instructors helping students to learn medical microbiology know how difficult it is sometimes to find that "typical" image to illustrate a point. We expect this book will provide instructors easy access to those images that are frequently used in the teaching of diagnostic microbiology. The book should be useful to medical students taking their microbiology courses, pathology residents and infectious diseases fellows, pathologists and physicians, and of course laboratory scientists and diagnostic microbiologists, both beginning and returning.

We have attempted to include in this atlas those images that cover the most frequently detected pathogens and most frequently performed tests in the clinical laboratory. We realize that this is an impossible task, because the numbers of pathogens and types of tests vary greatly from site to site. If we have failed to cover a subject that you think should be included, please feel free to send us your suggestions, or even better, specimens and samples, to be considered for the next edition of this book. Finally, we ask all of our readers, and particularly our critics, to bring to our attention any errors that they may find in the book. We thank all of you in advance for your help.

## Technical Note

The light and fluorescence microscopic pictures were taken with a Zeiss Universal microscope (Carl Zeiss Inc., West Germany) equipped with Zeiss and Olympus (Olympus Optical Corp., Ltd., Japan) lenses and a Nikon F5-FA camera (Nikon Corp., Tokyo, Japan). The microscopic images were obtained using a Nikon E camera with a Micro-Nikkor 95 mm F1.8 lens and a Contax RTS camera with a Carl Zeiss J-Plassar 40 mm F2.8 lens. Kodachrome 25 Professional film (Eastman Kodak Co., Rochester, New York) and Tri-X and Provia Fujichrome film (Fujifilm Photo Film Co., Ltd., Tokyo, Japan) were used for most of the images. The electron micrographs were taken with a Philips EM800 electron microscope (Philips Electronic Instruments Co., Eindhoven, The Netherlands) and Electron Image film 90-163 (Eastman Kodak Co.). The magnification stated in the legends corresponds to approximately that of the printed figure. However, for calculation, the measurement provided in the legend should be used.

## Acknowledgments

Special thanks for their extra effort go to all the staff of the Division of Medical Microbiology at UCL Medical Center. We are especially grateful to Sandra Armaso, Jessica Blanding, Miryam Cordero, Neil Entwistle, Kaye Evans, Jazel Shaghi, Sachin Schuman, Grace Tan, and Miriam Varkatas for their help in preparing the specimens and collecting the material for inclusion in this book.



**T**

*his book is dedicated to all past and present members of the staff of the Division of Medical Microbiology at UCI Medical Center.*

*Their commitment to patient care and teaching made this book possible.*

*And to Frank Perzle and James C. Taylor for their unending support.*

# Contents

- 1 Laboratory Safety, 1
- 2 Specimen Collection Containers 9
- 3 Specimen Processing, 14
- 4 Gram Stain 25
- 5 Micrococci 32
- 6 Streptococci 57
- 7 Anaerobic Gram-Positive Bacilli & Actinomyces spp. 46
- 8 Enterobacteriaceae 56
- 9 Other Gram-Negative Microorganisms 72
- 10 Anaerobic Bacteria 88
- 11 Mycobacteria 95
- 12 Microbial Pathogens Isolated and/or Identified by Tissue Culture or Other Special Methods, 102
- 13 Antimicrobial Susceptibility Testing, 107
- 14 Mycology, 113
- 15 Parasitology, 146
- 16 Virology, 176
- 17 Immunocytology, 194
- 18 Molecular Techniques, 200

## CHAPTER 1 Laboratory Safety

**D**uring the last decade there has been a dramatic increase in the number of guidelines, recommendations, regulations, and standards introduced for the safety of personnel working with potential pathogens in clinical laboratories. It is mandated by OSHA (Occupational Safety and Health Act created 1970) that guidelines be implemented and practiced and that protective clothing, containment devices, and decontamination equipment and materials be available in the workplace. Procedures regarding biosafety requirements and practices must be available and familiar to all laboratory personnel.

### Protective Clothing and Equipment

An important component of Universal Precautions and safe work practices is the effective use of protective gloves, eye and mouth protection, and laboratory coats all considered to be Personal Protective Devices. Protective gloves must be worn whenever there is a potential for direct contact with infectious materials. Disposable latex or vinyl gloves are recommended and should immediately be discarded in a designated waste container when removed.



**FIG. 1-1** Remove one glove by pulling the outside bottom of that glove toward the palm and fingers with the other gloved hand, turning the glove inside out as it is removed. Do not touch your skin with the outer glove surface.



**FIG. 1-2** Pull the glove completely off and hold this inside-out glove to the gloved hand. Do not touch the outer surface of either glove with the ungloved hand.





**1-11** Using the equipped hand, grip the inside of the bottom of the second glove and pull the other glove off.



**1-12** Pull the second glove over the contaminated glove that is held in the hand, starting to completely stretch out to ensure the first glove is not to avoid contamination of the ungloved hands.



**1-13** All gloves should be immediately discarded in a waste container designed for contaminated materials.



**1-14** Eye covering such as goggles should be worn whenever working with body fluids or other specimens that may splash or penetrate facial eye areas. Goggles with sideguards should be provided for such studies.



**1-15** A face mask is worn to prevent mucous membranes of the nose and mouth and to prevent air around particles that may enter the respiratory or digestive tracts. The mask should fit tightly around the face and provide the OSHA-mandated level of protection. Face masks should always be worn if there is a chance of splashing clinical material. In some settings waterproof face shields may be required for some tasks involving aerosols and blood.



**1-10** This photograph illustrates an improperly fitted face mask that does not provide adequate protection. There is a potential for microorganisms to enter through the gaps and come into contact with mucous membranes or to be inhaled.



**1-11** A full face shield protects the eyes, nose, and mouth.



**1-12** A full face shield is worn to protect the eyes, nose, and mucous membranes.



**1-13** Full contact laboratory coats with closed ends and wrists should be worn to cover all potential exposed skin and clothing. Full contact laboratory coats should be worn over street clothes and should be removed before leaving the laboratory.



**1-132** Gait suits protect the skin on the arms.



**1-133** Protective gloves should be worn even the least risk to provide a proper sterile barrier and biocontainment.

## Disinfection, Waste Disposal, and Sterilization

### Disinfection



**1-140** Decontamination of work surfaces at the end of each shift and cleaning of spills should be performed with a fresh 10% solution of bleach/bleach bleach.



**1-141** Paper towels should be placed on the contaminated work surface and wetted with diluted bleach/bleach bleach (approximately 1:10) to prevent spillage at the end of work shift.



**1-142** Bleach and paper towels should also be used for disinfecting the laboratory surface after each spill. Drop dead with waste records.

## Waste Disposal

Contaminated materials should be disposed of in proper containers that are properly labeled and fitted with lids. Separate containers are available for plastic and glass and needles. It is important to dispose of materials in proper containers to avoid potentially serious injury or illness to individuals coming in contact with these materials.



**9-17** Agar plates are discarded in a properly labeled container with a lid.



**9-18** Needles with attached syringes are discarded into hard-sided puncture-proof containers. Needleless syringes are to be capped unless a specially designed, commercially available holder for recapping needles is used for that purpose.



**9-19** Glass vials containing liquid samples are also discarded into hard-sided, puncture-proof containers.

## Sterilization

All materials contaminated with potentially infectious microorganisms must be decontaminated before disposal. An autoclave is the most common instrument used to perform this function.



**1-200** Leavely packed materials that have been aseptically sealed (disposed in designated containers) (e.g., spill bags) and liquid wastes, are placed in the autoclave. The bags, which should never be filled completely, are loosely tied at the top to allow the steam to circulate and effectively decontaminate the contents. Adding approximately 1 cup of water to the bag facilitates steam production.



**1-201** The autoclave door is closed tightly. In the completion of the cycle, the door is opened slowly by the operator wearing heat resistant gloves (tools that must be passed a patient before being washed).

## Biological Safety Cabinet

The biological safety cabinet (BSC) is one of the most common containment devices used in laboratories. Air is decontaminated by high efficiency particulate air (HEPA) filters.



**1-202** A commonly used Class II biological safety cabinet that filter circulation both the operating and exhausted air. A BSC should be kept free of clutter. Routine inspection of the cabinet's certified technician is required.



**1-233** Regular monitoring of air flow should be performed to ensure proper velocity. Mechanical problems (filter clogging, and other routine maintenance at lowered velocity) which may allow airborne particulates escape the cabinet.

## Other Safety Practices

Use of aerosol-containing safety buckets in centrifuges, venting tightly sealed tubes in a biological safety cabinet (BSC), and securing gas tanks to a holding device mounted on the wall are other additional safety practices that are required to ensure a safe work environment.



**1-234** Tightly sealed centrifuge tubes are placed into safety carrier, then are themselves tightly sealed before centrifugation. This will avoid mixing or aerosol if tubes are broken during centrifugation.



**1-235** Care should be taken when opening a sealed centrifuge rotor to avoid contamination from a possible broken container during centrifugation.



**1-126** Pipette-aided procedures should be performed in a BSC only. The tube should always be inserted to reach the centrifugal particles. A 10-second wait before opening the tube avoids aerosol contact, especially when capable materials are being handled.



**1-127** All use vials should be secured in place by a locking device.

## Mailing Containers

Potentially infectious material must be shipped according to the federal codes of Interstate Shipment of Biologic Agents.



**1-128** The vial(s) must be placed in an inner seal container which is then placed in an outer leak-proof tube. The outer container must have sufficient material to absorb the entire fluid contents in the event of a leak.



**1-129** An official United States label must be affixed to the outer container.

## CHAPTER 2 Specimen Collection

### Containers

**I**t is the responsibility of the diagnostic microbiology laboratory to select and provide transport devices for collection of specimens from a variety of anatomical sites. Most specimen collection containers incorporate transport media that support the viability of microorganisms encountered in clinical specimens. The interpretation of results is dependent upon the quality of the collected specimen and the transport conditions. Improperly transported specimens may result in failure to isolate the causative microorganism. One of the most common and widely used transport systems is a plastic tube containing a sterile polyester-tipped swab and medium to prevent the drying of microorganisms, maintain the pH, and minimize oxygenation.



**FIGURE 2-1** Sterile diagnostic swabs culture collection and transport system consisting of a plastic tube containing one swab-tipped swab and modified Stuart's transport medium to support aerobic organisms (BBL Microbiology Systems, Cockeysville, MD). The capsule is located at the base of the tube and is covered by a protective sleeve on the outside. The capsule must be cracked after the specimen is collected, and the swab is placed into the container to keep the specimen on the swab moist. The use of cotton swabs should be avoided because they contain fatty acids that may be inhibitory to some bacteria.



**FIGURE 2-2** Neopenterygonal-wrapped swab (Kingfisher Type IV, Spectra-Gen Diagnostics, Chesham, NJ). This swab is used for collection of specimens from specific sites, especially for samples requiring and using modified specimens. Culture slant media should not be used for collection of certain specimens, because they can be toxic to some strains of *Mycobacterium tuberculosis* complex. They are useful for collection of *Mycobacterium tuberculosis* cultures.






**2-3 Respiratory specimen collection system.** A similar collection system as described in Figure 2-2, however this unit is used for collection of respiratory specimens for recovery of *Moraxella parvula*. The tube contains *Acinetobacter baumannii* transport medium, which also supports the viability of *Pfeifferia parvula*.



**2-4 Sterile, disposable collection and transport system for recovery of viruses (BBL Microbiology Systems, Bockharyville, MD).** The cap and holder are shown separately as the top end is inserted at the bottom. The tube contains a virus-typed media and viral transport medium containing disinfectant and antibiotics. This prevents specimen drying, helps maintain viral stability, and restricts the growth of other microbial contaminants.



**2-5 Marpen collection kit** containing two discs of deoxyribonucleic acid (DNA) for preparation of a direct smear at the time of collection and a slide holder (Type M100010, Sanjour, CA). The slide will be stained later in the laboratory with a direct fluorescent antibody stain.



**2-6 Chlamydia trachomatis collection kit** containing two discs of deoxyribonucleic acid (DNA) for preparation of a direct smear at the time of collection and a slide holder (Type M100010, Sanjour, CA). The slide will be stained later in the laboratory with a direct fluorescent antibody stain.



**8-7** Viral-chemically transport medium. Specimens for detection of viruses or chlamydia should be placed into a tube of vireo-transport or buffer transport medium containing fetal calf serum, buffer, and penicillin (Baxter Diagnostic Inc., Deerfield, IL)



A



**8-8** Bio-Bag Environmental System (BBL Microbiology Systems, Cockeysville, MD). The specimen for anaerobic culture is collected in these leak detectable or leak less detectable as a seal, and then placed in the gas impermeable methylene violet bag that contains suspension of methylene violet and hydrogen  $\text{CO}_2$  generator. The bag is sealed and each separate incubated to produce anaerobic conditions.



B

**8-8** A. Transgrow media (BBL Microbiology Systems, Cockeysville, MD) for recovery of anaerobic organisms. For most, the surface of the bottle contains modified Thayer Martin agar that has been prepared under  $\text{CO}_2$ . B. Junction plate (BBL Microbiology Systems, Cockeysville, MD) for recovery of anaerobic organisms. Medium is prepared in a bag, plastic container with a capillary lid. The specimen is inoculated directly and the mixture from the medium diffuses a barbitone color to create the  $\text{CO}_2$  atmosphere.



**2-10 Para-Pak Parasitology Collection Kit (Meridian Bioscience, Cincinnati, Ohio).** The commercially prepared kit contains one vial of modified polyvinyl alcohol (PVA) and one vial of buffered neutral formalin. Placing the stool specimen into these vials immediately after collection preserves the morphology of the parasites.

There are a variety of media for the collection and recovery of protozoans from stool. Some of these media can be used alive, while others require incubation and detection in an automated system. They all contain between 0.25% to 0.5% sodium polyacrylate (SPA), and the recommended blood-to-medium ratio is 1:2.5 to 1:10, depending upon the system. Examples of these media appear in Figures 2-11 to 2-14.



**2-11 Becton-Dickinson Becton-Dickinson blood culture bottle (BDL) Microbiology System, Cockeysville, MD.)** represents a non-automated blood culture system. Bottle contains oxygen gas head and an open slide paddle with three separate incubation, MacConkey, and yeast. Blood is inoculated into the bottle within 1 hour of collection and transported to the laboratory. In a hematology cabinet, an open slide paddle is attached to a bottle, inverted, and incubated. Subcultures are performed by reinserting the bottle again to incubate the paddle.



**2-12 Isometer tube-anticoagulant blood culture system, 10 and 10-ml tubes (Whisper Laboratories, Greensburg, PA.)** are recommended for the collection of protozoans and fungi. The tube contains a blood cell layer coat of heparin, polyacrylate (PA) flocculant, and SPA. Specimens collected in the 10-ml tube are centrifuged at 1000  $\times$  g for 10 minutes. The supernatant is discarded, and the sediment is placed into colored media. Specimens collected in the 1-ml tube are placed directly into media.



**12-18 BACTEC blood culture bottles** are used to detect bacterial and fungal infections in 24-hour, 27-hour and 48-hour bottles. The bottles are used with the Vacutainer tubes containing BACTEC Oxoid (Becton Dickinson, Sparks, MD). These media are intended for recovery of aerobic and anaerobic microorganisms.



**12-19 BACTEC 500 aerobic and BACTEC 500 ANAEROBIC blood culture media** to be used with the Vacutainer blood culture tubes (Becton Dickinson, Sparks, MD).

**12-20 Blood collection vacutainer tubes.** These blood collection vacutainer tubes (purple-top with lithium heparin, white-top with citrate, and yellow-top with EDTA) are used for the serologic diagnosis of infectious diseases.



## CHAPTER 3 Specimen Processing

**T**here are a number of important steps involved in specimen processing that include (1) entry of patient information into a computer system or manual log; (2) visual examination of the specimen to determine if acceptance criteria are met; (3) specimen preparation, selection, and inoculation of media; and (4) preparation and microscopic examination of direct smears. Upon receipt into the laboratory the technologist should verify that the specimen is labeled with the appropriate patient information. Unlabeled and mislabeled specimens should be handled according to the written laboratory policy. Any additional patient information needed to process the specimen should be retrieved from the laboratory information system (LIS) or other source.



**FIG 3-1** Entry of specimen information into the computer system. Upon receipt into the laboratory the system should be entered into the LIS or manual logbook. Patient information including diagnosis and antimicrobial therapy, complete specimen information including source, collection site, time of collection, and test requests, and specimen accession number should be available before specimen processing.



**FIG 3-2** Examination of the specimen container to confirm that the specimen is properly labeled. Attention to acceptance criteria for handling unlabeled and mislabeled specimens should be available and followed.




**2-23** Verifying that the information on the specimen container label corresponds with the laboratory request form. Improperly labeled specimens should be processed only if correction action and documentation are complete.



**2-24** Carefully examining the specimen to determine if it matches the specimen description provided on the laboratory request form.



**2-25** Examining the specimen container for the presence of visible leaks and other potential problems. Universal precautions should be used when handling all specimens.



**2-26** Microscopic examination is an important laboratory test for all specimens. However, this laboratory test is requested from experienced users. Some microscopic appearances of the four specimen specimens shown in this figure are: (top left) thick smudged particles; (top right) many white particles; (bottom left) smudged fibers; (bottom left) thick, smudged areas surrounded by cells.



**2-7** A well-written procedure manual must be available in sterile surface positions for reviewing and inspecting specimens and guidelines for specimen processing. This manual should be consulted for complete information on specimen processing, inoculation of media, and incubation and atmospheric conditions.



**2-8** A collection of enriched, differential, and selective media should be available for culture. A daily supply of media should be stored at room temperature and easily accessible to the user.



**2-9** A daily supply of prepared laboratory media should be available and used only on selected specimens as outlined in the procedure manual. **2-10** The appropriate media and dyes should be selected on the basis of the specimen source and culture requirements. These media should be labeled with the patient identifying information number before they are incubated.





**3-103** A selection of irrigating loops including extruding loops should be among the supplies for processing.



**3-104** An excellent method for preparing a cast is to first fresh tissue is to press the tissue directly into the thin surface. Needle slides should be used to pressure the tissue flat until tissue malleable.



**3-105** Small specimens placed in covered Petri-dishes before inoculation into media like methyl methylene blue is useful. The specimen is placed into a sterile petri dish, and use a two-needle syringe and syringe and combine the specimen with a microorganism in laboratory.



**3-106** Another method for homogenizing tissue is using a disposable homogenizing kit (Biospec Products - Cary 22). This system minimizes contamination and splashing of the specimen.



**3-107** A portion of the specimen is placed into the form of a sterile tube, and a sterile 0.5 ml of medium broth is added to increase the specimen.



**3-108** The specimen is homogenized in the specimen bottle and pushes the probe down onto the tissue. The procedure should always be performed in a biological safety cabinet.






**FIGURE 10-100** Specimens must be examined carefully before being below of media. A representative portion that has an absent and appearance that lacks a typical structure, presence (colorless, abnormal color or other signs of an infectious process) should be observed.



**FIGURE 10-101** *a*, Values as shown in this figure must be distributed from specimen. The impression that a specimen contains cells is confirmed by preparing a smear and examining the specimen microscopically. If the microscopic examination confirms that the specimen is not a specimen, the specimen should be rejected and accompanied the routine bacterial culture. *b*, Gram stain (low power) of a good quality specimen specimen characterized by field of numerous epithelial cells. Most of the tissue cells seen are mononuclear. *c*, Gram stain (low power) of a poor quality specimen specimen, characterized by numerous squamous epithelial cells. The values of an in a specimen will not accurately reflect the clinical situation and may yield misleading results.



**2-20** The portion of the specimen used likely to harbor the infecting microorganisms is selected. A sample of this specimen specimen should be taken from the central (blood) portion.




**2-21** Selection of specimens that appear macroscopically to be from an infected source, such as the routine field (class) collection, being not acceptable unless numbers/contaminants should be avoided. The same state should be repeated class to microscopy between the microorganisms and the gross appearance of the specimen.



**2-22** The specimen should be used when preparing a dry smear. When transferring a specimen onto a slide with a wet, the wet should be with force stain of the surface using a stick. And the specimen, repeat dry smears and require a new specimen.



**2-23** When a wet specimen is used to transfer the slide into a slide, enough specimens must be placed on the slide to form a circle of approximately 1.5 mm in diameter. The slide should be dry in the biological safety cabinet before it is used and stored.




**Fig. 20-200** It may be necessary to use a forceps perpendicular to the slide to deposit a small amount of specimen directly on the slide to form a circle of approximately 1.5 cm in diameter. The bottle on the other hand, should be deposited in a temperature and atmosphere suitable for processing; success would be a development of body fluids.



**Fig. 20-201** A small drop of a liquid specimen is placed in a well on a slide and covered immediately with a cover slip. This preparation is examined under a microscope to detect bacterial seeding under various parameters and for other investigations.



**Fig. 20-202** A. If a specimen is in a vial, the open glass should be sterilized by rotating the vial across the quadrants of the plate. B. If only one vial is available, it should be rotated in a small amount of broth (0.75 to 1.0 ml) and the resulting suspension can be used to test this plate and prepare the rest. Sampling should be performed in a biological safety cabinet.



**2-25d, cont'd.** **C.** The primary inoculum is spread with a loop using a back-and-forth motion. A wire loop should be washed between each successive quadrant (cont'd.)



**2-26** Force forceps are used to process bacteria such as an herbicide. The forceps are held at an angle to avoid the outer surface of an agar plate, and then to place the specimen into a sterile container (cont'd.)




**2-27** Inoculum in syringes should be discarded, this one first inoculated at the lab into a sterile culture with one set to the study body fluid.



**2-28** If toxic body fluids are being tested by centrifuge then the supernatant should be separated off the sediment and a fresh pipette should be used to mix and inoculate the sediment into media and slides. Because this small to be centrifuged should be inoculated directly with a pipette.



**2-29** A variety of atmospheric conditions and incubation temperatures should be available for growth of pathogenic microorganisms. Incubation temperatures should include at least 30°C, 35°C and 42°C along with aerobic, anaerobic, microaerophilic and microaerophilic conditions. Most routine aerobic cultures should be incubated at 35° to 37°C in a humid atmosphere of 5% to 7% CO<sub>2</sub>. The temperature, CO<sub>2</sub> content and humidity should be monitored daily. The CO<sub>2</sub> content in the media with a gauge and a stopcock (cont'd.)



**3-22** The temperature, humidity, and ventilation must be monitored.



**3-23** Media should be distributed as soon as possible after formulation. An aliquot here formulated media is wrapped in styrofoam to maintain a cool and moist environment.



**3-24** Pre-sterilized and sealed plastic bags for bulk media are used. The bags can be a short time before and after inoculation.



**3-25** An incubator, although not an ideal for creating an anaerobic growth that holds only incubates until an incubator jar for three or more plates.



**3-26** The Gulfair incubator jar (which carries Diagnostic Microbiology Systems, Cockeville, TN) contains a hydrogen seal (hydrogen on edge, supposedly maintains blue color) and a venting feature in the lid.



**2-24** The CoPath assembly growth (Becton Dickinson Microbiology Systems, Cockeville, MD) is based on the same principle. The bags are impermeable and therefore contain its own gas-generating kit (carbon dioxide).



**2-25** The assembly of CoPath bags is in a anaerobic atmosphere and also increased in a pouch and jar.



**2-26** Recovery of some fastidious microorganisms like *Campylobacter* and *Nitrosomonas* requires a microaerophilic atmosphere (partial O<sub>2</sub>) that can be created by barilets candle in a jar (jar).



**2-27** The microorganism placed a blood culture bottle into the BACTEC 520i Fluorescence system and blood culture system (Becton Dickinson Diagnostic Instrument Systems, Sparks, MD) is automatically combined into a culture of the medium and analyzed every 30 minutes by a photoluminescence system to a computer.



**2-400** This technology is placed in the incubator into the UV instrument (BioLabs, Boston, MA). A connector is placed onto the top of the bottle before it is used. Each bottle is exposed to a coil of the transmitter and receiver every 30 minutes by a pressure gauge attached to a computer.



**2-401** Facility workers should be instructed that each area is to be secured as possible after the slide is prepared and fixed.



**2-402** Mailed smears should be read immediately after 24 hours. Abnormal results must be relayed immediately to a physician. The date and time of the call must be documented. The type of smear result that require a verbal report should be listed in the procedure manual.



**2-403** Special processing should also include a variety of rapid direct tests that can assist the physician with the patient's diagnosis. Examples include direct immunofluorescence as a cryoprecipitation test application test. Smear tests

## CHAPTER 4 Gram Stain

**T**he Gram stain is used to classify microorganisms on the basis of their Gram staining characteristics, size, shape, and arrangement of cells. It is one of the few tests in clinical microbiology that can assist in the rapid presumptive diagnosis of an infectious disease. It is also used to assess the quality of the clinical specimen based on the relative cellular content. Bacteria are first stain gram-positive, gram-negative, or gram-variable. The gram-variable appearance can be due to over- or under-decolorization age of the microorganisms, influence of antimicrobial treatment, and other factors. The staining reaction is dependent on the microbial cell wall composition. The tightly cross-linked peptidoglycan layer and the teichoic acid found in gram-positive microorganisms cause them to be resistant to acetone-alcohol decolorization. They retain the crystal violet stain and are purple in color. The hypoteichoic-acid-rich cell wall of gram-negative microorganisms is disrupted by acetone-alcohol and the crystal violet leaks out of the less tightly cross-linked cell wall structure. The subsequent counter stain can then be seen, resulting gram-negative organisms pink in color.

The size of most bacteria ranges from a large cell (10µm to 30µm in length) to a small cell (1µm to 2µm in length). The shapes of bacteria are usually described as cocci, bacillary, or coccobacillary, and the arrangements as pairs, clusters, chains, branching, filamentous or coryneform. Yeast can be seen as single cells, often with budding and/or hyphal filaments. Microorganisms may be seen intracellularly within somatic cells. Staining characteristics of individual bacterial cells include bipolar, beaded and irregular. Bacteria may appear rounded, pointed, flattened, or needle. Microorganisms other than bacteria and yeast observed in Gram stains, although they are better characterized with other stains, are *Trichomonas trophozoites*, *Pharyngocystis serial cysts*, *Trichomonas gonads trophozoites*, and *Strongyloides larvae*.


Somatic cells can also be seen in a Gram stain, best observed with methanol-fixed smears as shown here. White blood cells (WBC) and epithelial cells tend to stain a pink color whereas red blood cells (RBC) appear tan to buff colored. The presence and rough quantitation of somatic cells should be noted when interpreting a Gram stain.

Gram stains should be examined under low power magnification (10X objective) for the quality of overall staining, the thickness, and for evaluation of somatic cells; microorganisms should be observed under oil immersion (100X objective).






4-1 Polymorphous leukocytes (PMN) and gram (and trypan blue) stained diplococci suggestive of *Streptococcus pneumoniae* (x1,250). The clear area surrounding the bacterial cells is suggestive of a capsule. A rapid-level bacterial culture, too may be performed on the clinical specimen to confirm the organism identification.



4-2 PMN, red blood cells (RBC) and gram-positive cocci in pairs.




4-3 Gram-positive cocci in chains resembling streptococci (x1,250).



4-4 Gram-positive cocci in pairs, chains and clusters resembling streptococci (x1,250).



4-5 RBC and gram-positive cocci (x1,250). Note the somewhat large, irregular and may be also suggestive of streptococci (streptococci, streptococci, or Streptococcus spp.) diplococci in pairs with irregular shape and appearance (multicellular wall cocci) (x1,250).



4-6 Leukocytes, severely stained, gram-negative, round and gram-variable cocci (x1,250). In each case, the gram-variable appearance may be due to phagocytosis and partial degradation of the cell components by the PMN.




**4-7** PAS and iron-haematoxylin (in place) show red lung tissue and clusters remaining of mixed epithelioid and mesothelial (x1,200)



**4-8** PAS and iron-haematoxylin, gross granular nuclei with pleomorphic nuclei suggestive of carcinoma (x1,200) (some appear vacuolated, epithelioid and pleomorphic) which is suggestive of sarcomatous or *Utricularia* neoplasia



**4-9** Typical squamous epithelial cells with eosinophilic granules (eosinophilic) in nuclei (nuclei with rounded or flat ends) some in double suggestive of histiocytosis (x1200) Clear background and large epithelial cell edges typical of pleural cells with normal cytoplasmic structure



**4-10** Endothelial PAS and iron-haematoxylin, large granular blue nuclei resembling *Chromidium* (pleomorphic nuclei) (x1,200)



**4-11** PAS and large granular nuclei (x1,200) There appears to be few nuclei that are irregular in shape and overlapping, giving the appearance of coal in clusters. The basal appearance may be due to the presence of spines, particularly those by the PAS or because the structure is a small piece of haematoxylin which coats unevenly with the Gross stain



**4-12** Large, granular nuclei with eosinophilic (x1200) Some cells appear as dendritic, suggestive of *Chromidium* stain


**4-10** *Trich* and long, thin, beaded branching intensely stained green-pinkish filamentous structures of *A. nidulans* or *Aspergillus* spp. (CCL258) (the green may be confused with a modified Gram's) inter-larva stain. *Aspergillus* spp. are partially acid fast.



**4-11** Branching filamentous slender green-pinkish beaded filaments from an aerobic colony. This morphology is suggestive of *Trichomyces* (C 1274) from a colony of a pink filamentous aerobic fungus *Trichomyces* *Chet* *Chet* spp.




**4-12** *Trichomyces* spp. (segment like) curved beaded thin filamentous structures (C 1254)




## Gram-Negative Microorganisms


Gram-negative microorganisms can appear as cocci, coccobacilli, and bacilli. Cocci can appear as singles, as pairs with flattened adjacent sides, and as clusters. Coccobacilli are usually small to medium in size. Bacilli can vary in appearance from small, finely staining rods to large, plump rods with regular staining.




**4-16** Gram-negative cocci in singles and pairs. The adjacent sides of the diplococci appear flattened ( $\times 1,700$ ). This microorganism is *Streptococcus pneumoniae* stained from a blood culture broth.




**4-17** Gram-negative coccobacilli in pairs ( $\times 1,000$ ). This microorganism is *Klebsiella aerogenes* from a blood culture broth.




**4-18** Gram-negative coccobacilli in clusters ( $\times 1,700$ ). This microorganism is *Klebsiella pneumoniae* from a blood culture broth.




**4-19** Gram-negative coccobacilli in pairs. This microorganism is *Klebsiella pneumoniae* from a blood culture broth ( $\times 1,700$ ). The morphological similarities between *Klebsiella* and *Enterobacter* are usually indistinguishable to gram-negative diplococci. The HACEK group of fastidious gram-negative bacilli (*Haemophilus*, *Aerobacter*, *Citrobacter*, *Edwardsiella*, and *Kingella*) often do not stain as typical rod-shaped, but also appear as small, coccobacilli.



**M1156** BCC and plenty gram-negative bacilli with bipolar staining in a blood culture broth, suggestive of "rodentic" belonging to the Bacillaceae (M1156). The most dark bacilli are characteristic of *Listeria* spp.




**M1158** BCC and long gram-negative bacilli with rounded ends suggestive of *Proteus* or *Providencia* spp. (M1158) in a blood culture broth.




**M1705** Freshly differentiated PABs and double pleomorphic B100s together with flagellates of *Haemophilus* spp. or *Pasteurella* spp. (M1705).



**M1159** BCC and pleomorphic bacilli during gram-negative bacilli with penicillins suggestive of *L. monocytogenes* spp. (M1159).



**M1706** PABs (BCC) and slightly curved comma-shaped yellow and pinkish microorganisms suggest *Campylobacter* spp. (M1706). **M1708** Bright masses of *Streptococcus* in a prepared blood smear (M1708). These helical bacteria measure approximately 3 µm in diameter with a pole to pole distance made up to 10 base cells.




## Yeast

Yeast cells are larger than bacteria, usually ranging from 5 to 15  $\mu\text{m}$  in size. They can appear as single cells with buds, in clusters, and with typical filaments. Staining of cells can be somewhat speckled rather than confluent.



**4-44** Budding yeast cells with single blastoconidia and typical yeasts, staining surrounded by clear halos (note redness by the polysaccharide capsule) representative of *Candida* spp. (X 1150)




**4-45** Budding yeast cells and an irregularly stained polysaccharide capsule representative of *Candida* spp. (X 1250)



**4-46** Heavily stained 5 to 4  $\mu\text{m}$  wide branching mycelium hyphae with purple cells and 45° branching configurations representative of *Aspergillus* spp. (X1150)

## CHAPTER 5 *Micrococcaceae*

**T**hree genera of the family Micrococcaceae have been associated with infections in humans: *Staphylococcus*, *Micrococcus*, and *Sarcoscyces*. Of these, *Staphylococcus* spp. are by far the most common cause of human infections. Both *Micrococcus* and *Sarcoscyces* spp. can be found in the environment, as normal flora of the skin and respiratory tract, and as pathogens associated with catheter tips and other hardware.




**5-1** Gram stain showing BCG and gram-positive acid-fast bacilli, and gram-negative coccobacilli of *Micrococcus* (120 $\times$ ). In this figure, the coccobacilli are mostly loose.



**5-2** Culture of *Staphylococcus aureus* on 5% sheep blood agar. Typical appearance of a confluent bacterial "mat" culture grown overnight on a complete 5% animal.



**5-3** **A**, Culture of catalase-negative *Staphylococcus* on 5% sheep blood agar. Typical culture of catalase-negative *Staphylococcus* on 5% sheep blood agar and 5% **B**, Culture of *Micrococcus* spp. on 5% sheep blood agar. The colonies of *Micrococcus* appear brown, yellow. *Micrococcus* spp. usually grow consistently over 48 hours, appearing 48 hours of incubation. After 48 hours, the colonies are visible.



**8-4 Gram stain of *Streptococcus* species.** Microscopically indistinguishable (in shape) from staphylococci and appear to be more "colony" than gram-positive cocci ( $\times 1250$ ).



**8-5 Catalase test.** The test is performed by adding 1% hydrogen peroxide ( $H_2O_2$ ) in a test tube on a glass slide or by adding a cotton swab on a reaction slide to a drop of  $H_2O_2$  on a slide as shown here. The appearance of bubbles indicates that the organism oxidizes (has hydrolyzed  $H_2O_2$ ) into oxygen plus water. Staphylococci and micrococci are differentiated from other aerobic gram positive cocci by a positive catalase test. (left) No bubbles appear in a negative test result. (right)



**8-6 Fermentation of glucose is one of the methods used to differentiate staphylococci and streptococci.** When grown in an oxidation-fermentation (O/F) medium, staphylococci produce acid (yellow color) from glucose under anaerobic conditions, resulting in a yellow precipitate at the surface of the agar (tube on the left). The microorganism also produces gas (bubbles) (tube on the right).



**8-7 Glucose oxidation in O/F medium.** Micrococci are aerobic glucose oxidizers (yellow color in the open tube) but do not ferment the carbohydrate (left, no color change [green] in the sealed tube). The acid residues produced by oxidative organisms appear at the surface and gradually sink back throughout the medium.



**5-8 Bacitracin susceptibility test.** Susceptibility to 0.04 units of the antibiotic bacitracin is also used to differentiate *Staphylococcus* spp. from *Micrococcus* spp. The surface of a Mueller-Hinton agar plate is inoculated with the microorganism, the disk is applied, and the plate is incubated overnight. Staphylococci are resistant to 0.04 units bacitracin (zone of inhibition less than or equal to 9 mm, organism on the right) and micrococci are susceptible (zones of inhibition greater than or equal to 10 mm, organism on the left).




**5-9 Slide coagulase test performed on a glass slide.** Cells from a fresh colony are inoculated into saline to make a smooth suspension (left), and a drop of rabbit plasma is added. If the organism possesses bound coagulase ("clumping factor"), the enzyme acts on fibrinogen in the plasma and causes clumping of the bacteria, as shown on the right. *S. aureus* is the most common pathogen among the catalase-positive, gram-positive cocci, and it is easily differentiated from other staphylococci by the coagulase test. Coagulase is a thermostable enzyme found primarily in *S. aureus*. There are two forms of coagulase: bound and free.



**5-10 The tube coagulase test** detects free coagulase. Microorganisms are incubated in plasma for 2 to 4 hours, and the tubes are turned on their sides, as shown here. Free coagulase acts on prothrombin and fibrinogen in rabbit plasma and forms a fibrin clot (tube on the left).






**9-10 Methanol red agar.** **A,** Methanol red agar differentiates *S. aureus* from *S. enteritidis* based upon gas production. The medium contains 2.5% NaCl which inhibits the growth of many nonenterogens. *S. aureus* will grow red (ferment methanol) and gas production, which causes the slanted surface to change from pink to yellow (acid) on the left of the dish. The left half of the agar plate is inoculated with coagulase negative staphylococci. Coagulase-negative staphylococci will grow on the medium, but because they do not ferment methanol the slanted red half does not change to yellow and the colonies and surrounding medium appear pink. **B,** Colonies of *Staphylococcus aureus* on methanol red agar (left up). Methanococcus fermenting colonies are turned yellow. **C,** Streak of *S. enteritidis* on methanol red agar (left up).



**9-18 Antibiotic susceptibility testing of *S. aureus*.** Test sensitivity to several strains of *S. aureus* on antibiotic discs.



**9-19 *Staphylococcus aureus* growth on methylene blue agar.** Growth on a 0.5% (w/v) methylene blue agar used to detect methanococci. The colonies are white on 2.5% NaCl. *Staphylococcus enteritidis* is a better methanococcus. Growth that occurred after 14-hour incubation at 14°C from a spot of an inoculum of BSA is shown on the top and the lack of growth of a methanococcus viable if grown on the bottom.



**0-19** Colonial morphology of *Streptococcus agalactiae* on 5% sheep blood agar. Colonies resemble those of *Streptococcus pyogenes*.



**0-20** Novobiocin susceptibility test, used to differentiate *S. agalactiae* from all *Streptococcus* on the right from other species of *Streptococcus* (e.g., *Streptococcus agalactiae* susceptible; *Streptococcus* on the left) encountered in other specimens.

## CHAPTER 6 *Streptococcaceae*

**T**he family *Streptococcaceae* and related genera include some genera encountered in the clinical laboratory: *Streptococcus*, *Abstraxia*, *Enterococcus*, *Gemella*, *Chlorostreptococcus*, *Lactococcus*, *Leuconostoc*, *Pediococcus*, and *Streptococcus*. The more commonly isolated genera are illustrated here. The streptococci are catalase-negative, gram-positive, coccoid to coccobacillary in morphology (forming pairs and chains), colonial morphology and hemolysis on 5% sheep blood agar are very helpful characteristics used for preliminary identification. Hemolysis is classified as alpha, beta or gamma.



**6-1 Alpha hemolysis on 5% sheep blood agar plate.** Alpha hemolysis is an indicator test of the type of red blood cells (RBC) causing a given organism to grow on a medium. The colonies (small dots) surrounding the colony. Various streptococci of various and numerous strains are alpha-hemolytic.



**6-2 Colony of gamma-hemolytic streptococci on 5% RBC.** Microscopic view of gamma hemolytic streptococci (small dots) on 5% RBC medium (surrounding the colony).



**6-3 Beta hemolysis on 5% sheep blood agar plate.** Beta hemolysis is a complete lysis of RBC surrounding the colony that can result in a very mucous colony.




**6-4 Colony of beta-hemolytic streptococci on 5% RBC.** Microscopic view of beta hemolytic streptococci (small dots) on 5% RBC medium (surrounding the colony) that can result in a very mucous colony.



**8-3** *Staphylococcus hemolyticus* on 5% sheep blood agar. Leucocytes are present to denote lack of hemolysis of the 5% surrounding *Staphylococcus* colonies.



**8-4** Colony of *Staphylococcus streptococcus* ( $\beta$ 102). On 5% sheep blood it is apparent hemolytic activity of this organism produced by *Staphylococcus*.



**8-7** Gram stain of streptococci in broth culture ( $\beta$ 10204). Gram stain of a positive blood culture broth demonstrating gram-positive cocci arranged in chains. Streptococci are normal organisms in the upper respiratory tract and the gastrointestinal tract. (For this species, Lance field's of specimens from the site are not helpful in diagnosing infections caused by the pathogen; streptococci such as *S. pyogenes* and *S. pneumoniae*.)



**8-8** Strain of group A streptococci on 5% sheep blood agar. Colonies of group A streptococci on 5% sheep blood agar are small and surrounded by wide zones of beta hemolysis.



**8-8 A.** Colonies of group A streptococci on 5% sheep blood agar (left) and *Streptococcus selectivus* agar (right), a primary culture medium containing group A streptococci. Observe that the colonies on the SBA are smaller in size than those on the plate on the right. Groups C, F and G can also be isolated from the upper respiratory tract, and their colonies can be confused with group A streptococci if 5% sheep blood agar is used. Localization of a white, opaque mass on this medium may indicate (1) 1-7% agar and streptococci (1-2% agar), (2) beta-hemolytic growth of group A streptococci, (3) streptococci with alpha-hemolysis, (4) streptococci and group A streptococci. **B.** Colonies of *Streptococcus pyogenes* on 5% sheep blood agar. Colonies are small (1 mm) in diameter, and are surrounded by a zone of beta-hemolysis.



**8-10 PHS test.** The presence of an autoagglutination of units that agglutinates the minute PHS (*Staphylococcus aureus*) is a characteristic property of group A streptococci (beta-hemolytic and *Streptococcus pyogenes*) and *Streptococcus pneumoniae*. A positive result is indicated by a red color (plate on the right).



**8-11 Streptococcus susceptibility test.** The large white autoagglutination is an abnormal result on the PHS test for the presumptive identification of group A beta-hemolytic streptococci. A 0.04 L hydroxyethyl piperazine sulfoxide (HEPES) buffer solution is used as the suspension on sheep blood agar. After overnight incubation at 37°C, any mass of autoagglutination is later tested as a positive test, and the autoagglutination is presumptively identified as group A streptococci by serotyping.



**6-10 Gram stain of vaginal secretions.** Group II streptococci are purple-stained cocci pairs (aggregates of streptococci) (1,200 $\times$ ). Group II streptococci *S. agalactiae* are visible in the yellowish part of smears, occasionally colonizing and adhering to mucus of cervix and Clavicle in the major cases of neonatal sepsis and meningitis. Microscopically this occasionally occurs in pairs or chains (see figure).



**6-10 Colonies of group II streptococci on 5% sheep blood agar (100 $\times$ ).** Group II streptococci colonies are larger than other beta-hemolytic streptococci, but the hemolytic zone surrounding the colony is smaller, compare with size of beta-hemolytic with Figure 6-5 on p. 39.



**6-10 Hippurate hydrolysis test.** Group II streptococci and some strains of other bacteria hydrolyze hippurate resulting in the formation of glycine and ammonia/ammonium. A suspension of the microorganism is incubated for 2 hours at 37°C in a hippurate solution and then the indole test reagent (this is added) formation of glycine. If it is present it is tested by development of a purple color (color: 10 strains give color (100 $\times$ ). A negative reaction colorless (right) contains colorless.



**6-10 CAMP Test.** A positive CAMP test is Group II streptococci demonstrating the arrowhead-shaped enhancement of beta-hemolysis that occurs when the hemolytic beta toxin produced by *S. agalactiae* (the narrow point) is mixed bottom fully opposite streptolysin O (the plasma) (the pleurogram) is in synchronicity with the CAMP factor protein produced by Group II streptococci (narrow perpendicular to the streptolysin test test zone (arrow)). The CAMP test which is sensitive to the genus *Clavicle*, *Alloca* and *Morax* *Petrus* is an alternative to hippurate hydrolysis.



**0-16** Brown stains of respiratory secretions concentrated in the posterior, hook-shaped deflexion region of a *Dipterocarpaceae* seedling.



**0-17** **A.** Colonies of *Streptomyces penicillatus* on 0% sheep blood agar. These young colonies are rod-shaped with complete edges, somewhat rounded, and about 1 mm in diameter. They are surrounded by a zone of alpha hemolysis. **B.** Colonies of *Streptomyces penicillatus* on 0% sheep blood agar. These slightly-rod colonies display the typical structure as covered by the white, leathery, verrucose typical of T genus color. **C.** Colonies of *Streptomyces penicillatus* on 0% sheep blood agar. Characteristically, T genus-like colonies may simply flatten out as they age.








A



**9-18** **A, Colonies of *Streptococcus pneumoniae* on chocolate agar.** Medium was inoculated (100  $\mu$ l) by rolling in both large areas of alpha hemolysis. **B, Colonies of *Streptococcus pneumoniae* on chocolate agar (microscopy).** Colonies are quite flat on chocolate agar.



**9-19: Optochin susceptibility test.** Colonies of *Streptococcus pneumoniae* (pyridoxylated) on optochin agar (10% hydroxyurea hydrochloride) obtained on the paper disk applied to the surface of an inoculated 15-day blood agar plate. A zone of greater than 14 mm in diameter is presumptive evidence for *S. pneumoniae pneumoniae*. The size of white zone inoculated with culture suspension. Zones of less than 14 mm in diameter are questionable and should be confirmed with the bile solubility test.



A



**9-20** **A, Bile solubility test.** Green suspension of white colonies with diameter 1 cm or within 10 minutes at 37  $^{\circ}$ C in the presence of 1% sodium hydroxide. A drop of the suspension is applied to the surface of a rich double (left) media right if the organism is *Streptococcus pneumoniae*. Colonies on the blood agar plate on the left also showed a bile solubility test. In showing that the isolate is not susceptible to *Streptococcus pneumoniae pneumoniae*. **B, Bile solubility test (control test).** *Streptococcus pneumoniae pneumoniae* is not solubilized in the bile added to the plate surface. After the bile has is opened, only the zone of hemolysis remains, as shown here.



**8-21 Antimicrobial susceptibility testing of *Streptococcus pneumoniae*.** *Streptococcus pneumoniae* is inoculated by gentamicin resistance using a 100-µl penicillin disk + 1 µg streptomycin or gentamicin with the Kirby-Bauer disk diffusion method consisting of a plastic strip impregnated with varying concentrations of the antibiotic in this case, penicillin. Following overnight incubation, the point at which the first degree of inhibition of growth of the isolate appears (the 2 µm diameter scale on the strip is best point to the nearest inhibitory concentration point)



**8-22 Gram stain of a viridans streptococcus species in blood culture broth (1:10000).** There are numerous species of the viridans streptococci group which are normal mucosal microflora in mammals. Only a few species characterize the group and, each look like the majority they usually appear in long chains or pairs, microaerophilic from a blood culture broth medium.



**8-23 Optochin susceptibility test performed on a species of viridans streptococci.** Colonies of viridans streptococci are resistant to optochin (disc-inhibiting these organisms) concentrated organisms from *Streptococcus pneumoniae* (shown in figure 8-17)



**8-24 Identification of viridans streptococci with carbohydrate fermentation reactions.** Culture for identification requires every streptococci isolating 10-20 cells optochin fermentation, 6-10% NaCl serum ferment redness within antibiotic tests success and results (shown here). Other characteristics such as serum Voges-Proskauer reaction for lactate dehydrogenase release and tolerance to grow at 10% and 45°C will help identify the strain (significant results)



**8-26** Colonies of *Petrococcus* spp. on 8% sheep blood agar. *Petrococcus* spp. resemble *Staphylococcus* or other catalase-positive and coagulase-negative staphylococci in their colonial morphology and susceptibility for growth on many culture media. They are, in fact, isolated from other sites besides the nose, and recognized by their response to coagulase (NI) and DNA (see Figure 8-16) and gelatinase (A) test results. They are also indole-positive but PYR negative.



**8-27** *Streptococcus salivarius* on 5% sheep blood agar plate. In a family, *Streptococcus* is the most common organism isolated by the laboratory; some strains of various streptococci, *Lactococcus*, and *Lactobacillus* spp.



**8-28** MacConkey susceptibility test. *Proteus mirabilis* and *Proteus vulgaris* are distinguished from other *Staphylococcus* organisms by their resistance to macConkey. There is no growth of bacteria across the line, demonstrating that placed on an inoculation of *Proteus* on a 96-plate with a plate containing a range of isolates surrounding the various susceptibility test plate.



**8-29** 48-hour *Staphylococcus aureus* suspension broth and the addition of a toxin plate on the presence of toxin. The toxin is a dark, vertical band in a tube containing the agar which contains 0.2% (w/v) of a 10% (w/v) solution of sodium hydroxide. The toxin is a dark, vertical band in a tube containing the agar which contains 0.2% (w/v) of a 10% (w/v) solution of sodium hydroxide. The toxin is a dark, vertical band in a tube containing the agar which contains 0.2% (w/v) of a 10% (w/v) solution of sodium hydroxide.



**8-20** Colonies of *Enterococcus* spp. on 0% sheep blood agar plate. The colonies are round, white to gray-white, spherical (2 to 3  $\mu\text{m}$  in diameter), and are usually nonhemolytic.



**8-21** Colonies of *Lactococcus* spp. on 0% sheep blood agar plate. Usually, the colonies of *L. lactis* spp. (streptococcal subgroup 1-12) upon growth on blood agar (48 to 72 hours) reveal alpha hemolysis (opacity). However, these also compare in size to those produced from *Enterococcus* spp. by their own size (2 to 3  $\mu\text{m}$ ) and by the production of gamma and/or theta toxin (beta-hemolytic *Enterococcus faecalis*).



**8-22** 60% streptococci (left) and 0.1% (right) *Streptococcus*. The tube reaction of 60% (left) (dark) is indicating that the streptococci grow in the presence of the acid but do not produce gas. Growth in a broth containing a 0.1% concentration of 0.1% (left) after 24 hours is indicated by turbidity and a change in the indicator from pink to yellow following overnight incubation. The reactions show a test similar to the identification of *Enterococcus* spp. A positive PYR test (Figure 8-18) may also be used to confirm the identification of non-pathogenic *Enterococcus* spp. (enterococcal group II streptococci) (left) but unable to grow in the presence of 0.1% NaCl broth.

## CHAPTER 7 *Aerobic Gram-Positive Bacilli and Actinomyces spp.*

**C**ommon genera of aerobic gram-positive bacilli include *Bacillus*, *Corynebacterium*, *Erpobacterium*, *Clostridium* (gram variable), *Lactobacillus*, and *Listeria*. More clinical isolates include species of *Armadillibacterium*, *Citrobacter*, and *Phlebotomus*. Also included in this section are *Nocardia* and *Streptomyces* genera consisting of filamentous branching gram-positive bacilli.

Preliminary differentiation among these genera is based on microscopic morphology and the catalase reaction. *B. thur* spp., *Corynebacterium* spp., *Listeria*, *Clostridium*, and *Phlebotomus* spp. are catalase-positive; of these, *Bacillus* spp. is the only spore former. More examples of the macroscopic morphology of the gram-positive bacilli are presented in the Gram stain section.

**7-1 Gram stain of *Bacillus* spp.** Large gram-positive bacilli with squared-off ends measuring 0.8 × 7.0 μm occurring singly and in chains (FIG. 7-1) may appear gram-variable. Some species demonstrate subspore-like morphology located centrally or laterally along the axis. The microscopic morphology may resemble that of *Clostridium* spp.



**7-2 Colonies of *Bacillus* spp. on 5% sheep blood agar.** Colonies of *B. thur* are usually round and surrounded by a zone of beta hemolysis. Colonies on case agar are rough and dry with shaggy sides.





**FIG 2-5** Growth of *Staphylococcus aureus* on 5% sheep blood agar. Some species of bacteria are catalase/oxidase and fermenting bacteria, such as the *Staphylococcus aureus* and *Paraflocculus* spp. *Staphylococcus aureus* and other Gram + cocci are usually catalase/oxidase. They grow in air environments. An ability that helps create living colonies.



**FIG 2-6** Growth of *Staphylococcus aureus* on egg yolk agar. This test on egg yolk agar's ability to demonstrate rapid to identify *Staphylococcus aureus*. *Staphylococcus aureus* is one of many species that produce the enzyme lecithinase, demonstrated by a zone of opacity (opacity) in the agar extending away from the bacterial growth. In this case, the left. *Staphylococcus aureus* can also be identified on egg yolk agar demonstrated by a zone of opacity (opacity) in the center of the bacterial growth (opacity) surrounding the egg yolk, as shown on the right.



**FIG 2-7** Growth of *Staphylococcus aureus* on 5% sheep blood agar. Many *Staphylococcus aureus* proteolytic enzymes that can hydrolyze gelatin. In the center of the tube in this figure, exposed, undissolved, a very thin layer is the characteristic detection of gelatinase activity. By doing so, the tube is gelatin, leaving only the clear (transparent) that allows the ring was broken in the organism's expansion.



**FIG 2-8** Growth of *Listeria monocytogenes* on 5% sheep blood agar. Colonies of *Listeria monocytogenes* produce small (less than 1 mm in diameter) smooth, convex, but not translucent. The colonies are surrounded by a very characteristic narrow zone of beta hemolysis on 5% sheep blood agar and may be confused with group II beta-hemolytic streptococci. Differentiating characteristics include *Listeria*'s narrow zone of beta hemolytic, gelatinase, casein production and growth on bile and hydrolysis of novobiocin. Test species by doing a tug test.



**7-7** *Staphylococcus hyalei* catalase test. Most strains of *Citrobacter amara* (group 1) can hydrolyze urea within 2 hours. This photograph shows positive gas from the 100-minute test at 10°C and negative (right). These tests are described in Figure 6-14 and p. 46 in *Microbiology: An Atlas*, 14e.



**7-8** Gram stain of *Staphylococcus epidermidis* (strong, small gram-positive bacilli) and *Staphylococcus aureus* (F and L forms resembling Chinese letters). Most of the *Staphylococcus* spp. isolated in the clinical laboratory are normal flora of the skin and mucous membranes. It is difficult to differentiate the staphylococcal species from the pathogenic *S. aureus* and *S. epidermidis* on the basis of the Gram stain.




**7-9** Colonies of *S. aureus* on 5% sheep blood agar. Colonies are smooth and a little fluffy, showing enough adhesion to 5% to 10% L<sub>1</sub>. This species can cause life-threatening sepsis.




**7-10** Colonies of *C. difficile* on 5% sheep blood agar. *Clostridium difficile* can grow well on 5% sheep blood agar as relatively opaque colonies.



**7-14** Culture of *E. coli* on Tryptone agar. A white line medium should be used along with tryptic blood agar medium (propagation as *E. coli* is common as well as in different *E. coli* strains like other commensal ones). On the selective media in tryptic blood agar or on a slant culture (Tryptone agar) the colonies of *Escherichian coli* have a grey black appearance.



**7-15** Culture of *Corynebacterium* spp. from *Corynebacterium* spp. (old *Corynebacterium* spp. is called *Corynebacterium* spp.) on tryptic blood agar. Small colonies and can easily be confused with colonies of streptococci, especially when isolated from the respiratory tract. Some species appear white and appear resembling *Streptococcus* morphology.



**7-16** Lactifer-methylene blue stain. 1. *Lactobacillus* streptococci are made granules when stained with the Lactifer-methylene blue stain. The stain is best performed on reference strains on a Lactifer agar stain. *Methylenblau* deposits on various groups of *Lactobacillus* like stain.




**7-16** Gram stain of *Lactobacillus* spp. In weak and quantitative, *Corynebacterium* spp. appear as long non-spore forming gram positive bacilli. They stain blue-purple in short and photophobic cocci-like. When they appear microscopically *Corynebacterium*. They can be identified as streptococci. *Lactobacillus* are normal flora of the upper respiratory tract, intestinal tract, and vagina, although they often fail to grow in culture.

**7-18** Culture of *Lactobacillus* spp. Culture on alpha hemolytic, surrounding culture (streptococci) and *Corynebacterium*. It may be necessary to differentiate these groups depending on the species source. In most cases, both hemolytic and streptococci are normal flora and do not require identification. They are usually isolated from streptococci. Growth on TM agar will differentiate *Corynebacterium* like bacilli which resemble to gram.





**7-16** Colonies of *Corynebacterium rubeum*. Colonies are usually small (less than 1.0 mm in diameter) and inconspicuous. A general color reaction of brownish red prevails for the mature stage.



**7-18** Gram stain of *Corynebacterium rubeum* (Gram positive). *Corynebacterium rubeum* is a small, gram-positive to gram-variable rod-like microorganism that can be isolated from the genitourinary tract of humans. This organism was at one time included in the genus *Corynebacterium* and *Corynebacterium*.



**7-19** Colonies of *Corynebacterium rubeum* on phosphate agar. After 24 hours of incubation, colonies of *Corynebacterium rubeum* appear as tiny, pinpoint colonies on phosphate agar. Production of the characteristic pink to orange color for the diagnosis of bacterial vaginosis occurs in more than 90% of women who are colonized asymptotically with *C. rubeum*.



**7-17** *Corynebacterium rubeum* on triple sugar iron (TSI) agar. *Corynebacterium rubeum* produces H<sub>2</sub>S precipitates along the surface in a TSI agar that settles to the bottom of the tube at 37°C. This is a helpful feature that separates this organism from other gram-positive bacilli.



**7-20** *Staphylococcus hyalae* hydrolysis test. *S. hyalae* hydrolyzes gelatin. *Staphylococcus hyalae* (the microorganism hydrolyzes T<sub>2</sub> agar) is a rod-shaped bacterium that is isolated at 37°C as described in Figure 6-14 on p. 40. It produces a red color in the first segment of a gelatin tube within 10 minutes after the addition of catalyase. A negative reaction is colorless.



A





B

**FIGURE 1.** Colonies of *Aeromonas hydrophila* (left) and *Aeromonas caviae* (right). Colonies of *Aeromonas hydrophila* are on the left and *A. caviae* on the right; both fermentative after 48 hours, and can be confused with streptococci when isolated from a flow culture. As shown here on the left, the colonies are anaerobic after overnight incubation compared with the strong beta-hemolysis of group A streptococci. *Aeromonas* enter was previously included in the genus *Corynebacterium* but is a similar species. Microscopically, these protospiral bacilli resemble *Corynebacterium* spp. **B.** Colonies of *Aeromonas hydrophila* (left) and *Aeromonas caviae* (right) on a dark agar surface.



**FIGURE 2.** Colonies of *Rhodospirillum rubrum* on 8% sheep blood agar. *Rhodospirillum rubrum* is a red to pinkish "nonfermenter" genus. *Rhodospirillum rubrum* was also previously included in the genus *Corynebacterium* as well as the genus *Mycobacterium*. The microorganism grows well on sheep blood agar and reduces iron-sulfide (Sulfide) but it does not grow on MacConkey agar. Colonies usually develop a pink pigmentation of the prolonged incubation.



**FIG. 2-23** A, B, *Stenotrophomonas* appearance of *Stenotrophomonas* sp. (strain 10). Young cells from culture medium 24 hours after inoculation (A) usually first appear rodlike. B, Cellular morphology changes to a more rod-like form within 24 to 48 hours. Many forms may enter the filamentous, beaded filamentous stage fragment into rods and cocci during the growth cycle. C, Small heat-stable, glycogen-rich bodies. *Stenotrophomonas* sp. may be partially acid fast when stained with a modified Kinyoun stain using 2% NaOH as the decolorizing agent.




**FIG. 2-24** Colonies of *Stenotrophomonas* sp. on 5% sheep blood agar. Colonies are small, glistening filaments that spread into white films. This microorganism can be grossly identified by its characteristic morphology and by hydrolysis of casein and yellow pigment on sheep blood agar.



**FIG. 2-25** Very young colonies of *Stenotrophomonas* on 5% sheep blood agar (strain 10). Colonies are usually believed to be unusual corynebacteria because of their lack of Gram stain control (CCO) group A. They are characterized by colonies with a filamentous appearance that is observed with a light microscope under low power.



**FIG. 1001** Gram-stained *Mycobacterium* spp. (1000 $\times$ ). Most diposules are branching beaded. (Numerous granules per the beaded appearance. 4  $\mu$ m scale bar.) They can also appear as rounded or coccobacillary forms that also have beaded chains resembling the dimer because of the heat-resistant properties of their outermost cell wall outer membrane.



**FIG. 1002** Gram-stained *Mycobacterium neoaurum* (1000 $\times$ ). *Mycobacterium neoaurum* is a diposule as it has an acid-fast cell wall modified lipopolysaccharide (LPS) in the cell wall. The beaded appearance helps to distinguish the microorganism from other acid-fast organisms.



**FIG. 1003** Colonies of *Mycobacterium neoaurum* on 5% sheep blood agar. Colonies usually appear as confluent and as they mature they can develop a number of yellow-ochre to a chalky white to yellow-orange but not purple from a red tint. Most colonies have a white color.



**FIG. 1004** Colonies of *Mycobacterium neoaurum* on 50% egg yolk agar. Yellow-orange colonies of *Mycobacterium neoaurum*. The yellow color is more common with colonies of most *Mycobacterium* spp.



**FIG. 1005** Colonies of *Mycobacterium neoaurum* growing on 50% egg yolk (100 $\times$ ). Cells are as large as 40  $\mu$ m and frayed after 4 days incubation at 27°C. (Microphotography has been taken of a non-glucose energy form.)




**2-21 Identification of the genus *Bacillus* with biochemical reactions.** *Bacillus* spp. are anaerobic, sporulated rod-shaped bacteria. In this tube on left, hydrolysis was (dark red color), and gas in the medium of a fresh culture. Both *B. anthracis* and *B. thuringiensis* react in this way.




**2-22 Amino acid hydrolysis reactions of *Bacillus anthracis*.** Hydrolysis in the system level is usually made on the basis of hydrolysis of amino acids: cystine, tyrosine, methionine and tryptophane. Digestion of the agar around the colony is interpreted as a positive reaction. Both *B. anthracis* and *B. thuringiensis* give identical results on this plate. *B. anthracis* and *B. thuringiensis* (top and bottom quadrants) are negative.




**2-23 Amino acid hydrolysis reactions of *Bacillus anthracis*.** All amino acid reactions are negative on this amino acid hydrolysis plate. *B. anthracis* and *B. thuringiensis* are by themselves.




**2-24 Amino acid hydrolysis reactions of *Streptococcus* spp.** All reactions are positive on this amino acid hydrolysis plate. Unlike *B. anthracis*, *Streptococcus* spp. hydrolyze cystine, tyrosine, methionine and tryptophane. Catalase morphology resembles the double disc-wise colonies of *Bacillus* spp. Although *Streptococcus* spp. are Gram-positive they are more negative, distinguishing the genus from *Bacillus* spp.




**T-26 Green spore of *Streptomyces* spp. (1000x).** *Streptomyces* spp. are characterized as ellipsoidal, 1.5 by 2 micrometers long, beaded, gram positive filamentous bacilli indistinguishable from *Streptomyces* spp. Unlike *Streptomyces* spp. the parent and first spore is orange.



**T-28 Brownish multicellular silver spore of *Actinomyces* spp. (1000x).** The spore was prepared from a *Streptomyces* spore. In the center is a characteristic "silver granule" which consists of a granular, succinylated, uncolored to granular material. The *Actinomyces* spp. are gram-positive bacilli that vary in size from short, diplostriched forms to long, branching filaments. Unlike *Streptomyces* spp. they are not partially visible.



**T-27 Colonies of *Actinomyces* spp. (10x).** These are the characteristic "white to pink" colonies of *Actinomyces* spp. located at the edge of a 1 by 2 cm agar plate, seen, held under aseptic conditions. These colonies are large and raised with irregular margins.



**T-29 Chemical identification of *Actinomyces* spp. (10x).** With *Actinomyces* spp. culture strains tested. The 1st colored liquid is (1) including *Actinomyces* (small). Additionally, *Actinomyces* (small) hydrolyzes casein (2nd color, orange) but not urea (yellow color in center).

## CHAPTER 8 *Enterobacteriaceae*


**T**he *Enterobacteriaceae* are the most common bacterial isolates encountered in the diagnostic laboratory. They can be isolated from a wide variety of specimen sites as normal flora (predominantly in the intestinal tract) or pathogens. The gram-negative bacillary microscopic morphology is similar among all species. Colonies are usually large and gray on 5% sheep blood agar. Rapid lactose fermenters appear pink on MacConkey agar and have a greenish metallic sheen on Etest methylcellulose blue (EMB). The *Enterobacteriaceae* are very active biochemically and can be identified with a wide variety of tests. Four characteristics of the *Enterobacteriaceae* are growth on MacConkey agar, oxidase-negative, glucose fermentation, and reduction of nitrate to nitrite.



**FIGURE 8-1** *Enterobacteriaceae* on 5% sheep blood agar. *E. coli* colonies (sheep blood on 5% sheep blood agar showing large, full gray non-hemolytic colonies) form pink or colorless colonies instead of any color.



**FIGURE 8-2** *Proteus* species on 5% sheep blood agar. *Proteus* spp. are reported as a strain on the gram, from the large and colorless streak line often extending in waves. This characteristic of *Proteus* spp. is called an arrow and suggests that the microorganisms move by means of flagella.



A




B




C

**8-3 A, *E. coli* on MacConkey agar.** Rapid lactose fermenting strains of *E. coli* appear as shiny pink colonies on MacConkey agar. **B, *E. coli* on MacConkey agar (close-up).** **C, *E. coli* on 5% sheep blood agar (close-up).** Colonies are shiny, opaque, cream-colored, and attain 2 to 4 mm diameter overnight.




**8-4 *Klebsiella pneumoniae* on MacConkey agar.** Rapid lactose fermenting colonies of *Klebsiella pneumoniae* appears pink, large, glistening, and mucoid. This strain is probably encapsulated and therefore appears mucoid. Although this appearance is associated with *Klebsiella pneumoniae*, it is not unique for that species.




**8-5 Pigmented *Serratia* sp. on MacConkey agar.** These colonies appear red and should not be confused with the pink color due to lactose fermentation shown in 8-3. Rare strains of *Serratia* spp. produce pigment, which is seen on all solid media including the blood agar plate.






**8-6 Enterobacteriaceae on MacConkey agar.** All fast lactose fermenters and lactose nonfermenters of the Enterobacteriaceae appear as colonies like the isolated colonies on MacConkey. The colonies of the late lactose fermenting organisms appear slightly pink after 48 to 72 hours of incubation.



**8-7 Enterobacteriaceae on Hektoen enteric agar (HE).** Most species of Enterobacteriaceae except for the rapid lactose fermenters including some Enterobius, Klebsiella, Enterobacter and Citrobacter species appear colorless on either MacConkey or even weakly blue like (HE) agar. Food spoilage fungi should be incubated/checked over selective medium such as HE agar or yellow prime desferrioxime (YLD) agar to help differentiate among occasional pathogens. Fermentation of lactose occurs, and/or which results in yellow or pinkish-colored colonies due to the pH change on the bromothymol blue and acid fastness indicators.




**8-8 Salmonella spp and other Enterobacteriaceae on Hektoen enteric agar.** If an organism does not ferment lactose, acetoin, or indole the colonies will appear green color of the medium. If the microorganism produces H<sub>2</sub>S as important characteristic of the isolated agar, the colonies will have black centers and will appear completely black with prolonged incubation. These black colonies are characteristic of most Salmonella spp on HE agar.



**8-9 Identification of the Enterobacteriaceae by the spot-test level method.** Spot tests are now modified for the rapid, presumptive identification of many nonfermentation groups. For this method, blue paper or a paper disk is impregnated with the substrate. To perform the test one or more colonies are selected with a sterile wooden applicator stick from the described medium and applied to the spot surface.



**8-10 Gelatin spot test.** Some microorganisms possess other cytochrome oxidase or cytochrome oxidase which catalyze the transport of electrons from donor compounds (NADH) to electron acceptors ( $O_2$ ). In this test, the substrate phosphorylase/kinase dihydroxybenzoyl serine as artificial electron acceptor for the enzyme oxidase. The dye is reduced and forms the reduced cytochrome oxidase blue. For this test, filter paper is saturated with hydroxybenzoyl serine. Colonies of the microorganisms to be tested are placed onto the filter paper with a sterile swab, replicate streak. An immediate color change from deep blue indicates positive test result (as in the left). The color change from white to purple and thereafter the transition should not change the color of the paper (as in the right).



**8-11 Methyl red test ( $pH < 4.5$ ).** *Escherichia coli* produces formic acid. The enzyme  $\beta$ -galactosidase hydrolyzes lactose from lactose. Rapid lactose fermenter (however possess a permease) enzyme that quickly the substrate yielding pink color in Methyl red test. Production of  $\beta$ -galactosidase only results in delayed lactose fermentation (2 to 12 days). The permease ferments perme without enzyme. The  $pH < 4.5$  is typical for  $\beta$ -galactosidase.  $\beta$ -galactosidase hydrolyzes from the substrate-DMFC into water-soluble D-galactose compound and glucose. A yellow color is formed in methyl red ( $pH < 4.5$ ) and no color change (negative) indicates no color change from white ( $pH > 6.8$ ).



**8-12 Indole spot test.** Indole is one of the degradation products of the amino acid tryptophan. Indole produced by a microorganism grows by reaction that is catalyzed by tryptophanase can be detected by its ability to combine with certain reagents to form colored compound. For this spot test, filter paper is coated with other reagent such as indoleacetylacetate/indole or paraformaldehyde/indole/ethylglyoxal (Kovacs reagent). Colonies of the microorganisms to be tested are rubbed onto the filter paper with a sterile swab or applicator stick. If paraformaldehyde/indole/ethylglyoxal is used, a fluorescent color appears immediately in the presence of indole. If Kovacs reagent is used, a bright pink color develops. In this figure using filter paper coated with paraformaldehyde/indole/ethylglyoxal the colonies on the right is positive while the one on the left is negative. The pink color is the result of a lactose-positive color referred from Methyl red test.



**8-13 Methyl red-Voges-Proskauer (Methyl red test).** Both tests are performed from the same inoculum were prepared, which is divided by testing. The methyl red test is used to determine the pH of the end products of glucose fermentation. In converting glucose to pyruvate acid, some enterobacteria produce acid and produce a pH less than 4.4 while others produce the methyl red color red, while others produce colorless pH greater than 6.0, which is in the indicator pH low. The Voges-Proskauer test determines the presence of acetyl and ketone is added to the broth suspension. Ketone is oxidized to diacetyl, yielding a red color pH greater than 6.0. The left set of tubes shows a positive MR and a negative VP (no color change), and the right set shows a negative MR, and a positive VP.



**8-14 Nitrate utilization test.** Success of the nitrate reduction test is the ability to reduce nitrate to a color change of carbon. In this test, a colony is inoculated into the surface of a nitrate broth in the medium is incubated overnight at 37°C. The products contain nitrate (blue), nitrite (red), and ferrous sulfate (blue). The presence of a blue color indicates the presence of nitrate and products and a positive nitrate test (left tube). If the test is negative (blue is no color change right tube).

**8-15 Phosphatase activity test (PAP) test.** This test is based on the principle that some microorganisms possess an enzyme phosphatase that breaks down phosphatase resulting in the production of phosphoric acid. A substrate medium growing on a phosphatase agar plate is tested by adding a few drops of a 10% ferric chloride solution. A green color (left tube) indicates a positive test result. This color remains colorless with a negative result (right tube). This test is helpful in differentiating the coliform (PAP-positive) from the other Enterobacteriaceae (PAP-negative).






**Q18: Grams test.** Microorganisms that possess strong gram negative hydrophobic tails (lipopolysaccharide and O<sub>2</sub>) do not take color better (they are pink) due to the presence of the alcohol in the dye mixture. This type of stain is called Gram stain. Gram stain is used to identify bacteria. The left and middle tubes are positive. The left tube was inoculated with *Proteus mirabilis* (a gram negative bacterium) and the middle tube was inoculated with *Staphylococcus aureus* (a gram positive bacterium) in 24 hours. The right tube was inoculated with *Staphylococcus aureus* (a gram-positive microorganism).

**Q19: Positive Fall and gram tests.** Four species of the *Escherichia coli* group are EHEC and were positive. *Proteus mirabilis*, *Proteus vulgaris*, *Morganella morganii* and *Providencia stuartii* (fall) were tested and produced false reactions and also caused antibiotic resistance or could be a *Staphylococcus*. A spot culture test could identify the species. *Proteus mirabilis* is gram-negative and *Proteus vulgaris* is gram-positive. The clinical test would differentiate the microorganisms. *Morganella morganii* is gram-negative and *Providencia stuartii* is gram-positive.





**FIGURE 10-10** *Bacillus cereus* decarboxylates lysine. The decarboxylase enzyme removes  $\alpha$ -amino groups from the amino acid lysine and urethane. The decarboxylase enzyme removes carbonic group from urethane di-carboxyl base within the amino acid (left tube) and urea containing (right tube) plus the amino acid substrate are liberated. Decarboxylation and decarboxylation are reversible reactions. In the case decarboxylation is involved with urea, several amino acids are, notably all lysine turn yellow during acidification of the test tube (dominant purple) by the acid end-product of glucose fermentation. If the amino acid decarboxylated the  $\alpha$ -amino acid product is available substrate to convert to an alkaline pH (purple). In this assay, negative (yellow) result from left and positive (red) color are positive and lysine is negative.



**FIGURE 10-11** Triple Sugar Iron (TSI) agar slant and butt tubes. The agar slants fermentations of glucose, lactose and sucrose. H<sub>2</sub>S production and gas formation. The medium contains one part glucose in 10 parts lactose and sucrose for an anaerobic slant and placed in tubes if the red cross-section ferments glucose since the slant tube (slant and butt) will turn red (yellow) because fermentation will begin catabolic process on the agar surface resulting in alkaline and products and a pink slant if lactose and/or sucrose are fermented, the entire tube remains yellow, because of the greater concentration of these carbohydrates. H<sub>2</sub>S reacts with ferrous sulfamate to form a positive black and product, ferrous sulfide. Gas formation is detected by the presence of bubbles or cracks in the agar. The *left tube* ferred *Salmonella typhi* (Gross) (K1) is a non-fermenter and identification system. Fermentations can occur in the upper portion of the first tube, glucose and lactose fermentation, H<sub>2</sub>S production, gas formation, and phenylalanine decarboxylation. Lysine decarboxylation occurs in the lower (anaerobic) portion. In this example, glucose is fermented, and H<sub>2</sub>S is produced in the TSI (left tube). The yellow color (acidic pH) are demonstrated, is indicated by the black end-product of H<sub>2</sub>S production. The same reaction appears in the upper portion of the *right tube* ferred a positive in the lower portion. This reaction suggests *Salmonella typhi*.



**8-20 TSI slant reactions.** The 10 TSI tubes have the following reactions\* (left to right):

	Slant	Butt	H <sub>2</sub> S	Gas formation
#1	acid	acid	0	0
#2	acid	acid	0	+
#3	alkaline	acid	0	+
#4	alkaline	acid	0	0
#5	alkaline	acid	0	0
#6	alkaline	acid	0	0
#7	alkaline	acid	0	+
#8	acid	acid	+	0
#9	alkaline	acid	slight	0
#10	alkaline	acid	+	+


\*Acid = yellow; alkaline = red; H<sub>2</sub>S-positive = black; gas = cracks, bubbles.



**8-21 r/b<sub>1</sub> reactions.** The eight r/b<sub>1</sub> tubes have the following reactions\*:

	Slant	Butt	H <sub>2</sub> S	Gas formation	PAD	Lysine
#1	acid	acid	0	+	0	0
#2	alkaline	acid	0	0	0	0
#3	alkaline	acid	0	+	0	0
#4	alkaline	acid	0	0	0	+
#5	alkaline	acid	0	0	+	0
#6	alkaline	acid	+	0	+	0
#7	alkaline	acid	slight	0	0	+
#8	alkaline	acid	+	+	0	+

\*Acid = yellow; alkaline = red; H<sub>2</sub>S-positive = black; gas = cracks, bubbles; PAD positive = brown; lysine-positive = red.



**8-22 r/b<sub>2</sub> medium and reactions.** This tube has an upper (aerobic) and lower (anaerobic) portion. Indole is interpreted in the upper portion and ornithine decarboxylase and motility are interpreted in the lower portion. This tube is similar to the conventional MIO (motility-indole-ornithine) medium.

The four r/b<sub>2</sub> tubes shown have the following reactions\*:

	Motility	Indole	Ornithine
#1	+	+	+
#2	0	+	+
#3	+	0	+
#4	0	+	0

\*Indole-positive = red (+) at the surface of the top portion, ornithine positive = purple red, motility positive = turbid medium with a hazy stab line

**8-23 CIT/RHAM medium and reactions.** The upper (aerobic) portion of this medium is for demonstration of citrate utilization, and the lower (anaerobic) portion is to detect rhamnose fermentation.

The four CIT/RHAM tubes have the following reactions\*:

	Citrate	Rhamnose
#1	0	0
#2	0	+
#3	+	+
#4	+	0

\*Citrate-positive = blue; citrate-negative = green; rhamnose-positive = yellow; rhamnose-negative = green.





**8-24 Characteristic reactions of *Salmonella* spp. (except *Salmonella typhi*)**

TSI	r/b <sub>1</sub>	r/b <sub>2</sub>	CIT/RHAM
alk/acid	alk/acid	ornithine +	+/+
H <sub>2</sub> S+	H <sub>2</sub> S + lysine + PAD 0	motile indole 0	

**8-25 Characteristic reactions of *Salmonella typhi*.**

TSI	r/b <sub>1</sub>	r/b <sub>2</sub>	CIT/RHAM
alk/acid slight H <sub>2</sub> S +	alk/acid slight H <sub>2</sub> S + lysine + PAD 0	ornithine 0 motile indole 0	0/0







**8-26** Characteristic reactions of *Proteus vulgaris*.

TSI	r/b <sub>1</sub>	r/b <sub>2</sub>	CIT/RHAM
acid/acid	alk/acid	ornithine 0	+ / 0
H <sub>2</sub> S +	H <sub>2</sub> S +	motile	
	PAD +	Indole +	
	lysine 0		

**8-27** Characteristic reactions of *Yersinia enterocolitica*.

TSI	r/b <sub>1</sub>	r/b <sub>2</sub>	CIT/RHAM
alk/acid	alk/acid	ornithine +	0 / +
H <sub>2</sub> S 0	lysine 0	non motile	
	H <sub>2</sub> S 0	Indole +	
	PAD 0		





**8-28** Characteristic reactions of *Shigella* spp.

TSI	r/b <sub>1</sub>	r/b <sub>2</sub>	CIT/RHAM
alk/acid	alk/acid	ornithine 0	0/0
H <sub>2</sub> S 0	lysine 0	nonmotile	
	H <sub>2</sub> S 0	indole 0	
	PAD 0		



**8-29** Characteristic reactions of *E. coli* using MICRO-ID strip. The MICRO-ID (Organon Teknika, Durham, N.C.) strip contains 15 biochemical tests for the rapid identification (4 hours) of the *Enterobacteriaceae*. The system is based on the principle that the heavy inoculum suspension of the organism to be tested contains high levels of preformed enzymes. The tests included in the system are Voges-Proskauer (VP), nitrate (N), phenylalanine deaminase (PD), H<sub>2</sub>S, indole (I), ornithine (OD) and lysine decarboxylase (LD), malonate (M), urea (U), esculin (E), ONPG, arabinose (ARAB), adonitol (ADON), inositol (INOS), and sorbitol (SORB). The test results are interpreted visually, coded numerically, and compared to identification tables or the MICRO-ID profile code book to make the identification. The interpretation of this example is:

VP	N	PD	H <sub>2</sub> S	I	OD	LD	M	U	E	ONPG	ARAB	ADON	INOS	SORB
0	+	0	0	+	+	+	0	0	0	+	+	0	0	0

The identification is *E. coli* and the profile code number is 23430. The code was derived by dividing the tests into 5 groups of three and assigning a score to each. The first test of each set is given a score of 4, the second test is assigned a 2 and the third test is assigned a 1. If any of the tests are positive, it is given its score, but it receives a 0 if it is negative. The numbers in a set can range from 0 to 7.



**B-30 Characteristic reactions of *E. coli* using the API 20E strip.** The API 20E strip (bioMérieux Vitek, Hazelwood, Mo.) is a self-contained system of 20 microtubes of dehydrated substrates, a miniaturized version of conventional procedures, designed for overnight incubation. Identification is made by adding necessary reagents and then visually interpreting the results. The tests included in the system are: ONPG, arginine dihydrolase (ADC), lysine (LDC) and ornithine (ODC) decarboxylase, citrate (CIT),  $H_2S$ , urea (URE), tryptophan deaminase (TDA), indole (IND), Voges-Proskauer (VP), gelatin (GEL), glucose (GLU), mannitol (MAN), inositol (INO), sorbitol (SOR), rhamnose (RHA), sucrose (SAC), melibiose (MEL), amygdalin (AMY), and arabinose (ARA). An oxidase test must be performed separately. Numerical coding of results allows computerized interpretation of patterns, lists of which are available in a codebook or in computerized form. The interpretation of this example is:

ONPG	ADC	LDC	ODC	CIT	$H_2S$	URE	TDA	IND	VP	GEL	GLU	MAN	INO	SOR	RHA	SAC	MEL	AMY	ARA
+	0	+	+	0	0	0	0	+	0	0	+	+	0	+	+	0	+	0	+

The identification is *E. coli* and the profile code number is 5144552. The code was derived by dividing the tests into 7 groups of three (oxidase reaction is the third test in the last set) and assigning a score to each. The first test of each set is given a 1, the second test is assigned a 2 and the third test is assigned a 4. If any of the tests are positive, it is given its score and receives a 0 if it is negative. The total scores in a set can range from 0 to 7.



**B-31 Characteristic reactions of *Enterobacter aerogenes* (top) and *Enterobacter cloacae* (bottom) using API 20E strips.** The interpretation of these tests in this example is:

	ONPG	ADC	LDC	ODC	CIT	$H_2S$	URE	TDA	IND	VP	GEL	GLU	MAN	INO	SOR	RHA	SAC	MEL	AMY	ARA
E clo	+	+	0	+	+	0	0	0	0	+	0	+	+	0	+	+	+	+	+	+
E aer	+	0	+	+	+	0	0	0	0	+	0	+	+	+	+	+	+	+	+	+

The profile code number for *E. cloacae* (E clo) is 3305573.

The profile code number for *E. aerogenes* (E aer) is 5305773.

The major differences between the two species are the reactions with arginine, lysine, and inositol.



**B-32** Characteristic reactions of *Proteus mirabilis* (top) and *Proteus vulgaris* (bottom) using API 20 E strips. The interpretation of these tests in this example is:

	ONPG	ADC	LDC	ODC	CIT	H <sub>2</sub> S	URE	TDA	IND	V P	GEL	GLU	MAN	INO	SOR	RHA	SAC	MEL	AMY	ARA		
P vul	0	0	0	0	+	+	+	+	+	0	P vul	0	+	0	0	0	0	+	0	+	0	
P mir	0	0	0	+	+	+	+	+	0	0	P mir	0	+	0	0	0	0	0	0	0	0	0

The profile code number for *P. vulgaris* (P vul) is 0674021.

The profile code number for *P. mirabilis* (P mir) is 0734000.

The differences between the two species are the reactions with ornithine, indole, sucrose, and amygdalin.



**B-33** Characteristic reactions of *Serratia marcescens* (top), *Enterobacter aerogenes* (middle), and *Klebsiella oxytoca* (bottom) using API 20 E strips.

	ONPG	ADC	LDC	ODC	CIT	H <sub>2</sub> S	URE	TDA	IND	V P	GEL	GLU	MAN	INO	SOR	RHA	SAC	MEL	AMY	ARA
S mar	+	0	+	+	+	0	0	0	0	+	S mar	+	+	+	+	0	+	0	+	+
E aer	+	0	+	+	+	0	0	0	0	+	E aer	0	+	+	+	+	+	+	+	+
K oxy	+	0	+	0	0	0	0	0	+	+	K oxy	0	+	+	+	+	+	+	+	+

The profile code number for *Serratia marcescens* (S mar) is 5307723.

The profile code number for *E. aerogenes* (E aer) is 5305773.

The profile code number for *Klebsiella oxytoca* (K oxy) is 5245773.

The differences between *S. marcescens* and *E. aerogenes* are the reactions with gelatin, rhamnose, and melibiose.

The differences between *K. oxytoca* and *E. aerogenes* are the reactions with ornithine and indole.



**8-34** Characteristic reactions of *Yersinia enterocolitica* and *Shigella sonnei* using API 20 E strips.

	ONPG	ADC	LDC	ODC	CIT	H <sub>2</sub> S	URE	TDA	IND	V P	GEL	GLU	MAN	INO	SOR	RHA	SAC	MEL	AMY	ARA	
Y ent	w+/0	0	0	0	0	0	+	0	0	0	Y ent	0	+	+	0	+	0	+	0	0	+
S son	w+/0	0	0	+	0	0	0	0	0	0	S son	0	+	+	0	0	+	0	0	0	+

The profile code number for *Yersinia enterocolitica* (Y ent) is 1014522.

The profile code number for *Shigella sonnei* (S son) is 1104112.

The differences between the two species are the reactions with ornithine, urea, sorbitol, rhamnose, and sucrose.

w = weak positive.



**8-35** Characteristic reactions of *Proteus mirabilis* (top), *Morganella morganii* (middle), and *Proteus vulgaris* (bottom) using API 20 E strips.

	ONPG	ADC	LDC	ODC	CIT	H <sub>2</sub> S	URE	TDA	IND	V P	GEL	GLU	MAN	INO	SOR	RHA	SAC	MEL	AMY	ARA
P mir	0	0	0	+	+	+	+	+	0	0	P mir	0	+	0	0	0	0	0	0	0
M mor	0	0	0	+	0	0	+	+	+	0	M mor	0	+	0	0	0	0	0	0	0
P vul	0	0	0	0	+	+	+	+	+	0	P vul	0	+	0	0	0	+	0	+	0

The profile code number for *P. mirabilis* (P mir) is 0734000.

The profile code number for *M. morganii* (M mor) is 0174000.

The profile code number for *P. vulgaris* (P vul) is 0674021.

The differences between *M. morganii* and *P. mirabilis* are the reactions with citrate, H<sub>2</sub>S, and indole.

The differences between *M. morganii* and *P. vulgaris* are the reactions with ornithine, citrate, H<sub>2</sub>S, sucrose, and amygdalin.




**10-38** **Correlating colonial morphology with biochemical identification.** An important step in the final identification of an isolate is to compare the biochemical reactions with the colonial morphology and position on the primary plating media. In this example, the colonies of *Salmonella* spp. are 1)  $\beta$ -positive on the 1% agar and the 1% D10 also positive on the 1% D10 E.


## CHAPTER 9 Other Gram-Negative Microorganisms

**T**he gram-negative microorganisms presented in this section are those that are not in the family *Enterobacteriaceae*. Included are the *Nisseria* and *Moraxella* species, miscellaneous fastidious gram-negative bacilli (*Moraxella*, *Bordetella*, *Breuerella*, *Pasteurella*, *Pasteurella*, *Kingella*, and *Capnocytophaga* spp.), the curved or comma-shaped gram-negative bacilli (*Klebsiella* and *Vibrio* spp.), and the nonfermentative gram-negative bacilli (*Pseudomonas*, *Syngomonas*, *Stenotrophomonas*, *Serratia*, *Acinetobacter*, and *Flavobacterium* spp.).

### *Nisseria* and *Moraxella* spp.




**FIG. 1** Gram stain of *Neisseria* spp. (small, paired spheres) from a blood culture showing gram-negative diplococci characteristic of *Neisseria* spp. The adjacent sides of the cell pairs lack lateral attachment; they are tapered. The diplococci are gram-negative and, as a result, they stain gram-variable in gram-positive ( $\times 1,200$ ).



**FIG. 2** *Moraxella meningitidis* on chocolate agar. *Neisseria meningitidis* colonies on a chocolate agar plate were found at 37°C in 18, 113, and 49% to 60% humidity. The colonies are more gray than yellow and nonhemolytic; characteristics that distinguish this pathogenic species of *Neisseria* from most of the nonpathogenic species. In other cultures, 148 hours' age, white and/or gray, filamentous colonies in the area of heavy growth.



**FIGURE 9-12** *Moraxella meningitidis* on 8% sheep blood agar. (Also increased opacity, opaque glaucous, and 1 to 1.5 mm in diameter. The colonies are grayish to color and not repeat pinkish as shown in this figure.)



**FIGURE 9-14** *Moraxella meningitidis* on enriched Thayer-Martin (MTM) agar. MTM is a selective agar medium for the isolation of *N. gonorrhoeae*. *Neisseria meningitidis* grows (nonoptimal), viable, opaque, and translucent. Use to inhibit gram-positive and other gram-negative microorganisms, mainly nonerythroid and for serotyping of *Neisseria* spp. After overnight incubation at 35°C in 5% CO<sub>2</sub>, colonies are small (0.5 to 1.0 mm in diameter), gray/glistening, and opaque. They increase in size with prolonged incubation. Culture should be held for 72 hours before reporting as negative.



**FIGURE 9-13** *Moraxella (Branhamella) catarrhalis* on 8% sheep blood agar. *Moraxella catarrhalis* grows well on routine laboratory enriched media including blood and chocolate agar media. Colonies are evenly opaque, round, and white. Microscopically they may be confused with nonpathogenic strains of *Neisseria* spp. but they are so differentiated by catalyphase utilization tests. They reduce NO<sub>2</sub> and unlike *Neisseria* spp. they produce DNase.



**FIGURE 9-15** Identification of *Moraxella* spp. by catalyphase utilization. The standard method for identifying *Neisseria* spp. is to determine oxidase and catalyphase (optimal) agar medium. A nonoxidative agar with 2% of each of the catalyphase glucose, lactose, sucrose, and sucrose. In this example, oxidase is produced in both glucose and maltose identifying the strain as *N. meningitidis*. This is an unreliable method because *Neisseria* spp. are oxidative, and the medium was designed to detect fermentative microorganisms.





**9-7 Identification of *Haemophilus* spp. by API 20E-NEPH-2.** API 20E-NEPH-2 (bioMérieux) is a 20-well microtiter for determining catalytic tests: catalase, oxidase, nitrate reductase, urease, and lysine decarboxylase. It is included in carbohydrate control (CTC), glucose (G), indole (I), malonate (M), lysine (L), urea (U), and acetone (A) 250 µl, and beta lactamase (B-LAC). Placed red in the wells are the reaction pattern. Glucose fermentation is characteristic of *H. meningocoli*, *H. parvulus* without sedimentation in glucose only.



**9-8 Identification of *Haemophilus* spp. by Gram-negative catalase-negative microorganism (GNM) panel.** The 24-well microtiter strip (bioMérieux) is a panel microtiter for conventional single-substrate biochemical tests to identify *Haemophilus meningocoli* spp. and the GNCE (Gram-negative catalase-negative) and *H. parvulus*. *Haemophilus* (family Pasteurellales) *Haemophilus* spp. group of organisms. Tests include: acetate-p-nitroprusside (A), phosphorylase (PH), protein-p-aminosalicylic (PA), gamma-glutamyl-p-nitroprusside (GPN), glycerol-p-nitroprusside (GLN), lysine-p-nitroprusside (LN), p-nitroprusside phosphorylase (PN), glucose (GL), acetone (AC), malonate (MT), hydrolytic nitrocellulose (HTC), urease (UR), ureastatic (US), urea (U), and urea-urea (UU) tests (GNM) catalase (CAT) and nitroprusside (NP) tests. The organism in this photograph is *H. meningocoli*.

### ***Haemophilus* spp.**




**9-9 *Haemophilus influenzae* on chocolate agar.** *Haemophilus influenzae* colonies on chocolate agar are smooth round, white, flat, umbonate (upside) and purple in color. *Haemophilus influenzae* (bioMérieux) and *H. influenzae* grown both before and after on chocolate agar.



**9-10 *Haemophilus influenzae* on chocolate agar (pinax up-view).** Colonies of smooth, encapsulated organism *Haemophilus influenzae* (type) were cultured on chocolate agar on chocolate agar or on chocolate agar. In this example, *Haemophilus influenzae* (type) which has been described as *Haemophilus influenzae* (type).



**9-11 X, Identification of *Haemophilus influenzae* by X and V strips.** The growth requires both X to use the iron (hemin) and V factor (nicotinic acid) (see also the [NAE] abbreviation) or nicotinic acid or nicotinamide phosphate (NADP) also known as *co* factor. The two important features for identifying this pathogen are: 1) the method like paper strips or disks are impregnated with X factor, V factor, and both factors V response (1) requires a spread over a surface that is devoid of the factors, and the strips are then applied. The plate is incubated overnight at 37°C, and observed for the presence of growth surrounding the strips. This is why growth only in the presence of both X and V factors confirmed by the fact that growth surrounding the X strip. **9-12 *Haemophilus parainfluenzae* multiplying around *Staphylococcus aureus*.** The blood in the agar provides the factor (X factor) and the ability to synthesize the NAD (V factor) allowing growth of the *Haemophilus influenzae* near the streak.



**9-13 Identifier of *Haemophilus influenzae* by the Porphyria Production Test (PPT).** This test differentiates *Haemophilus parainfluenzae* species from *Haemophilus influenzae* (type) and *H. parainfluenzae* from an enzyme porphobilinogen synthase that converts hematin to uroporphyrin III and porphobilinogen. *Haemophilus parainfluenzae* and *Haemophilus influenzae* detect porphobilinogen. Porphobilinogen exhibits fluorescence under a Wood's lamp (UV light at 365nm). With Wood's lamp: *H. parainfluenzae* appears red and *H. influenzae* fluoresces colorless. Under a Wood's lamp, a fresh culture of *H. parainfluenzae* fluoresces and *H. influenzae* does not on their base.

**9-14 *Haemophilus parainfluenzae* on a chocolate agar.** Colonies of *H. parainfluenzae* are approximately 1 mm in diameter, round, somewhat opaque and glistening. This species is commonly found in the upper respiratory tract as part of the normal flora. It is important to distinguish a beta *Haemophilus influenzae* a probable pathogen in the lower respiratory tract. Both species can cause significant pleuropneumonia and meningitis. Differentiation of the species is performed with the X and V strips (Figure 9-11) or the porphyria production test (Figure 9-13).




**8-14 Identification of *B. abortus* and *Brucella melitensis* with metabolic systems.** The metabolic biochemical differentiation system (BD) Microtiter assay system (Calspan Co.) Microtiter assay paper disks impregnated with substrates to test for the differentiation of proteolytic and gas-negative microorganisms incubated 48 hr. There are approximately 90 different substrates that may be used optionally with systems 3, 4, 5 and 6. Each plate is used for testing. In this example the following metabolisms are tested: glucose (G), melitin (M), urease (U), lysine (L), ornithin decarboxylase and cytochrome (C).


Organism	G	M	U	L	MD	C
<i>B. abortus</i>	+	-	-	+	+	+
<i>B. melitensis</i>	+	-	+	+	+	+

The distinguishing feature is that *Brucella abortus* does not urease lysine.

## ***Bordetella*, *Brucella*, and *Pasteurella* spp.**



**8-16 Gram stain of *Brucella* with G-1000G.** 70% long primary axis bacilli stained into a typical picture of a pairset with fourfold. Both *Bordetella* and *Brucella* spp. may have similar morphology.



**8-18 *Bordetella bronchiseptica* on 1% sheep blood agar.** Colonies on base (1) can be observed. Be, and that they grow well on 1% sheep blood agar incubated overnight. Their appearance and color on blood agar is similar to many of the aerobic catalase gas-negative bacilli but they do not grow on MacConkey agar. Rapid urea hydrolysis distinguishes this species from other microorganisms.



**9-17** Identification of *Brucella abortus* by urea hydrolysis test. Urea hydrolysis test rapidly. A pink color change re-develops within 15 minutes after inoculation. The right test turned strongly pink during overnight incubation.




**9-18** Colonies of *Brucella* sp. on chocolate agar. Growth of this along on chocolate agar after 4 days of incubation. Colonies require prolonged incubation to reach sufficient growth for further biochemical testing.

**9-20** Identification of *Brucella* spp. by nitrate reductase test. Some *Brucella* species hydrolyze urea in less than 10 minutes and/or nitrate to nitrite or nitroprusside. This holds a certain positive and reduced NO<sub>2</sub> to NO<sub>3</sub>. After 48 hours the addition of a trimethylamine sulfide and indole alpha-naphthylamine are added. The presence of nitrite is indicated by a red color change. If nitrite was not reduced or was re-oxidized nitrite storage gas, the bright orange coloration. A negative appearing as a red-brown is confirmed by adding more fluid, which reduces the oxidized nitrite to nitrite and covers the red color to form red, as shown here.



**9-21** Colonies of *Brucella* sp. on 6% sheep blood agar. These colonies of this after 48 hours old. There are no complete aerobic growing, but colonies may still grow on blood agar media without added extra blood. Growth may take 48 hours to become visible and even longer for a subculture from a single colony to grow. Colonies are small point, white to cream in color, non-glancing. The growth on media is incubated 15°C to 25°C O<sub>2</sub> for 7 days for best decaying as required.





**9-21 *Pasteurella multocida* on 5% sheep blood agar.** *P. multocida* is a somewhat fastidious organism that grows slowly on 5% sheep blood agar as a round grayish, nonhemolytic colony. It is oxidase- and catalase-positive and produces indole and ornithine decarboxylase.

### ***Eikenella*, *Kingella*, and *Capnocytophaga* spp.**




**9-22 *Eikenella corrodens* on 5% sheep blood agar.** Microscopically, *Eikenella corrodens* are tiny gram-negative bacilli. They are fastidious microorganisms and their growth is enhanced by CO<sub>2</sub> and the X factor. The microorganism grows on 5% sheep blood agar and pits or corrodes the agar surface; colony edges tend to spread. The colonies are distinguished by their bleachlike odor. The species is oxidase-positive and catalase-negative. This microorganism does not grow on MacConkey agar and is indole- and urea-negative. It does reduce nitrate to nitrite and produces lysine and ornithine decarboxylases.




**9-23 *Kingella kingae* on 5% sheep blood agar.** Microscopically, *Kingella* spp. are plump, gram-negative coccobacilli and can appear as gram-positive because of their tendency to retain crystal violet. Colonies are oxidase-positive and catalase-negative. Colonies of *Kingella kingae* are beta-hemolytic, and the species is indole- and nitrate-negative. They are fastidious microorganisms; growth on MacConkey agar is variable.



**9-24** Colonies of *Kingella kingae* on 5% sheep blood agar showing  $\beta$ -hemolysis. This figure shows the very distinct beta hemolysis displayed when the agar plate is held up to a light source.




**9-25** Colonies of *Capnocytophaga* on 5% sheep blood agar. They are fastidious, facultative anaerobes and may appear as tiny yellow colonies after overnight growth on blood agar. Growth is enhanced by  $\text{CO}_2$ . Colonies are slow-growing and often require 48 hours to detect distinct colonies. Microscopically, *Capnocytophaga* spp. are thin, gram-negative bacilli with pointed ends. These microorganisms were previously classified as *Bacteroides* spp. and as CDC group DF-1.



**9-26** *Capnocytophaga* on chocolate agar. Growth on chocolate agar shows the effect of gliding motility on the colonies produced by these microorganisms. A film surrounding the colonies is characteristic for growth of this microorganism. These microorganisms do not grow on MacConkey agar and are catalase-, oxidase-, indole-, and urea-negative.

## ***Campylobacter* and *Vibrio* spp. (curved gram-negative bacilli)**




**9-27 Gram stain of *Campylobacter* spp. ( $\times 1250$ ).** Microscopically, *Campylobacter* spp. appear as curved, comma-shaped, S-shaped, and gull-winged gram-negative bacilli. These shapes result from two cells remaining attached after division.




**9-28 Identification of *Campylobacter jejuni* by nalidixic acid susceptibility.** This microorganism grows best at 42°C, but will not grow at 25°C. It is microaerophilic, requiring 3% to 15% O<sub>2</sub> and 3% to 5% CO<sub>2</sub>. Selective media (i.e., CAMPY blood agar) is required to isolate this organism from mixed flora. It is oxidase- and nitrate-positive and urea-negative. A distinguishing characteristic is its susceptibility to nalidixic acid (*top disk*) and resistance to cephalothin (*bottom disk*).



**9-29 *Campylobacter jejuni* on TSI slant with lead acetate strip.** H<sub>2</sub>S is detected with lead acetate paper, but not in the TSI medium. *Campylobacter jejuni* also hydrolyzes hippurate (Figure 6-14 on p. 40).



**19-19** **A**, *Vibrio cholerae* on 0% sheep blood agar. Colonies of *Vibrio cholerae* grow well on sheep blood (strong tryptic media). Growth is enhanced by addition of 1% NaCl. *Vibrio cholerae* (strain 76 NaCl) is sodium-positive, DNase 0-1%, indole (+), urea- and catalase-positive. **B**, *Vibrio cholerae* on Trypticase-yeast-bleb water-sucrose (TYBS) agar. Colonies of *Vibrio cholerae* appear yellow due to excessive fermentation of TBS. Other agar (agarose NaCl) are in label.



**20-20** **A**, *Vibrio parahaemolyticus* on 0% sheep blood agar. *Vibrio parahaemolyticus* does require NaCl for growth and can tolerate 10% NaCl. Colonies on 0% sheep blood agar are approximately 1 to 5 mm in diameter, dull and can resemble a pale blue-green color. Colonies grow well on 0% blood agar containing catalase, indole- and urea/urea, decarboxylase positive. This species ferments glucose, malic acid, and sucrose. This organism has been associated with wound infections and gastroenteritis. **B**, *Vibrio parahaemolyticus* grows colorless and *V. cholerae* (profeta scintillans) on TYBS. *V. parahaemolyticus* can ferment almost all sugars (but not lactose and glycerol).



## Nonfermentative Gram-Negative Bacilli



**FIGURE 8-22** Oxidative-fermentative (OF) medium. In a tube, many Gram-negative bacilli from oxidative and fermentative pathways are detected in OF medium (right) and often a color transition like pH values in ThiO. The protein-cysteine/serine ratio in OF media prevents the acidification of weak acid by cellular end products. The color transition described in this section on oxidative and fermentative gram-negative bacilli. The tube on the left contains the basal medium alone; the other two tubes contain bacteria. A positive, but inconsistent, is indicated by a yellow color under tube with a negative test (oxidative media substrate) shows an color change (green right tube) fermentative substrate produces acid in both the oxidase (left) tube and the fermentative tube (colored with different of Figure 8-4 on p. 17)

## Pseudomonadaceae



**FIGURE 8-23** **A**, *Pseudomonas aeruginosa* on 0% sheep blood agar. The culture shows a typical growth on 0% and agar on typical (yellow-green and later brownish). However, the pattern—no growth of its own species. But, a common diagnostic pigment is produced by *P. aeruginosa* strains. Most colonies have a distinct grape-like color (purple) and fluoresce. *P. aeruginosa* may be distinguished from other species by its ability to grow at 42°C. **B**, *Pseudomonas aeruginosa* on 0% sheep blood agar. Oblique lighting demonstrates multiple colors and the surface of colonies.



**9-24** *Pseudomonas aeruginosa* on MacConkey agar. (Lab. MacConkey, cdk, see note on about cultured agar) growing on MacConkey agar. Colonies grow well at 37°C aerobically. They will also grow at 25°C.



**9-25** Pigment production by *Pseudomonas aeruginosa*. The strains of *Pseudomonas aeruginosa* growing on the Congo red agar. The red pigment is due to the production of pyoverdins, a water-soluble pigment, and the yellow-green color is due to pyoverdins, a fluorescent water-soluble pigment.




**9-26** *Pseudomonas aeruginosa* on TSI slant. On a TSI agar slant *Pseudomonas aeruginosa* appears as a MacConkey, nonfermenter, *Serratia* on a layer of pyoverdins.




**9-27** Identification of *Pseudomonas aeruginosa* by characterization in G20 tubes. The G20 system (Difco Laboratories, Lincoln, Neb.) is useful for identifying colonies of gram-negative bacilli that do not belong to the *Enterobacteriaceae*. The system consists of the G20 glucose- $\beta$ -D-glucosaminidase G20 tube (Figure 9-14) and the Ind G20 TSI plate (Figure 5-54). The correspondence between a fermenting organism (under a 5-hour temp.) to red color under the conditions indicate  $H_2$  gas production. *P. aeruginosa* does not ferment glucose.



**9-38 Identification of *Pseudomonas aeruginosa* by GNF, uree, and Pseudosel (42P) agar.** In this figure the GNF slant described in Figure 9-37 is on the left, a positive urea slant (Figure 8-16 on p. 61) is in the center, and 42P, the second tube of the N/F screen, is on the right. The 42P medium, also known as pseudosel agar, is incubated at 42°C for 18 to 24 hours. This microorganism grew at 42°C and produced the blue-green pyocyanin pigment. Based on the reactions in this figure, the microorganism can be definitively identified as *Pseudomonas aeruginosa*.




**9-40 *Pseudomonas stutzeri* on 5% sheep blood agar.** Characteristic colonies of *Pseudomonas stutzeri* on 5% sheep blood agar are buff-colored, dry, wrinkled, and adhere to the agar surface, making it difficult to remove them from the plate. Most strains grow at 42°C and reduce nitrates to nitrogen gas.




**9-39 Identification of *Pseudomonas aeruginosa* by Uni-N/F-Tek Plate.** The Uni-N/F-Tek plate (Remel Laboratories, Lenexa, Kan.) is a 13 test unit consisting of H<sub>2</sub>S, indole, a carbohydrate control, glucose, xylose, mannitol, lactose, maltose, acetamide, escullin, urea, DNase, and ONPG (Figure 8-12 on p. 59). Reactions for *Pseudomonas aeruginosa* are:


GLU	XYL	MAN	LAC	MAL	ACET
+	+	+	+	0	0
ESC	Urea	DNase	ONPG	H <sub>2</sub> S	IND
0	+	+	+	0	0



**9-41 *Sphingomonas* (formerly *Pseudomonas*) *paucimobils* on 5% sheep blood agar.** Colonies of *Sphingomonas paucimobils* on 5% sheep blood agar are small and yellow-pigmented after 24 hours of incubation. This species was previously known as *Pseudomonas paucimobils* and CDC group 11k-1. It is a slow-growing microorganism at 35°C, but optimal growth occurs at 30°C. This microorganism does not grow at 42°C or on MacConkey agar. It is oxidase-, escullin-, ONPG-, and DNase-positive and nitrate- and indole-negative.




**9-42 *Stenotrophomonas (Xanthomonas) maltophilia* on 5% sheep blood agar.** Colonies of *Stenotrophomonas maltophilia* on 5% sheep blood agar are chartreuse to lavender-green and have a characteristic strong ammonia odor. Unlike other pseudomonads, this species is oxidase-negative. It is resistant to many antimicrobials, although trimethoprim-sulfamethoxazole is usually effective.



**9-43 Identification of *Stenotrophomonas (Xanthomonas) maltophilia* by the Uni-N/F-Tek.**

Reactions for *Stenotrophomonas maltophilia* are:

GLU	XYL	MAN	LAC	MAL	ACET
weak	0	0	0	+	0
ESC	Urea	DNase	ONPG	H <sub>2</sub> S	IND
+	+	+	+	0	0



**9-44 *Shewanella (Pseudomonas) putrefaciens* on 5% sheep blood agar.** *Shewanella* is the new genus name for the microorganism previously classified as *Pseudomonas putrefaciens* and CDC group 1b. Colonies of *Shewanella putrefaciens* on 5% sheep blood agar are slightly viscous and mucoid and pinkish to red-brown or orange-tan in color.




**9-45 TSI reaction of *Shewanella (Pseudomonas) putrefaciens*.** The key characteristic that differentiates *Pseudomonas putrefaciens* from frequently encountered nonfermenters and other related microorganisms is its ability to produce large amounts of H<sub>2</sub>S in TSI or KIA slants.




**9-46 DNase agar.** Production of the enzyme DNase, which hydrolyzes DNA, may be used to differentiate nonfermenting gram-negative bacteria as well as *Staphylococcus aureus* and *Serratia marcescens*. This DNase test medium contains toluidine blue complexed with DNA. Hydrolysis of DNA by the inoculated microorganism causes changes of structure of the dye to yield a pink color. Toluidine blue may inhibit growth of some microorganisms, so equivocal results should be retested with another method.

## Other Nonfermentative Bacilli




**9-47 *Acinetobacter* (*Acinetobacter calcoaceticus* var. *anitratus*) *baumannii* on MacConkey agar.** *Acinetobacter baumannii* was formerly designated as *Acinetobacter calcoaceticus* var. *anitratus*, *Moraxella lwoffii*, *Herellea vaginicola*, and *Mima polymorpha*. They do not produce a pigment on blood agar, but appear faint pink on MacConkey agar, as shown in this figure. They are oxidase- and nitrate-negative. Microscopically, these microorganisms appear as coccobacilli, predominantly in pairs, and for this reason they have been confused with *Neisseria* and *Moraxella* spp.




**9-48 *Flavobacterium meningosepticum* on 5% sheep blood agar.** Colonies of *Flavobacterium meningosepticum* on 5% sheep blood agar after 48 hours of incubation are approximately 3 mm in diameter. Growth on MacConkey agar is variable, and most strains are oxidase-positive. This species is indole-, esculin-, ONPG-, and DNase-positive, urea-negative, and nonmotile.

**9-49 *Burkholderia (Pseudomonas) cepacia* on MacConkey agar.** This species, formerly known as *Pseudomonas cepacia*, grows slowly, especially when recovered from cystic fibrosis patients, in whom it is a significant pathogen. Colonies on MacConkey agar are often bright pink or red after prolonged incubation (as shown here) due to lactose oxidation. The species is oxidase-positive, although many strains display a weak oxidase reactivity. Most strains are lysine decarboxylase-positive and oxidize a number of sugars, including lactose.



**9-50 *Legionella pneumophila* on buffered charcoal-yeast extract (BCYE) agar.** BCYE agar is a selective medium for the recovery of *Legionella* spp. This buffered medium is the agar of choice for the isolation of *Legionella* spp. because it contains the requirements for optimal growth of the microorganism: L-cysteine, iron salts, and a pH of 6.9. Antibiotics are added to inhibit the growth of other bacteria. Growth appears in 2 to 3 days, and the colonies are circular, glistening, entire, and measure up to 4 mm, as shown here.




## CHAPTER 10 Anaerobic Bacteria

**A**naerobic bacteria are unable to multiply in the presence of atmospheric oxygen, although they have varying degrees of oxygen tolerance. The pathogenic strains are divided into two groups: exogenous pathogens, primarily the spore-forming clostridia that live in the soil and cause traumatic wounds or sporadic infections to lead to tissue disease; and the endogenous pathogens, normal inhabitants of human mucosa or lumen that cause disease when they enter a normally sterile site after traumatic breakdown of the normal mucosal integrity.

Specimens must be protected from air during collection, transport, and laboratory handling to maximize recovery of these fastidious microorganisms. Because of the time and resources necessary to perform good anaerobic microbiology, most laboratories limit anaerobic processing to specimens likely to harbor anaerobes such as abscess aspirates, deep tissue biopsies, and specimens collected through a protected catheter system.



**10-1** Colonies of *Bacillus fragilis* on anaerobic blood agar. *Bacteroides* is not as color as are clostridia, gram-negative bacilli grow less on this anaerobic medium. Agar must be supplemented with hemin and vitamin K<sub>1</sub> to support the growth of many species of anaerobes. These colonies appear to resemble yeast, not mold.




**10-2** Colonies of *Bacillus fragilis* on anaerobic bile esculin agar. This medium can also grow clostridia to inhibit. *Bacteroides* gram-negative bacilli like this isolate does not resemble other than those of the *B. fragilis* group, which grow as large, gray to black colonies due to their hydrolysis of esculin in the medium.




**10-3 Growth of clostridia on egg yolk agar.**

The production of the iron-sulfide element (paper-like precipitate) extending from the surface into the medium (pink color) and lysine interference (clear on surface of colony that lifts) is used to differentiate among Clostridium species and to help identify some Lactobacillus species.



**10-4 Nitrate stick test.** In test for reduction of nitrate to nitrite, a nitrate reductase enzyme paper-disk is placed over the base of an invertible culture of the microorganism. After 48 hours incubation, drops of the nitrate reagent solution and *α*-naphthylamine are added to the disk. A reduction indicates the presence of nitrate.




**10-5 Spot inside stick test.** To test for production of indole, phenylalanine paper-disk is placed on an area of growth of a microorganism on a medium containing tryptophan (such as casein blood agar base). After the disk has become moist to contact, a drop of indole-pyrazole-sulfide reagent is placed on the disk. A greenish-brown color indicates the presence of indole, a breakdown product of tryptophan. Indole is used to help differentiate among some *Enterobacteriaceae*, *Proteobacteriaceae*, *Actinobacteriaceae*, *Penicillium*, and *Aspergillus* spp.



**10-6 Medium polymerized sulfonate (EPS) stick test.** Polystyrene or one other of the sticky substance growth substrates that is inhibited by EPS.







**10-7 Special potency antimicrobial disks for presumptive identification of anaerobes.** The pattern of susceptibility to the three antibiotic disks: kanamycin (1000 µg), colistin (10 µg), and vancomycin (5 µg) can help differentiate among anaerobic genera using these criteria:

Microorganism type	Kanamycin	Vancomycin	Colistin
<i>Clostridium</i>	S (1, 2)	S	R
<i>Bacteroides fragilis</i> group	R (1)	R	R
<i>Bacteroides ureolyticus</i> group	S	R	S
<i>Fusobacterium</i> species	S	R	S
<i>Veillonella</i> species	S	R	S
<i>Porphyromonas</i> species	R	S	R
<i>Prevotella</i> species	R	R	V
<i>Peptostreptococcus anaerobius</i>	R (3)	S	R
Other gram-positive cocci	S	S	R

- (1) Some strains are kanamycin-resistant
- (2) V. variable; R, resistant; S, susceptible
- (3) Rare strains are susceptible




**10-8 A, *Bacteroides fragilis* on Brucella blood agar.** The large gray, mucoid colonies and resistance to all three potency disks described in Figure 10-7 are typical of members of the *B. fragilis* group. The extra filter paper disk is for performance of the nitrate test. **B, Gram stain of *Bacteroides fragilis*.** Pale-staining, gram-negative bacilli with rounded ends. Both pleomorphism and irregularity in staining can be observed in this smear ( $\times 1250$ ).



A



B




C

**10-9 A, Gram stain of *Prevotella melaninogenica* ( $\times 1250$ ).** Small, pleomorphic, gram-negative coccobacilli and bacilli can be observed in this smear, a characteristic appearance of *Prevotella melaninogenica*. **B, *Prevotella melaninogenica* on *Brucella* blood agar.** The small, dark-appearing colonies, resistant to kanamycin (less than 10 mm zone of inhibition) and vancomycin, will fluoresce brick-red under ultraviolet light (Woods lamp). **C, *Prevotella melaninogenica* on *Brucella* blood agar.** After several more days of incubation, these colonies will take on a dark brown to black color, due to assimilation of heme in the medium.



**10-10 *Bacteroides fragilis* and *Prevotella melaninogenica* on *Bacteroides* bile esculin (BBE) agar.** The *Prevotella* species was streaked on the upper half of the plate, but growth was inhibited. The *Bacteroides fragilis* grows profusely, displaying a typical gray precipitate in the agar surrounding the colonies.



**10-11** *A. Pseudomonas aeruginosa* on *Serratia* blood agar. The dry white breakdown colonies form colony in center of the petri dish and the antibiotic disk pattern (Ampicillin) and colonies susceptible to ampicillin (antibiotic susceptibility of *P. aeruginosa*). *P. aeruginosa* hydrolyzes the casein that is present in the plate. On the disk, will turn blue because the species is ampicillin-resistant. **B. Dried stain of Pseudomonas aeruginosa (x1000). Pale yellow long, slender gram-negative, barrel-shaped tapered pointed ends drawn showing the characteristic appearance of *Pseudomonas aeruginosa*.**



**10-12** *Paenibacillus amarus* on *Serratia* blood agar. This gram-negative results in plain medium, weak gas, catalase, oxidase and hemolytic activities and susceptibility to vancomycin.

**10-13** *Paenibacillus amarus* on *Serratia* blood agar. This species forms very colonies but small diameter to large size colonies appearing pink-purple color. It is susceptible to kanamycin and vancomycin and resistant to clavulanic.




**10-14** *Bacillus pumilus* (biochemical) (chemo) resistance of *Paenibacillus amarus* (antigen). *B. pumilus* *P. amarus* (antigen) (chemo).



**10-15 Colonies of *Propionibacterium acnes* on *Brucella* blood agar.** These young colonies are small and white to gray-white; however, as they age, they can appear yellowish. They have been referred to as anaerobic diphtheroids because, when stained and examined microscopically, they resemble diphtheroids. *Propionibacterium acnes* are indol-positive, as seen on this plate, and catalase-positive. When an anaerobic diphtheroid is catalase- and indol-positive, it can be presumptively identified as *Propionibacterium acnes*.




**10-18 Reverse CAMP test for presumptive identification of *Clostridium perfringens*.** *C. perfringens* is streaked vertically, and *Streptococcus agalactiae* is streaked horizontally in this test. The hemolysis of the clostridium is synergistically enhanced by the hemolysis of the streptococcus in an arrowhead-shaped pattern.



**10-16 Gram stain of *Clostridium* species ( $\times 1250$ ).** Cells are parallel-sided, long, thin, gram-variable, and some show swollen ends indicative of spore formation.



**10-17 *Clostridium perfringens* on *Brucella* blood agar.** Colonies display the double-zone of beta hemolysis typical of this species. Colonies are large with peaked centers and irregular edges after 48 hours incubation.



**10-19 *Clostridium perfringens* on egg yolk agar.** Colonies and surrounding medium display the expanding white precipitate stimulated by lecithinase production, typical of *C. perfringens*.



**10-20** *Cinetobacter affinis* on Brucella blood agar. A distinct blue-black spiral pattern and positive indole test are a given positive diagnosis. This species has a distinct bitter taste, and color and smell become characteristic only of incubated liquid. Laboratory diagnosis of *Cinetobacter affinis* starts usually depends on first look of the taste or smell rather than on isolation of the microorganism.



**10-21** *Cinetobacter affinis* and *C. pasteurianus* on egg yolk agar. The *C. pasteurianus* (upper half of plate) produces abundant lactinase, while the *C. affinis* (lower half) grows but shows no reaction.



**10-22** *Cinetobacter ferrireducens* on Brucella blood agar. This species also shows the spiral pattern and positive indole test of a given positive organism. *C. ferrireducens* is one of several species of *Cinetobacter* that is able to grow anaerobically. Spores are formed more readily however, during anaerobic growth.



**10-23** *Cinetobacter ferrireducens* on Brucella blood agar. Growth on egg yolk is characterized by its ability, such that growth will cover the entire plate in a white film within a few days.



**10-24** Brown stain of *Cinetobacter pasteurianus* spores. This microorganism displays unusual, wavy spores and the green rhizoids staining appear of *Cinetobacter* species. Spores may not take up the stain, so they may appear as clear zones.


## CHAPTER 11 *Mycobacteria*

**T**he *Mycobacterium* species include the *Mycobacterium tuberculosis* complex and the nontuberculous mycobacteria. They are slender rod-shaped aerobic bacteria that require special media for growth. Acid-fastness is a key characteristic of these bacteria. Mycobacteria may be divided into slow and rapid growers. Slow growers require at least 7 days of incubation to produce visible colonies, while colonies of rapid growers usually appear in less than 7 days. The rapid growers may be confused with other related genera, *Corynebacterium*, *Nocardia*, and *Rhodococcus* because they may be partially acid-fast.

**11-1** *Mycobacterium smegmatis* strain JH1000. Mycobacteria are easily stained with carbol fuchsin, which binds the mycolic acids in their lipid-rich cell walls. This stain cannot be removed (decolorized) with acid alcohol and, therefore, the microorganisms are referred to as acid fast. There are two common acid-fast stains: Ziehl-Neelsen (ZN) and Kinyoun's. The difference is that the ZN stain requires heat during the staining process because the phenol component goes into effect in the Kinyoun's method. The decolorizer, acid alcohol, and the counterstain, methylene blue, are the same, for both methods. When stained, acid-fast bacteria stain red and the background takes on a blue hue.



**11-2** *Mycobacterium fluororescens* strain JH1000. In this method, a fluorescent dye is used to visualize mycolic acid in the cell walls of mycobacteria. Both auramine and rhodamine are dyes in fluorescein that can be used to stain mycobacteria. When stained with a fluorescent mycobacteria are bright yellow (auramine) or red (rhodamine) or magenta (rhodamine) with a black background due to the potassium permanganate counterstain. The advantage of the fluorescent stain technique is that it allows for a more rapid diagnosis of acid fast bacteria stained smears. This lower magnification increases the viewing field and decreases the reading time.






**11-3** *M. tuberculosis* on Middlebrook 751 J agar. Colonies appear white, uncolored dry and wrinkled. Middlebrook 751 J agar is widely used medium for the isolation and enumeration of mycobacteria, and is of special value. Although not the gold standard, 751 J agar can detect acid fastness that supports the accuracy of fast and routine assays of *M. tuberculosis* complex strains. The area for detection of nonmycobacterial isolates. One here is a rough colony of *M. tuberculosis* growing on Middlebrook 751 J agar.



**11-4** *M. tuberculosis* on Löwenstein-Jensen (LJ) agar slant. LJ agar for both good culture standard for isolating mycobacteria. It contains a regulated whole egg, glycerol, potassium and zinc. *M. tuberculosis* growth is to be like the growth of a stationary form. The characteristic of LJ medium is that it becomes hydrolyzed when culture starts to grow. It is a red color, most for detection. Limited here are rough, half-colored colonies that appeared within 3 weeks typical of *M. tuberculosis*.



**11-5** Colonies of *M. tuberculosis* (TBE). Very young colonies less than 24 days old, growing on Middlebrook 751 J agar. Typical characteristics of the very beginning of the growth characteristic of *M. tuberculosis* can be observed.



**11-6** Colonies of *M. tuberculosis* (TBE). Older colonies (1 to 4 weeks old) growing on Middlebrook 751 J agar. The colonies have a rough appearance and color, compared to the 6 days ago.



**11-27** *M. avium* on Middlebrook 7H9 agar. The densest culture layer of *M. avium* (10<sup>8</sup> CFU/ml) on the Middlebrook 7H9 agar plate is smaller overall colony than a slightly younger and a larger overall colony. These cultures containing both colony types are characteristic of *M. avium*. In most cases, *M. avium* strains can be distinguished from other *M. avium*-complex members (such as *M. abscessus*) by their ability to grow at 45°C.



**11-28** *M. avium* on Lefsonstein-Jensen agar slant. *M. avium* does not grow on Lefson-Jensen agar if it appears as a thin layer on a media of incubation. Colonies are buff colored or cream to red, but much smaller than *M. tuberculosis* (see page with Figures 11-44).




**11-29** Smooth colony of *M. avium* (ATCC). Microscopically, the smooth colony type on agar is a 10<sup>8</sup> CFU/ml with smooth or striated morphology (Figure 11-30).



**11-30** Rough colony of *M. avium* (ATCC). The face of the rough colony has a characteristic "fuzzy" appearance.






**18-10** *M. luteus* on Middlebrook 7H11 agar. *M. luteus* is a photochromogen group 1 of Runyon's classification—a chemical method for differentiating the mycobacteria. Growth appears on Middlebrook 7H11 agar in approximately 1 week. A characteristic feature of the photochromogens is their dependence on exposure to visible light for pigment production. Colonies here are the rough, wrinkled colonies of *M. luteus* before exposure to light; the colonies are dull-colored.




**18-11** *M. luteus* on Middlebrook 7H11 agar, following exposure to light for a few hours. Colonies of *M. luteus* have a strong yellow color. Another characteristic of *M. luteus* includes a positive catalase reaction, utilization of oleic acid, and acid hydrolysis of Tween 80—a fatty acid.



**18-12** *M. luteus* on *Lathrobium-like* substrate agar. The slant on the left was exposed to light, resulting in the strong yellow-colored colonies, while the slant on the right was not exposed to light, and the colonies remained dull-colored.



**18-13** Culture of *M. luteus* (CFU). Microscopically, this brown-red colony, growing on Middlebrook 7H11 agar, shows a dark, dense center and a wrinkled protuberance—a key attribute of *M. luteus*.




**11-15** *M. chelonae* on Middlebrook 7H11 agar. This rapid grower can appear on culture media within 2 to 4 days. *M. chelonae* belongs to the *M. fortuitum-chelonae* complex. *M. chelonae* can be distinguished from *M. fortuitum* because it does not reduce nitrates or assimilate iron, but it is susceptible to polymyxin B and resistant to ciprofloxacin.



**11-16** *M. chelonae* on chocolate agar. *M. chelonae* can grow on chocolate agar and appear as smooth, opaque colonies resembling staphylococci or yeast. Because these microorganisms can be associated with skin infections and can be confused morphologically with other isolates also associated with skin infections, definitive identification of such colonies is needed. Although they are not readily stained by the Gram stain, they may appear as weakly staining, beaded, gram-positive bacilli, suggesting the possibility of a rapidly growing mycobacterium.



**11-17** Colony of *M. chelonae* ( $\times 20$ ). Examined microscopically, the colonies appear dark and dense with smooth edges and a somewhat lighter center, as shown here.



**11-18** *M. fortuitum* on Middlebrook 7H11 agar. The rapidly growing, rough colonies are shown after 3 days of incubation. *M. fortuitum* can also grow on modified MacConkey agar, without crystal violet, at 37°C, at 43°C on 7H11 and LJ, and in 5% NaCl at 37°C. These characteristics help to distinguish this species from *M. chelonae*.



**14-120** *M. fortuitum* on L-lysine-desoxycholate agar. The rough colonies are shown growing on L agar after 48-hour incubation.



**14-121** Detail of *M. fortuitum* (14-120). The rough colony of *M. fortuitum* observed microscopically with transmitted light demonstrates a somewhat smooth edge with a very dense core.



**14-122** *M. goodii* on 20% agar. Yellow-orange pigmented colonies of *M. goodii* are belong to the furrow group (*M. goodii* group). Colonies are laid by pigmented colonies on the absence of exposure to light. Colony given after 48 hours of 1 week.



**14-123** *M. goodii* on L-lysine-desoxycholate agar. The yellow-orange pigment is also produced on L agar.




**11-023** Colony of *B. anthracis* on methylene blue agar (100x). The entire colony is convex and it has characteristic non-specific antimicrobial lysis.



**11-024** *Bacillus anthracis* tolerance test. Some species of anaerobic bacteria, including the rapid growers found for *B. anthracis*, can grow at 20°C on an agar-based medium containing a 0.04% *M. anthracis* desoxyribose (D) while *B. anthracis* grows (right).



**11-025** Iron uptake tests. The iron uptake test is used as a way to determine the taxonomic and other rapid growers from *B. anthracis*. *B. anthracis* and other rapid growers have the ability to convert ferric iron/iron(III) chloride to an iron(II) form. The iron(II) is soluble in a reddish brown or rust color at the bottom. In this illustration, the tube on the left is the negative control (the entire tube was clear), while the culture is negative and the tube on the right, inoculated with *B. anthracis* in the positive control.



**11-026** Nitrate reductase test. Most strains of *B. anthracis*, *B. cereus*, and some strains of *Bacillus subtilis* convert the nitrocellulose barrier inside of *M* into minute amounts into the culture medium. Drops of paper strips impregnated with all substances are used to test the source, which reacts with a reagent to build in the presence of a primary oxidant to produce a yellow color. In this example, in which the paper strips were removed after color development, the tube on the left is positive (yellow) and the tube on the right is negative (colorless).

## CHAPTER 12 *Microbial Pathogens Isolated and/or Identified by Tissue Culture or Other Special Methods*




**FIGURE 12-1** *Cytomegalovirus*. Cytochrome cell culture assay using MRC-5 cells. Phase contrast microscopy ( $\times 400$ ). Culture of C cells at low density allows one to be quite formal very specific since the technique can be proven to accept passive neutralization. For this reason, the direct plaque assay for the detection of the *C. cytomegalovirus* is not recommended. Laboratory diagnostic method. Figure 12-1 shows the basic effect (usually by a representative cell using tests added to the MRC-5 monolayer) while Figure 12-2 shows a representative neutralizing antibody (the test) that has been added, preventing development of the cytopathic effect.




**FIGURE 12-2** *Cytomegalovirus*. Direct fluorescent antibody assay (DFA) of cervical cells ( $\times 400$ ). To perform the assay cells from the cervix (or other) are collected with a swab placed on a glass slide, and fixed with acetone. A fluorescent labeled specific antibody is then added and the specimen is examined using a fluorescence microscope. In this figure, 1 *Cytomegalovirus* elementary bodies have been stained with a monoclonal antibody to the outer envelope protein.



**12-4** *Chlamydia trachomatis*. Direct antigen detection. Three dot blots are arranged in a dot blot filter for antibody immunoscreening systems for the direct detection of *C. trachomatis*. For each of the two assays depicted here, the specimens is placed in the lower window of the card. The cell-specific antigen through the filter and binds to a *Chlamydia*-specific antibody located on a line in the filter at the end of the second window. A positive reaction results in the formation of a color line across and on the left side of each system. Another antibody located in the filter at a greater distance from the specimen window is used as a control to ensure that the specimen has migrated through the filter.




**12-5, a, b** *Chlamydia trachomatis*. Shell vial culture. Fluorescent and immunofluorescence micrographs. Shell vial with 960,000 to 1.6 million cells are frequently used for the isolation of *C. trachomatis* from clinical specimens. Following inoculation of the sample onto the cells, the cell monolayers are checked for 48 hours. If viable *Chlamydia* were present in the sample they will multiply in the cells of the monolayer and produce intracellular inclusions. The final inclusions can be stained with a fluorescent-labeled monoclonal antibody (Figure 12-5a) or with immunofluorescence (Figure 12-5b) without antibody labeled with immunofluorescence (Figure 12-5a). A similar approach can be used for the detection of *Chlamydia pneumoniae* and *Chlamydia psittaci*.




133

**13-7-9** *Chlamydia trachomatis* and *Chlamydia pneumoniae*, infected tissue culture cells viewed by transmission electron microscopy. Left, *Chlamydia trachomatis* inclusion (in focus) by TEM (magnification 21 to 73 times of in-focus growth). The elementary bodies are electron-dense inclusions of the organism and are dense compact spherical bodies approximately 500nm in diameter. The outer (cell) bodies are large, non dense, have dense structure and measure up to 1,000 nm in diameter. Clearly, normally the inclusion of *C. trachomatis* surround the inclusion of the host cell (Chlamydia from the cytosol's membrane) thus they their material (Figure 13-7, 010710). *C. pneumoniae* may form multiple inclusions (Figure 13-9, 010710). Some of the elementary bodies of *C. pneumoniae* may form a peritrichous shape (Figure 13-9, 010710).



134



135


12-10



12-11



12-12




#### 12-10-12 Mycoplasma and Ureaplasma species


**Mycoplasmas** are the smallest freely living organisms. Unlike other bacteria they lack cell walls and are so gelatinous on the surface of the host cells that they are intracellular. Complex media and tissue culture techniques are used for their isolation. *Ureaplasma parvum* was originally named *T-titer mycoplasma* due to the small size of their colonies (usually 15 to 30 µm in diameter) (Figure 12-12) viewed through a light microscope. *U. urealyticum* takes an alkaline color change (purple) without turbidity in a medium that contains urea (Figure 12-11, left side, control tube plus the right). The *pro-fal* mycoplasmas grow better in 1% CO<sub>2</sub> and produce the typical Friedberg's appearance morphology (Figure 12-12 about 14, dense colonies viewed through a light microscope).



**12-13 *Ureaplasma urealyticum*. Transmission electron microscopy ( $\times 77,000$ ).** The ultrastructure of the mycoplasmas and ureaplasmas is pleomorphic and fairly indistinct. They measure 200 to 300 nm in diameter and have a cell membrane, but they do not have a cell wall.



**12-14 Mycoplasmas and ureaplasmas. Antibody detection by microimmunofluorescence (MIF) test.** There are several approaches, including complement fixation, metabolic inhibition, and fluorescence assays, for detecting antibodies to these microorganisms. For the MIF test, the serum of the patient is reacted with colonies of mycoplasma present on the slide, and antihuman fluorescein-labeled globulin is then added. This test is fairly specific and sensitive and allows for the detection of IgM as well as IgG antibodies.



## CHAPTER 13 Antimicrobial Susceptibility Testing (AST)

**N**umerous antimicrobial susceptibility test procedures have been described. Qualitative tests, such as the agar disk diffusion method, categorize the bacterial isolate as susceptible, intermediate, and resistant. Quantitative methods are used to determine a specific endpoint, referred to as the minimal inhibitory concentration (MIC; the minimal amount of antimicrobial agent that inhibits visible growth).

In the disk diffusion method, used for rapidly growing aerobic bacteria, a standardized pure culture inoculum is swabbed onto the surface of a Mueller-Hinton agar plate. Filter paper disks impregnated with antimicrobial agents are placed on the agar plate. After overnight incubation, the zone of inhibited growth surrounding each disk is measured and results are compared to guidelines published by the National Committee for Clinical Laboratory Standards (NCCLS, Wayne, Pa.).



**13-1 Preparation of standardized inoculum.** *Staphylococcus aureus* colonies of similar morphology are selected on a routine culture from a broth, 1:5 in 10% saline culture and inoculated into 5 to 9 ml Mueller-Hinton broth or trypticase soy broth. After the broth is incubated for 2 to 4 hours or the selected colonies are traced back from 0.45% sterile saline or broth without its alcohol, a turbidity suspension of the desired turbidity. The direct method is preferred for those concentrations that are turbid and have unappreciable growth to broth.



**13-2 Using a turbidity meter to adjust the inoculum.** The suspension is vortex mixed for 25 to 30 seconds before determining its turbidity. A meter is used to adjust the inoculum to a 0.5 turbidimetric turbidity standard, which is equal to approximately  $1.5 \times 10^8$  colony-forming units (CFU) per ml. The turbidity is adjusted to match the standard by diluting with sterile saline or broth.



**13-3 Disk diffusion method—Agar plate inoculation.** A sterile cotton swab is dipped into the inoculum and the entire surface of the Mueller-Hinton agar plate is swabbed three times by rotating the plate approximately 60° between streaking to ensure even distribution. The plate is allowed to stand for at least 3 minutes, but no longer than 15 minutes, before applying the disks.



**13-4 Disk diffusion method—application of disks.** Disks are applied to the agar surface using a dispenser. No more than 12 disks should be placed on a 150 mm plate. A disk should not be relocated once it has made contact with the agar surface because antimicrobial diffusion begins immediately. Agar plates must be incubated within 15 minutes of disk application.



**13-5 Disk diffusion method—*Staphylococcus aureus* American Type Culture Collection (ATCC) 25923.** The zone of growth inhibition is measured using a millimeter ruler or calipers. Results are interpreted as susceptible (S), intermediate (I), or resistant (R). The size of the disk (6 mm) is included in the measurement; therefore, those agents with no zone are measured as 6 mm and are always interpreted as resistant. The microorganism tested on this Mueller-Hinton agar plate is susceptible to all antimicrobial agents.



**13-6 Disk diffusion method—*Escherichia coli* ATCC 25922.** The isolate tested on this Mueller-Hinton agar plate is interpreted as susceptible (S) to all antimicrobial agents. Reading clockwise from the top, they are mezlocillin (MZ), amikacin (AN), ampicillin (AM), cefazolin (CZ), cefotaxime (CTX), cefuroxime (CXM), cephalothin (CF), gentamicin (GM), and tobramycin (NN); the three disks in the center of the plate are trimethoprim/sulfamethoxazole (SXT), ceftolixan (FOX), and ticarcillin/clavulanic acid (TIM).




**10-7 Disk diffusion method—*Pseudomonas aeruginosa* antibiotic screen.** In the healthy lung, the isolate is sensitive to tetracycline, fusidic acid (FT) and susceptible to all others. Day 14 disk diffusion (DIF) results (FT) sensitive (S), tetracycline (TET) sensitive (S), fusidic acid (FT) sensitive (S), fusidic acid (FT) sensitive (S), rifampin (RF) sensitive (S), gentamicin (GM) sensitive (S), streptomycin (SM) sensitive (S), chloramphenicol (CM) sensitive (S), trimethoprim-sulfamethoxazole (TS) sensitive (S), and piperacillin (PI) sensitive (S).




**10-8 Disk diffusion method—*Pseudomonas aeruginosa*, a resistant strain.** 1 month in the healthy lung, the isolate shows that the isolate is resistant to doxycycline (DO), rifampin (RF) and susceptible to the remaining. The isolate is sensitive to tetracycline (TET), fusidic acid (FT), gentamicin (GM), streptomycin (SM), fusidic acid (FT), rifampin (RF) and erythromycin (EM). The isolate is susceptible to chloramphenicol (CM), fusidic acid (FT), streptomycin (SM), rifampin (RF), gentamicin (GM), erythromycin (EM), trimethoprim-sulfamethoxazole (TS) and piperacillin (PI).




**10-9 Disk diffusion method—*Streptococcus pneumoniae* secondary level battery when susceptible to resistant (less than 50%) to vancomycin.** A positive test on Mucic Acid agar with (M) drug disk. Clindamycin (CL) sensitive (S), rifampin (RF) sensitive (S), fusidic acid (FT) sensitive (S), gentamicin (GM) sensitive (S), streptomycin (SM) sensitive (S), chloramphenicol (CM) sensitive (S), trimethoprim-sulfamethoxazole (TS) sensitive (S), erythromycin (EM) sensitive (S), fusidic acid (FT) sensitive (S), rifampin (RF) sensitive (S), gentamicin (GM) sensitive (S), streptomycin (SM) sensitive (S), chloramphenicol (CM) sensitive (S), trimethoprim-sulfamethoxazole (TS) sensitive (S).



**10-10 E test method—*Streptococcus pneumoniae*.** The first MIC is 0.06 µg/ml (0.06 µg/ml) to determine the MIC of an organism. Microdilution minimum inhibitory concentration (MIC) is 0.06 µg/ml. The MIC of piperacillin (PI) is 0.06 µg/ml (0.06 µg/ml) and the MIC of gentamicin (GM) is 0.06 µg/ml (0.06 µg/ml).




**13-11 E test method—*Streptococcus pneumoniae*.** In this figure, the MIC of cefotaxime (CT) is between 1.5  $\mu\text{g}/\text{ml}$  and 2.0  $\mu\text{g}/\text{ml}$ . The MIC should be reported as 2.0  $\mu\text{g}/\text{ml}$ , which is interpreted as resistant.



**13-13 Broth Microdilution Method—Inoculated microdilution panel.** This Pasco microdilution panel (Difco Laboratories, Detroit, Mich.) includes wells for both identification biochemicals and quantitative antimicrobial susceptibility testing. If all antimicrobial dilutions show inhibition of growth, the endpoint is less than or equal to the lowest concentration tested, assuming that the positive growth control exhibits growth. Interpretation is based on NCCLS approved standards.



**13-12 Broth Microdilution Method—Inoculation of microdilution panel.** Microdilution panels are used to quantitatively measure the susceptibility of microorganisms to a battery of antimicrobial agents. The inoculum is mixed and poured into the inoculum/seed tray. The inoculator tips are lowered into the tray and the tips filled by capillary action. The inoculator tips are then lowered into the wells of the panel, as shown here.



**13-14 Broth Microdilution Method—Interpretation of microdilution panel.** Close-up view of 20 wells of the microdilution panel. Note the drug name and concentration in blue ink directly below the corresponding well. Row 1, left to right, the antimicrobials and interpretations are as follows: ampicillin/sulbactam (A/S) - no growth (MIC $\leq$ 2/1  $\mu\text{g}/\text{ml}$ ); clindamycin (CD) - no growth (MIC $\leq$ 2  $\mu\text{g}/\text{ml}$ ); erythromycin (E) - growth (MIC  $>$ 0.5  $\mu\text{g}/\text{ml}$ ); oxacillin (OX) - growth (MIC  $>$ 1  $\mu\text{g}/\text{ml}$ ); vancomycin (VA) - no growth (MIC  $\leq$ 2  $\mu\text{g}/\text{ml}$ ).

Row 2, gentamicin (GM) - interpretation: growth in all wells, MIC is  $>$ 8  $\mu\text{g}/\text{ml}$ .

Row 3, tobramycin (TO) - interpretation: growth in all wells, MIC is  $>$ 8  $\mu\text{g}/\text{ml}$ .

Row 4, wells 1 to 3, ceftriaxone (FRX) - dilutions, left to right, 32, 16, and 8  $\mu\text{g}/\text{ml}$ . Interpretation: growth in 8  $\mu\text{g}/\text{ml}$ , no growth in 16 and 32  $\mu\text{g}/\text{ml}$ ; MIC is 16  $\mu\text{g}/\text{ml}$ .

Row 4, wells 4 and 5, chloramphenicol (C) - dilutions, left to right, 16 and 8  $\mu\text{g}/\text{ml}$ . Interpretation: no growth in either



**18-10 Automated method—Preparation of standardized inoculum.** Several instruments for automated performance of MICs are available. Preparation of inoculum for the Vitek 2 (Becton Dickinson, Inc.) Automated MIC system is performed using the AutoMicros 360, which uses a single-board processor with a 1-MHz test filter. The inoculum suspension concentration and the turbidity measurement modulus are the same.



**18-11 Automated method—setting the Vitek card.** The standardized inoculum and the card containing the diluted antimicrobials are first placed into the filling module of the Vitek instrument. Once the inoculum is transferred into the readily reusable module through the plastic transfer tube, the card is sealed in the working module, as shown here.



**18-12 A, Automated method—Vitek card in the reader incubator.** The card has been sealed. Inward is required for a proper fit and then it is placed into the reader incubator as shown here. The results are determined by the electro-optical detection of turbidity growth in the wells of the test module. MICs are determined by visual analysis of bacterial growth rates in the presence of antimicrobial agents. **B, Automated Vitek card.** The Vitek Card capsule susceptibility card (Vitek 2) is used for performing MICs against any aerobic bacteria from water specimens. Once the growth-inhibited well reaches the turbidity threshold, the instrument reads the remaining wells, which contain various dilutions of concentrations in the tube. Turbidity (presence) are shown as by of the success, while only clear wells (susceptible) are shown as by the failure.

**10-09 Detection of beta-lactamase production by coliforms.** Beta-lactamase enzymes produced by bacteria break down penicillins and render the antibiotic unable to penetrate. The method chemically detects all known beta-lactamases. The colorless reagent (from Oxoid's Microbiology Systems, Catalogue No. M1) is impregnated with nitrocellulose on a chromatography strip. Dilutions of the microorganisms are added onto the disk. If the enzyme breaks the nitrocellulose beta-lactamase ring is hydrolyzed to a beta-lactamase which causes the color of the disk to rapidly change from white to red.



**10-10 Positive beta-lactamase test.** The pink color on the surface of the disk, that was inoculated with the test microorganism, indicates the presence of beta-lactamase. This test is used to determine beta-lactamase activity in *Staphylococcus*, *Haemophilus influenzae* (strain not clinically *Serratia marcescens*, *Neisseria meningitidis*, enterococci, and certain anaerobic bacteria including *Bacteroides* spp).



**10-08 Inhibition test/serum inhibitory test.** Determining the activity of an antimicrobial agent in the patient's serum may be useful in certain infections requiring prolonged antimicrobial therapy, such as endocarditis and osteomyelitis. For the inhibition test, the patient's serum containing an antimicrobial agent was diluted with the patient's own infusing plasma/serum. After inoculation, the test inhibitory (IT) is interpreted as the highest dilution of the serum showing no visible growth as shown on the far left tube, while the other three tubes illustrate growth of the microorganism.



## CHAPTER 14 Mycology

**I**dentification of fungi is still based largely on gross and microscopic morphological characteristics. In addition to visual inspection, commercially available biochemical, manual or automated systems, nucleic acid probes, and serological reagents are used. The fungi are broadly divided into yeasts, characterized by aqueous, creamy colonies, and molds that produce cottony or powdery growth. Some species, particularly those that cause systemic infections, are dimorphic and can assume either morphology.

This section presents the agents of fungal infections in a commonly used clinical classification. Systemic or deep-seated mycoses often begin in the respiratory tract and disseminate. Subcutaneous mycoses are frequently secondary to traumatic implantation. Systemic spread of these infections is less likely. Superficial infections are confined to the keratinous layer of skin, nails, and hair and in general there is minimal host inflammatory response. The "opportunistic mycoses" usually cause infections only in immunocompromised hosts. The sensitive fungi are common environmental inhabitants. The "true pathogens" category includes those fungal infections that do not fit easily into the other headings.

### Deep-seated Mycoses




**14-1** *Rhizoglyphus dermatitidis* colonies. *R. dermatitidis* growing on bread (left) releases spores with 10% desiccation potential and (micrographed), at 25°C the mold form of *R. dermatitidis* produces white to tan, cottony colonies that grow fairly rapidly in a week. With long-term incubation, at 12°C the mold form produces more or less unbranched, wavy looking colonies (not shown).



**14-2** *Blastobotrys dermatitidis* spores. *Blastobotrys dermatitidis* spores resemble *Aspergillus fumigatus*, at 25°C in the mold form. Spores of *R. dermatitidis* are macroconical, with a beaded, smooth wall, and are attached to short lateral or terminal hyphal branches giving the appearance of "claypegs." The spores measure 1.0-1.5  $\mu$ m in diameter.








**14-3, 4** Yeast form of *Blastomyces dermatitidis*. Lung tissue. Gomori methenamine silver stain (Figure 14-3,  $\times 1850$ ; Figure 14-4,  $\times 1250$ ). At  $37^{\circ}\text{C}$  the yeast form of *B. dermatitidis* appears as thick-walled, spherical cells with a wide-based single bud (Figure 14-3). Small forms may occasionally emerge, particularly in tissue, which cannot be easily differentiated from *H. capsulatum* (Figure 14-4).



**14-5** *Coccidioides immitis* growing on Sabouraud's dextrose agar. *C. immitis* grows well on Sabouraud's medium at  $25^{\circ}\text{C}$ . Within 1 or 2 weeks, the gray-tan colony has a powdery, cottony appearance resulting from the formation of arthroconidia from the hyphae. The colony becomes tan to brown with age, and the reverse side is white. Other colony colors are seen occasionally.




**14-6, 7** *Coccidioides immitis* wet mount (KOH; 14-6,  $\times 500$ ), and calcofluor white (14-7,  $\times 1200$ ) preparations. Mature spherules of *C. immitis* are round or oval, measure 15 to 60  $\mu\text{m}$  in diameter, have a well demarcated wall, and produce endospores by undergoing progressive cleavage. The endospores measure 2 to 5  $\mu\text{m}$  in diameter, do not bud, and are extruded from the mature spherules. Once they are extruded, the endospores initiate maturation, becoming spherules that eventually undergo cleavages that result in new endospores. It is important to recognize that immature spherules that lack endospores can be mistaken for other fungi, such as *Blastomyces dermatitidis*, and for artifacts.



**14-8 Lactophenol cotton blue preparation of arthroconidia of *Coccidioides immitis* ( $\times 500$ ).**

The arthroconidia of *C. immitis* appear in the branches of the hyphae as thick-walled, barrel-shaped structures measuring 4 to 6 by 2 to 4  $\mu\text{m}$ . Alternating with the arthroconidia are weakly stained empty cells, a characteristic that differentiates this microorganism from *Geotrichum* spp. The arthroconidia mature, break off, and following aspiration they can produce an infection in susceptible individuals.




**14-9 Tissue section of the lung stained with Gomori's methenamine silver showing a *Coccidioides immitis* spherule ( $\times 1250$ ).**

The wall of the spherule has ruptured and the endospores are released into the surrounding tissues. In general, the walls of the endospores stain darker than the walls of the spherules. Spherules that lack endospores need to be differentiated from other fungi including *Histoplasma capsulatum*, *Torulopsis (Candida) glabrata*, *Cryptococcus neoformans*, and *Paracoccidioides brasiliensis*.



**14-10 *Cryptococcus neoformans* colonies growing on chocolate agar.**


The colonies of *C. neoformans* are typically flat or slightly raised, shiny, mucoid as a result of the presence of a mucopolysaccharide capsule, and may have a wide variation in color ranging from cream to tan and pink. With age, the colonies tend to become drier and darker in color.




**14-11, 12 India ink (14-11) and Gram stain (14-12) preparations of *Cryptococcus neoformans* ( $\times 2250$ ).** The cells shown here are surrounded by a broad mucopolysaccharide capsule. The thickness of the capsule can vary widely from preparation to preparation, and in certain instances the capsule may be undetectable using India ink. The cells are round to oval, with a diameter ranging from 3 to 10  $\mu\text{m}$ , and have a single bud with a narrow neck. The India ink preparation is not as sensitive as the latex agglutination antigen detection tests for initial detection of *C. neoformans* in cerebrospinal fluid. Specimens negative by India ink should be tested for antigen; all cerebrospinal fluid specimens sent for India ink preparation should also be cultured.




**14-13 Urease test for *Cryptococcus neoformans*.** Over 99% of the *C. neoformans* isolates give a positive urease test, pink-purple color, within 15 minutes, in contrast to other urease-positive species of yeast that require more than 3 hours to give a positive reaction.




**14-14 Colonies of *Histoplasma capsulatum*.** Note the yeast form of *H. capsulatum* growing on chocolate agar at 37°C (right side), and the mycelial form growing on Sabouraud's dextrose agar at 25°C (left side). The yeast form appears waxy and moist with a yellow-tan color, while the mycelial form has a cottony appearance with a white-brown or pinkish color.




**14-15 Lactophenol cotton blue stained *Histoplasma capsulatum* ( $\times 500$ ).** After several days in culture at 25°C, thick-walled, tuberculate (knobby), and nontuberculate macroconidia appear, measuring 10 to 15  $\mu\text{m}$  in diameter. The microconidia of *H. capsulatum* emerge during the early stages of colony growth as spherical or oval structures, 3 to 5  $\mu\text{m}$  in diameter.




**14-16 Direct bone marrow preparation showing *Histoplasma capsulatum* stained with calcofluor white ( $\times 2500$ ).** Round to oval, small, 2  $\times$  5  $\mu\text{m}$ , narrow-based, budding yeast cells can be observed on this direct bone marrow preparation.




**14-17 Buffy coat preparation showing *Histoplasma capsulatum* ( $\times 1250$ ).** White blood cells containing blastoconidia of *H. capsulatum*. The budding yeast seems to be surrounded by a capsule. This pseudocapsule is thought to be an artifact resulting from shrinking tissues during fixation.




**14-18 Lung section showing *Histoplasma capsulatum* stained with Gomori's methenamine silver ( $\times 1250$ ).** The yeast form of *H. capsulatum* is usually found intracellularly in the cytoplasm of professional phagocytes. The microorganism is elongated, measures approximately 2 to 5  $\mu\text{m}$ , and multiplies by narrow-based, unequal budding.



**14-19 Mixed culture of *Histoplasma capsulatum* and *Cryptococcus neoformans*.** This patient had a double infection, and thus colonies of both *H. capsulatum* and *C. neoformans* are growing on this plate. This possibility should be taken into consideration, particularly in specimens from patients that are immunocompromised. Plates should be incubated for at least 4 weeks and should not be discarded earlier, even when a rapid grower is already present.




**14-20 *Paracoccidioides brasiliensis* stained with lactophenol cotton blue ( $\times 325$ ).** In this preparation, cultured at 37°C, large, thick-walled cells with multiple buddings are attached to the mother cells by narrow connections, giving the characteristic appearance of a "ship's wheel" or "steering wheel." At 25°C most strains grow for a long time without producing conidia. There are different types of conidia that may form, none of which is characteristic of this species.




**14-21 *Paracoccidioides brasiliensis*, Adrenal tissue. Gomori's methenamine silver stain ( $\times 500$ ).** This preparation demonstrates numerous thick-walled yeast cells, spherical to oval, measuring up to 60 to 70  $\mu\text{m}$  in diameter with multiple narrow-based buds that give the typical appearance of a "steering wheel" and "Mickey Mouse" type forms.

## Opportunistic Mycoses




**14-22 Lactophenol cotton blue preparation of *Alternaria* spp. ( $\times 500$ ).** The mycelia are septate and dark. The conidiophores produce tan-brown conidia that are large, approximately  $10 \times 30 \mu\text{m}$ , round at the end near the conidiophore, and narrow at the far end giving a clavate (clublike) shape, with a smooth or rough wall, and transverse and longitudinal septations.




**14-23 Diagram of the morphological structure of *Aspergillus* spp.** *Aspergillus fumigatus* (A) is uniseriate (e.g., has only one layer of phialides that cover the upper two-thirds of the vesicle). The phialides bear the conidia, which are extruded from the end of the unstricted phialide. *Aspergillus versicolor* (B) and *A. niger* (C) are biseriata, in which the vesicle is covered with a layer of short hyphal structures called metulae (the structures bearing the phialides) and another layer consisting of the phialides. The metulae and phialides of *A. niger* form a radiate arrangement (C) but those of *A. versicolor* do not (B). *Aspergillus flavus* may be both uniseriate and biseriata with the metulae and phialides covering the entire vesicle.




**14-24 Colonies of *Aspergillus flavus* growing on a sheep blood agar plate.** The colonies of *A. flavus* have a distinct velvety yellow to yellow-green or brown color. The green-brown color is more prominent in older cultures. The reverse is white to red-brown. This microorganism grows better at 37°C than at room temperature.




**14-25 Lactophenol cotton blue preparation of *Aspergillus flavus* (×225).** The conical heads measure approximately 300 to 400 μm in diameter and may have uniseriate and biseriate rows of phialides that cover the entire vesicle. Usually the proximal row of sterigmata is twice the length of the outer row. The conidia do not usually chain and tend to accumulate over the vesicle.




**14-26 *Aspergillus fumigatus* growing on a sheep blood agar plate.** The colonies of *A. fumigatus* grow well at 20 to 43°C. Incubating the specimens of *A. fumigatus* at temperatures above 43°C helps to inhibit the growth of contaminants. The color of the colonies ranges from white to green and with age they tend to turn gray, brown, or black.




**14-27 Gram stain of *Aspergillus fumigatus* (×1250).** *A. fumigatus* typically appears as uniform septate hyphae with dichotomous (two-way) acute (45°) branching. The hyphae have parallel walls without constrictions, and are fairly uniform in width, measuring 3 to 6 μm in diameter.



**14-28, 29** *Aspergillus fumigatus* stained with lactophenol cotton blue (14-28 [ $\times 225$ ]) or with calcofluor white (14-29 [ $\times 500$ ]). The conidiophore of *A. fumigatus* measures up to 300 to 500  $\mu\text{m}$  in length and 4 to 8  $\mu\text{m}$  in width. The vesicle is dome-shaped, 20 to 30  $\mu\text{m}$  in diameter, and merges with the conidiophore. Phialides are present only in the upper half or two thirds of the vesicle and extend in the same direction as the conidiophore. The small, oval-round, green-tan conidia are borne in chains from the tips of the phialides or sterigmata.



**14-30** Colony of *Aspergillus niger* growing on Sabouraud's dextrose agar. *A. niger* colonies are fairly distinct, with a woolly appearance, originally displaying a white-yellow color that turns to dark brown-black due to the conidial heads, with a white to tan basal layer.



**14-31** *Aspergillus niger* stained with lactophenol cotton blue ( $\times 500$ ). The conidiophores of *A. niger* are large, measuring up to 2 to 3 mm in length by 15 to 20  $\mu\text{m}$  in width. The vesicle is globose, 40 to 80  $\mu\text{m}$  in diameter, brown-black, with biserial phialides that cover the entire surface of the vesicle forming a radiate head.






**14-32 *Aspergillus terreus* colony on Sabouraud's dextrose agar.** This colony of *A. terreus* has a velvety white, cinnamon-buff to brown color on the surface, and is white to brown on the reverse side.




**14-33 *Aspergillus terreus* stained with lactophenol cotton blue ( $\times 500$ ).** The conidial heads of *A. terreus* measure up to 500  $\mu\text{m}$  in diameter and are biseriate. The conidia are spherical to oval with a hyaline smooth wall and, as shown in this figure, can form long chains.



**14-34 *Aspergillus versicolor* colony on Sabouraud's dextrose agar.** The colonies of *A. versicolor* usually display a wide range of colors including white-tan to yellow, green, and brown.



**14-35 *Aspergillus versicolor* stained with lactophenol cotton blue ( $\times 1250$ ).** This figure shows Hülle cells produced by *A. versicolor*. These structures are also produced by other species of aspergilli, and it is not clear what their function is. They are thick-walled, spherical or pear-shaped, and are produced as terminal or intercalary cells on the hyphae.




**14-36 *Candida albicans* growing on sheep blood agar plates.** The colonies of *C. albicans* are white or tan, opaque, smooth, and as shown in this figure, may have small extensions or "feet" that increase with the age of the colony. The colonies grow rapidly and can usually be detected in 24 to 48 hours when grown aerobically at 25 to 30°C.




**14-37 *Candida albicans* grown on a corn meal agar with trypan blue.** As shown in this figure, *C. albicans* grown at 25°C for 48 to 72 hours forms septate pseudohyphae with groups of blastospores at the septa, and large thick-walled terminal chlamydospores that are typical of this species.



**14-38 *Candida albicans* germ tubes stained with calcofluor white (×625).** A simple test for the identification of *C. albicans* is the germ tube test. Part of a colony is emulsified in fetal bovine serum and incubated at 37°C for 2 to 3 hours. The germ tubes are approximately half the width and two to four times the length of the yeast cells. The germ tubes produced by *C. albicans* are not constricted at the union with the blastoconidium as is the case with *Candida tropicalis*.




**14-39 *Candida albicans*. Blood culture, Gram stain (×1350).** Blood culture from a patient with a *C. albicans* fungemia. Blastoconidia can be observed budding and producing the pseudohyphae. The blastoconidia stain gram-positive, are round or oval, and measure approximately  $3 \times 5 \mu\text{m}$ .



**14-40 *Candida tropicalis*. Blood culture, Gram stain (×1250).** The blastoconidia of *C. tropicalis* tend to be more barrel-shaped and irregular than those of *C. albicans*.




**14-41 *Geotrichum* sp. grown on corn meal agar with trypan blue.** Chains of smooth, hyaline, cylindrical arthroconidia with round corners, are characteristic of this species. The arthroconidia originate by segmentation of the hyphae, and typically they germinate from one corner, giving the appearance of a "hockey stick." The formation of consecutive arthroconidia differentiates this genus from *Cordiloides immitis*.



**14-42 Diagram of *Zygomycetes* structure.** Typical microscopic features of the *Zygomycetes*, a group with generally aseptate hyphae, are shown here. The organisms of the *Rhizopus* spp. (A) have nodal sporangiophores, e.g., rhizoids (root-like structures) are formed at the point where the sporangiophores meet the stolon (horizontal section of hyphae from which sporangiophores and rhizoids arise). No rhizoids are found in the *Mucor* spp. (B). *Absidia* spp. (C) form internodal sporangiophores. When the wall of the sporangium breaks, a collarette remains at the top of the sporangiophore and the pointed columella becomes visible.




**14-43 *Absidia corymbifera* colony in Sabouraud's dextrose agar.** The colonies of *A. corymbifera* grow rapidly, filling the Petri dish, and have a woolly, white color that turns gray with age. The reverse side is white. In general, all colonies of *Zygomycetes* have similar colony characteristics.




**14-44, 45 Lactophenol cotton blue preparations of *Absidia corymbifera* (14-44,  $\times 80$ ; 14-45,  $\times 550$ ).** *A. corymbifera* produces few rhizoids and, in contrast to those of the genus *Rhizopus*, the sporangiophores grow at a point on the stolon that is between the rhizoids and not opposite them. The sporangiophores are characteristically highly branched (Figure 14-44). The sporangium of *A. corymbifera* has a pyriform shape, small size, measuring 20 to 40  $\mu\text{m}$  in diameter, with a conical columella and a distinct apophysis (Figure 14-45). The sporangiospores are oval in shape, measure approximately  $3 \times 5 \mu\text{m}$  and have a smooth wall.


**14-46 *Mucor* spp. stained with lactophenol cotton blue ( $\times 500$ ).** The sporangiophores of *Mucor* spp. have terminal spherical sporangia that measure up to  $300\ \mu\text{m}$  in diameter, containing numerous round sporangiospores. The sporangiophores are hyaline, smooth, with a gray-tan color, and they do not have apophyses. The columella is well developed and has a prominent collarette at the junction with the sporangiophore. No rhizoids are found in this species.



**14-47 *Mucor* spp. Lung. Gomori's methenamine silver stain ( $\times 1250$ ).** In tissues, the hyphae of *Mucor* spp. are broad, irregular, with bulbous lateral protrusions, and do not have septae; branching occurs preferentially at wide angles ( $90^\circ$ ).



**14-48 *Rhizopus* spp. Lactophenol cotton blue preparation ( $\times 80$ ).** The members of the genus *Rhizopus* are characterized by the presence of sporangiophores that grow opposite the rhizoids along the stolon. The columella is round or slightly elongated, and the apophysis is not obvious.



## Subcutaneous Mycoses



**14-49 Lactophenol cotton blue preparation of *Bipolaris* spp. ( $\times 550$ ).** The conidia of *Bipolaris* spp. are fusiform, approximately  $10 \times 20 \mu\text{m}$  in size, rounded at both ends, with the central cells similar in size and color to the distal ones. 3 to 5 light pseudoseptae, and a nonprotuberant hilum. The hilum is a "scar" at the point of attachment of a conidium to the conidiophore. The conidiophores bend at each point where conidia are formed, giving a zigzag appearance. The hyphae are septate.



**14-51 Colony of *Cladosporium carrionii* growing on Sabouraud's dextrose agar.** The colonies of *C. carrionii* are slow growers, with a flat or slightly raised center, and a velvety gray or green-black color.




**14-50 *Bipolaris* spp. Germ tube preparation ( $\times 575$ ).** *Bipolaris* spp. typically display orientation of the germ tubes along the axis of the conidium. This characteristic can be used to differentiate these organisms from members of the genus *Drechslera*, which form germ tubes perpendicular to the conidial axis. Members of the genus *Exserohilum* also form germ tubes along the axis of the conidium. However, the size and structure of the conidium is significantly different.



**14-52 *Cladosporium carrionii* preparation stained with lactophenol cotton blue ( $\times 550$ ).** The hyphae are dark, septate, and have branches. The conidiophores are elongated and produce chains of ellipsoid, smooth-walled conidia that measure approximately  $2 \times 5 \mu\text{m}$  in size. The conidia have a characteristic dark area at the ends called a disjuncter. The conidia closest to the conidiophore may have a "shield" shape.




**14-53 *Exophiala (Wangiella) dermatitidis*.** Lactophenol cotton blue preparation ( $\times 500$ ). In new cultures, oval and round budding yeastlike cells are formed. Subsequently, these cells produce septate hyphae with flask-shaped to cylindrical phialides. The conidia are round or oval, measure approximately  $3 \times 5 \mu\text{m}$ , and are found at the tip of the phialide and also along the hyphae. Growth and biochemical characteristics are used to differentiate these organisms from *Exophiala jeanselmei* and *Phaeoanellomyces werneckii*.





**14-54 *Exophiala jeanselmei* growing on Sabouraud's dextrose agar.** The colonies are brown or green-black, moist, and glistening. With age they become covered with velvety-grayish hyphae. The reverse is black.




**14-55 Lactophenol cotton blue preparation of *Exophiala jeanselmei* ( $\times 500$ ).** The conidiophores are elongated, tubular, and with a tapered, narrow end. The conidia of *E. jeanselmei* are smooth, thin-walled, and ellipsoid, measuring  $2 \times 3 \mu\text{m}$ , and can gather in clusters around the conidiophores and at points along the septate hyphae.




**14-56 Lactophenol cotton blue preparation of *Exserohilum* spp. ( $\times 500$ ).** The conidia are fusiform, measure approximately  $100 \times 15 \mu\text{m}$ , have a prominent truncated hilum, seven to eleven septa, and the septum next to the hilum is often darkly stained. The conidiophores are elongated and bend at the point of attachment of the conidia (geniculate), giving a zigzag formation. The hyphae are septate and dark.



**14-57, 58** *Fonsecaea pedrosoi* stained with lactophenol cotton blue ( $\times 550$ ). The hyphae of *F. pedrosoi* are septate, branched, and have a dark brown pigment. Three types of conidial formation may be observed. The *Phialophora* type (see *Phialophora* sp.), the *Cladosporium* type, and the *Rhinochladella* type. In the *Cladosporium* type (Figure 14-57) the conidiophores give rise to large shield-shaped cells that produce branching chains of oval conidia with dark scars in the hila. In the *Rhinochladella* type of conidiation (Figure 14-58), the conidiophores, arising terminally or laterally on hyphae, have denticles that produce ovoid or elongated conidia along the sides or the tips.



**14-59** *Phialophora verrucosa* stained with lactophenol cotton blue ( $\times 500$ ). This microorganism has typical vase- or flask-shaped or elliptical phialides, or conidiophores, with wide, flared, pigmented collarettes. The conidia are round or ellipsoid, hyaline, and measure approximately  $2 \times 4 \mu\text{m}$ . Hyphae are septate, brown, and branched.




**14-60** Gomori's methenamine silver stain of chromoblastomycosis ( $\times 1250$ ). Skin lesion from a patient with chromoblastomycosis. Thick-walled dark brown cells (called sclerotic bodies and the hallmark of chromoblastomycosis) that divide by septation can be observed. The following microorganisms are associated with chromoblastomycosis: *Botryomyces caespitosus*, *Cladosporium carrionii*, *Fonsecaea pedrosoi*, *Fonsecaea compacta*, *Phialophora verrucosa*, and *Rhinochladella aquaspersa*.






**14-21** Colony of *Pseudococconeis haptii* on *Suberostrea exilis* eggs. The colonies of *P. haptii* have a sticky surface that is white to gray-brown in color and gets darker with the age of the culture. The reverse is also white (white/light brownish eggs).



**14-22.** **22** **Left:** hyaline cotton-like state of *Pseudococconeis haptii* and *Arthropus penalis*, several states of *Pseudococconeis haptii* (juvenile stages) (1-200x). The *Arthropus* type of spores of *P. haptii* (Figure 14-23) may not develop from the spore to give a fan-like type of arthropods; spores translated at the base, and sometimes branched, cords of *Pseudococconeis haptii*. The hyaline are long and flexible, branch at acute angles and they may include septula. To describe the structure an organism (partial structure) of the *Arthropus* state of *P. haptii* (Figure 14-23) have involved hyaline cords (arthropods) or cylindrical (approximately 6-8  $\mu$ m in the basal end of *haptii* length) (3-220  $\mu$ m in diameter round) forms structures on these consisting of complex structures.


**14-64 Lactophenol cotton blue preparation of *Scedosporium prolificans* ( $\times 1250$ ).** In the annellation type of conidiation, the conidia are extruded from the tapered tip of the conidiophore, which bears a ring scar for each conidium. The conidiophores, or annellides, are short and do not branch. The conidiogenous cells are flask-shaped, with annellation and terminal conidia. The conidia are ovoid to pyriform (pear-shaped), truncated at the base, with a smooth and thin wall, approximately  $3 \times 10 \mu\text{m}$ .



**14-65 Colony of *Sporothrix schenckii* growing on potato dextrose agar.** At  $25^{\circ}\text{C}$  these colonies grow slowly, are moist, and have a wrinkled surface. The color of the surface initially is white-tan and turns dark brown or black with age. At  $37^{\circ}\text{C}$  the colonies on brain heart infusion agar (BHI) are cream-tan, smooth, and yeastlike (not shown).




**14-66 Lactophenol cotton blue stained preparation of *S. schenckii* ( $\times 500$ ).** Growth at  $25^{\circ}\text{C}$  produces round, oval, and pear-shaped conidia that measure  $2 \times 5 \mu\text{m}$  and are attached to both sides of the septate hyphae by short and thin denticles. Clusters of conidia can also be observed in a "daisy" pattern attached by denticles to conidiophores. At  $37^{\circ}\text{C}$ , oval or fusiform budding cells ("cigarlike") can be observed (not shown).



## Superficial Mycoses




**14-67** *Epidermophyton floccosum* on Sabouraud's dextrose agar. Colonies of *E. floccosum* grow slowly at 25°C and often appear white to tan, although darker colors including olive and khaki can appear with age. The surface is fluffy and powdery and can be flat or radially folded. The reverse side may have an orange-brown pigmentation.



**14-68** Preparation of *E. floccosum* stained with lactophenol cotton blue ( $\times 600$ ). *E. floccosum* have club-shaped, smooth, thin-walled macroconidia that usually grow in clusters directly from the septate hyphae. The macroconidia have less than five septa and measure 20 to 40 by 5 to 8  $\mu\text{m}$ . Chlamydospores, both terminal and intercalary, can be observed in old cultures.




**14-69** Colonies of *Malassezia furfur* growing on Sabouraud's dextrose agar. Olive oil was added to the right side of the plate. As shown in this figure, the colonies are only growing on the right side of the plate because of the growth requirement of long-chain fatty acids, provided by the olive oil. The colonies grow better at 37°C and are yellow-tan, smooth, and dry.



**14-70, 71** Gram stain (14-70;  $\times 1250$ ) and calcofluor white (14-71;  $\times 625$ ) preparations of *Malassezia furfur*. These yeastlike cells are actually phialides that have a spherical or ellipsoid shape and measure approximately  $3 \times 5 \mu\text{m}$ . The cells are round at one end and cut off at the other end, with an indistinct collarette, where the budlike structure forms singly on a broad base. In the calcofluor white fluorescence stain, the typical "bowling-pin" morphology can be clearly observed.



**14-72** *Microsporum audouinii* colony on Sabouraud's dextrose agar. The colonies of *M. audouinii* are usually white to tan, flat, and have a suede-like surface. The underside is frequently yellow-red-brown.



**14-73** *Microsporum audouinii* stained with lactophenol cotton blue ( $\times 1800$ ). Cultures of *M. audouinii* are often sterile. The hyphae are septate with intercalary and terminal chlamydospores. On high magnification, the terminal chlamydospores can be shown to have a pointed end.




**14-74, 75** *Microsporum canis* var. *canis* colony growing on Sabouraud's dextrose agar. *M. canis* grows rapidly at 25°C, producing colonies that are usually white-tan with a yellow-green lemon color at the periphery (Figure 14-74). The surface of the colony often has a radiate woolly appearance. The reverse side of the colonies is frequently golden yellow or brown (Figure 14-75).

**14-76** Lactophenol cotton blue preparation of *Microsporum canis* var. *canis* ( $\times 625$ ). The macroconidia of *M. canis* are spindle or fusiform in shape with a thick irregular rough wall containing between 5 to 15 cells. They measure 5 to 20 by 10 to 100  $\mu\text{m}$ , and characteristically have a knoblike end. The hyphae are septate.




**14-77** *Microsporum gypseum* colonies on Sabouraud's dextrose agar. The colonies of *M. gypseum* have a white-tan surface, a white, starburstlike border, and a suedelike appearance. Areas with red-brown color are common on the reverse side. The colonies grow rapidly.






**14-78** *Microsporium gypseum* stained with lactophenol cotton blue ( $\times 625$ ). The macroconidia of *M. gypseum* have an ellipsoidal to fusiform shape with a thin irregular wall, contain four to six cells, and measure 7 to 15 by 30 to 60  $\mu\text{m}$ . The site of attachment to the hyphae is usually flattened, while the distal end is more rounded. The hyphae are septate and the microconidia have a clavate or club shape, but this characteristic is not helpful for differentiation from other fungi that produce similar structures.




**14-79** Colony of *Trichophyton mentagrophytes* growing on Sabouraud's dextrose agar. The colonies grow well at 25°C and are usually cream-tan in color, although they may turn darker with age. The surface of the colony appears fluffy and powdery. The reverse of the colony may have a color ranging from tan to brown and dark red.



**14-80** *Trichophyton mentagrophytes* stained with lactophenol cotton blue ( $\times 500$ ). This fungus typically displays spherical or pyriform microconidia growing singly or in clusters on branched conidiophores. In addition, there are often characteristic septate spiral or coiled hyphae. Cigar-shaped macroconidia with three to six cells and measuring  $7 \times 40 \mu\text{m}$  can sometimes be found (not shown).



**14-81** *Trichophyton rubrum* on Sabouraud's dextrose agar slant. White, velvety, or fluffy colonies of *T. rubrum* grown at 25°C. Note the wine-red pigment produced by this organism. Pigments ranging from yellow to orange to red can be observed in different isolates of *T. rubrum*. This type of pigment, however, is not unique to *T. rubrum*. Other dermatophytes, including *T. ajelloi* and *T. mentagrophytes*, can produce red-orange pigments.




**14-B2** *Trichophyton rubrum* stained with lactophenol cotton blue ( $\times 500$ ). *T. rubrum* has septate hyphae with lateral, tear-shaped microconidia that measure  $3 \times 4 \mu\text{m}$ . The macroconidia are long with thin parallel walls and two to eight cells (not shown). Microconidia may form directly from the macroconidia.




**14-B3** Christiansen urea agar slant for demonstrating urease production by *Trichophyton mentagrophytes* and *Trichophyton rubrum*. The tube on the left side of the image has been inoculated with *T. mentagrophytes*, and the two on the right with *T. rubrum*. The tubes have been incubated for 3, 7, and 3 days, respectively. Note that the *T. mentagrophytes* yields a positive test after 3 days, while the *T. rubrum* is only weakly positive at 7 days.





**14-B4** Colony of *Trichophyton schoenleinii* growing on Sabouraud's dextrose agar. Young cultures of *T. schoenleinii* have a waxy appearance with a tan to brown color and honeycomblike thallus (mat of hyphae). This microorganism grows slowly, and with age the colonies become irregular with folded surfaces, as shown in this figure.



**14-B5** *Trichophyton schoenleinii* stained with lactophenol cotton blue ( $\times 500$ ). The hyphal ends are swollen, giving a "nailhead" morphology, and as a result of branching, often give the appearance of "favic chandeliers" or "antherlike" hyphae. Macroconidia and microconidia are rare, although chlamydoconidia are frequently observed.




**14-86, 87** *Trichophyton tonsurans* colony on Sabouraud's dextrose agar. *T. tonsurans* grows slowly at 25°C on Sabouraud's dextrose agar. The colonies can display a wide variety of colors ranging from white, creamy, yellow, tan, and pink. The colonies can be flat or raised with a velvety or powdery appearance. Note rugae (ruglike folds) cutting across the colony (Figure 14-86). The underside of the colony can range from yellow-brown to red-brown color (Figure 14-87).



**14-88, 89** *Trichophyton tonsurans*. Lactophenol cotton blue preparation ( $\times 500$ ). The septate hyphae have microconidia that may vary significantly in size and shape. Some of them are round or oval with a ballooned appearance, and others are clavate and pear-shaped. The microconidia are attached to branched conidiophores by a short stalk (Figure 14-88). Macroconidia are smooth, thin-walled, irregular, and clavate (Figure 14-89), although frequently they cannot be observed. Intercalary and terminal chlamydospores are common.







**14-90 Nutritional test for *Trichophyton tonsurans*.** For this test, the microorganism is grown on vitamin-free casein agar (T1), and on casein agar with thiamine added (T4). As shown in this figure, *T. tonsurans* requires thiamine for growth.




**14-91 *Trichophyton verrucosum* on Sabouraud's dextrose agar.** This microorganism has slow-growing colonies with a white, creamy color. The surface is velvety, and the center may be raised.




**14-92, 93 Lactophenol cotton blue preparations of *Trichophyton verrucosum*** (Figure 14-92,  $\times 550$ ; Figure 14-93,  $\times 225$ ). On Sabouraud's dextrose agar at 37°C, *T. verrucosum* produces many chlamydospores that are sometimes referred to as "chains of pearls" (Figure 14-92). Occasionally "antlerlike" branches at the ends of the hyphae can be observed, but they are rare compared with those of *T. schoenleinii*. On enriched media with thiamine, small, delicate, single microconidia, and long thin macroconidia shaped like a rat's tail can be found (not shown). In some preparations, hyphae produce terminal vesicles (Figure 14-93).



**14-94 Colonies of *Trichophyton violaceum* on Sabouraud's dextrose agar.** *T. violaceum* produces a waxy, irregular shaped, raised colony with areas that have a dark violet color and others that are white. Occasionally, isolates produce only white-tan colonies. The undersurface is purple. In general, this microorganism is a slow grower.




**14-95, 96 Lactophenol cotton blue (14-95; [ $\times 225$ ]), and calcofluor white (14-96; [ $\times 625$ ]) preparations of *Trichophyton violaceum*.** These figures demonstrate the large, branched, irregular shaped hyphae containing numerous intercalated chlamydoconidia. Macroconidia and microconidia are rarely produced.



## Miscellaneous Mycoses



**14-97** *Fusarium* spp. colony on potato-dextrose agar. *Fusarium* spp. grow rapidly in culture, and on potato-dextrose agar the colonies are cottonlike, usually white, turning pink-violet or brown at the center with age.




**14-98, 99** *Fusarium* spp. stained with lactophenol cotton blue (Figure 14-98,  $\times 700$ ; Figure 14-99,  $\times 500$ ). Typical *Fusarium* spp. microconidia with a fusiform or oval shape extending from delicate lateral phialides (Figure 14-98). The macroconidia of *Fusarium* spp. are produced on conidiophores after 4 to 7 days. The macroconidia are fusiform, usually curved, giving the appearance of a sickle, and have three to five septae (Figure 14-99).




**14-100 Colony of *Paecilomyces variotti* on Sabouraud's dextrose agar.** This organism is a fast grower that produces flat colonies with a tan-brown color and a powdery or suedelike surface.



**14-101 Lactophenol cotton blue preparation of *Paecilomyces variotti* ( $\times 550$ ).** The phialides, or sterigmata, bend away from the axis of the conidiophore, and thus they are called "ten pins." The conidia are elliptical or oblong, measure approximately  $2 \times 3 \mu\text{m}$ , and the chains do not branch.




**14-102 Diagram of *Penicillium* spp. "brush" or "penicillus."** The septate hyphae have branched and unbranched conidiophores. These form metulae (short, hyphal structures below the phialides) that give rise to flask-shaped phialides. The conidia are round, smooth, or rough and unbranched.




**14-103 Colony of *Penicillium* spp. on Sabouraud's dextrose agar.** The colonies of *Penicillium* spp. usually grow fast and have a powdery white, gray, or green surface color.




**14-104 *Penicillium* spp. stained with calcofluor white ( $\times 250$ ).** Fructing head of *Penicillium* sp. showing septate mycelia with conidiophores. The conidiophores measure 100 to 250  $\mu\text{m}$  and consist of phialides, or sterigmata, that extend directly from the conidiophore. Alternatively, as in this figure, they originate from metulae, giving a brushlike appearance, also known as a "penicillus." Extending from the tapered tip of the phialides are short, unbranched chains of conidia, measuring  $3 \times 5 \mu\text{m}$ , that can be spherical or fusiform and smooth or rough-walled.




**14-105 Colonies of *Penicillium marneffei* growing on Sabouraud's dextrose agar.** The *P. marneffei* colonies have a white-gray color that turns green with age and produces a characteristic dark red pigment when grown at 25°C.



**14-106, 107** Lactophenol cotton blue (14-106;  $\times 500$ ) and calcofluor white (14-107;  $\times 750$ ) preparations of *Penicillium marneffei*. At 25°C the conidiophores of *P. marneffei* are smooth and have three to five metulae, each of them with several phialides, producing smooth, spherical conidia in chains (Figure 14-106). The conidia measure approximately  $3 \times 6 \mu\text{m}$ . At 37°C, oval and elliptical yeastlike cells that measure  $3 \times 7 \mu\text{m}$  are produced. The cells replicate by fission and, as shown here, a distinct cross-wall is formed (Figure 14-107).



**14-108** Colony of *Scopulariopsis* spp. on Sabouraud's dextrose agar. This organism grows fast, producing colonies that vary in color from white to tan or brown and black. In most instances however, as in this case, the colony is tan. The surface at first is glabrous (smooth) and with time usually becomes powdery.




**14-109** *Scopulariopsis* spp. stained with lactophenol cotton blue ( $\times 550$ ). This organism typically produces chains of single-celled conidia originating from a conidiogenous annellide. The conidia are round or pyriform, measuring 7 to 8  $\mu\text{m}$  in diameter, with a thick wall that may be smooth or rough, giving the impression of "light globes," and a color that ranges from tan to brown. The hyphae are septate.

## Specimen Preparation and Identification Systems

Colonies of molds isolated in culture can be examined microscopically using several techniques, including the tease mount, the tape mount technique, and the slide culture method. These procedures should always be performed in a laminar flow biosafety cabinet.

**14-110 Tease mount technique.** For the tease mount, a fragment of the colony is collected using a wire or a loop and transferred to a glass microscope slide. A drop of lactophenol cotton blue is then added, and the specimen is teased using dissecting needles as shown in this figure. The teasing of the specimen needs to be done carefully so that on the one hand, isolated elements can be observed, while at the same time preserving the integrity of the overall structure of the microorganism.



**14-111 Adhesive tape technique.** A piece of transparent adhesive tape is used to collect the specimen by pressing the adhesive side against the surface of the fungal colony. Aerial elements will adhere to the tape, which is subsequently placed on a microscope slide containing a drop of lactophenol cotton blue. This technique is good for preserving the original relationship between spores and aerial hyphae. However, it usually cannot be applied to mold specimens that have few aerial mycelia or to yeast with a moist consistency.



**14-112 Slide culture method.** Two sterile rods are placed at the bottom of a Petri dish on which a sterile glass microscope slide is placed. Blocks or circles of agar are transferred aseptically to the microscope slide. A fragment of the fungal colony to be studied is inoculated onto the sides of the agar, which is coverslipped aseptically and incubated at 25°C. When the colony is mature, the coverslip is removed and both the coverslip and the growth on the slide surface below the agar block are mounted in lactophenol cotton blue and sealed with nail polish or mounting medium.





**14-113 API 20C clinical yeast system.** The API 20C system (bioMérieux Vitek, Inc., Hazelwood, Mo.) is a micromethod for the identification of most yeasts and yeastlike microorganisms. Microcupules containing carbohydrate substrates are inoculated with a suspension of the microorganism, and the strip is incubated at 30°C. Once growth occurs, cupules showing turbidity heavier than the O control cupule are considered to be positive. The results are converted to a seven-digit biotype, which is matched to the analytical profile index published by the manufacturer.



**14-114 VITEK YBC (yeast biochemical card).** The VITEK YBC (bioMérieux Vitek, Inc., Hazelwood, Mo.) card is a component of a semiautomated identification system for yeast and yeastlike organisms. The card contains 30 wells, of which four are negative controls. The biochemical tests include several conventional tests such as carbohydrate assimilation, urea hydrolysis, resistance to cycloheximide, and nitrate reduction. The YBC is incubated at 30°C in the VITEK reader/incubator module, and read at 24 and 48 hours. The automated system provides an identification based on the results of the different reactions.



**14-115 DNA probes.** Nucleic acid probes are now available commercially for the identification of *Blastomyces dermatitidis*, *Coccidioides immitis*, *Cryptococcus neoformans*, and *Histoplasma capsulatum*. To perform this test, growth from a colony is treated to extract the DNA, which is then hybridized to a probe provided in the kit (Gen-Probe Inc., San Diego, Calif.). One of the main advantages of this system is that the test can be performed before sporulation has occurred.




## CHAPTER 15 Parasitology


**P**articularly in areas of poor hygienic practices and sanitation, parasites are still among the major causes of disease and death around the world. Citizens of developed nations often acquire the microorganisms during travel. Direct examination of diluted material is usually the best method for laboratory diagnosis. Careful observation by skilled macrobiologists is essential. Serological requests and tests for some parasites are available commercially and at reference laboratories. This section includes images of protozoa (single-celled microorganisms), followed by helminths (multi-celled macroorganisms including worms). Helminths are divided into roundworms, called nematodes; flatworms including flukes, called trematodes; and tapeworms, called cestodes. Average size of the microorganisms is indicated in the legends.

### Protozoans


#### Intestinal protozoan: *Amoeba*.




**15-1** *Entamoeba histolytica*, Trophozoites, Parasit. *Trichrome stain* ( $\times 800$ ). The trophozoites of *E. histolytica* are round and measure from 12 to 20  $\mu\text{m}$  in diameter. With the trichrome stain the cytoplasm stains pink-orange, and the karyosome a central large mass. The organism has a delicate cytoplasm and a distributed uniformity in organization along the entire periphery.



**15-2** *Entamoeba histolytica*, Trophozoites, Parasit. *Trichrome stain* ( $\times 800$ ). This *E. histolytica* trophozoite shows an ingested red blood cell within cytoplasm. The nuclear structure is still fairly typical although the karyosome is eccentrically located. Ingestion of red blood cells by *E. histolytica* is very characteristic with lamellar bodies visible inside, but only a few red cells are still visible within the organism from the other species of *Entamoeba* in nature. In red cells karyome and cell debris are visible in the cytoplasm of this parasite. Macrophages that have ingested white material may sometimes be confused with *E. histolytica*.




**15-3 *Entamoeba histolytica*. Trophozoites. Liver. PAS stain ( $\times 1250$ ).** Amebic liver abscesses are relatively common complications in cases of invasive disease. The chromatin of this trophozoite is well defined and located along the nuclear membrane. The karyosome is not clearly visible. The morphology of the parasite is usually not well preserved in this type of specimen due to the autolysis of the tissues.




**15-4 *Entamoeba histolytica*. Cyst. Feces. Trichrome stain ( $\times 1500$ ).** Immature cysts of *E. histolytica* typically have a large nucleus that can be laterally displaced by a glycogen vacuole. The karyosome may be, as in this case, large and eccentrically located, although in general, they are small and centrally located. Several chromatoid bodies can be found surrounding the glycogen vacuole. The cysts of *E. histolytica* range in size from 12 to 16  $\mu\text{m}$ .




**15-5 *Entamoeba histolytica*. Cyst. Feces. Trichrome stain ( $\times 1250$ ).** Mature cysts of *E. histolytica* have four nuclei while immature cysts contain one or two nuclei. The cysts range in size from 10 to 18  $\mu\text{m}$  and frequently have a diffuse glycogen vacuole, as in this case. Red chromatoid bodies with rounded ends are frequently found. The karyosomes may be slightly eccentrically located, and the peripheral chromatin is fine and evenly distributed.




**15-6 *Entamoeba hartmanni*. Trophozoite. Feces. Trichrome stain ( $\times 1250$ ).** The trophozoite of *E. hartmanni* is round or ovoid and measures approximately 7 to 10  $\mu\text{m}$ . With the trichrome stain, the cytoplasm appears green-blue, and the single nucleus demonstrates a karyosome that is usually centrally located, although it can be laterally displaced. The chromatin is peripheral and has a fine appearance. This microorganism does not phagocytize red cells, a characteristic that is useful to distinguish it from *E. histolytica*.




**15-7 *Entamoeba hartmanni*. Cyst. Feces. Trichrome stain ( $\times 1250$ ).** The cysts of *E. hartmanni* appear similar to those of *E. histolytica*. Although they are generally smaller in size, 5 to 8  $\mu\text{m}$ , there is a significant amount of overlap, so they are difficult to differentiate. The mature cysts have four nuclei, but the immature cysts with one or two nuclei are more frequently found. The nucleus has a small centrally located karyosome with a uniformly distributed peripheral chromatin. In this particular figure, several typical red-staining chromatoid bodies can be observed.




**15-8 *Entamoeba coli*. Trophozoite. Feces. Trichrome stain ( $\times 1250$ ).** The trophozoites of *E. coli* are big, measuring 20 to 25  $\mu\text{m}$  in diameter and are round or amoeboid in shape. The single nucleus has a large, eccentrically located karyosome, and the chromatin has a clumpy, irregular appearance that stains green to purple with the trichrome stain. The cytoplasm of this microorganism usually contains bacteria, yeast, and other cell debris. In this figure, a glycogen-like vacuole can be observed in the cytoplasm.




**15-9 *Entamoeba coli*. Cyst. Feces. Trichrome stain ( $\times 2500$ ).** The cysts of this species are large, usually around 15 to 25  $\mu\text{m}$  in diameter, and while the immature cyst has one to two nuclei, the mature stage contains eight nuclei with a distinct, eccentric karyosome. This figure shows four of the eight nuclei in the same plane of focus. The peripheral chromatin is usually coarse, although it may have a smooth appearance. The chromatoid bodies are usually splintered with sharp ends, as shown here.




**15-10 *Endolimax nana*. Trophozoite. Feces. Trichrome stain ( $\times 2500$ ).** The trophozoites of *E. nana* are small, measuring approximately 8 to 10  $\mu\text{m}$  in diameter. The microorganism is spherical to amoeboid in shape and has a cytoplasm that frequently contains granules, vacuoles, and bacteria or other cellular debris. The nucleus is small with a relatively large, red-purple, centrally located karyosome. Trophozoites with irregular karyosomes are frequently found. Typically no marginated chromatin can be observed in the nucleus. This combination of a large karyosome surrounded by a clear halo gives an "owl's-eye-like" appearance to the nucleus of this species.




**15-11** *Endolimax nana*. Cyst. Feces. Trichrome stain ( $\times 1250$ ). The mature cysts measure 7 to 8  $\mu\text{m}$  in diameter, have four nuclei with a prominent karyosome and no peripheral chromatin. No chromatoid bodies are found in the cytoplasm of these microorganisms, although ingested debris and bacteria can be observed.



**15-12** *Iodamoeba bütschlii*. Trophozoite. Feces. Trichrome stain ( $\times 2500$ ). The trophozoites of this species are round or ovoid and measure approximately 12 to 15  $\mu\text{m}$ . The nucleus has a large, reddish, centrally located karyosome with no peripheral chromatin, giving these microorganisms a structural similarity to *E. nana*. Achromatic granules can sometimes be observed between the karyosome and the nuclear membrane.




**15-13** *Iodamoeba bütschlii*. Cyst. Feces. Trichrome stain ( $\times 2500$ ). The cysts of *I. bütschlii* have a single nucleus with a large eccentric karyosome and no peripheral chromatin. Typically the cysts contain large glycogen vacuoles that stain brown with iodine. In this figure, the glycogen is displacing the nucleus to the periphery of the cytoplasm. Cysts range in size from 10 to 15  $\mu\text{m}$  in diameter.




**15-14** *Blastocystis hominis*. Feces. Trichrome stain ( $\times 2000$ ). This microorganism is oval or spherical and ranges in size from 5 to 35  $\mu\text{m}$ . The central area resembles a vacuole that takes a green color with the trichrome stain. The cytoplasm is located in the periphery and usually contains one nucleus, although two to four can be found. Large granules with a dark red color can also be found in the periphery of the cytoplasm.


## Flagellates and Ciliates




**15-15 *Dientamoeba fragilis*. Trophozoite. Feces. Trichrome stain ( $\times 2500$ ).** *D. fragilis* trophozoites resemble amoeba and have a spherical or oval shape. They measure 8 to 12  $\mu\text{m}$ , and the cytoplasm frequently contains bacteria, yeast, and other types of debris, giving it a granular appearance. These microorganisms have one or two nuclei with a characteristic karyosome composed of four to eight lobules. No peripheral chromatin can be identified. The small size of *D. fragilis* and the faint structural staining characteristics make its detection difficult. There is no cyst stage in this species.




**15-16 *Giardia lamblia*. Trophozoite. Feces. Trichrome stain ( $\times 2500$ ).** The trophozoites of *G. lamblia* have a characteristic pear shape, measuring 10 to 20  $\mu\text{m}$  in length by 10 to 15  $\mu\text{m}$  in diameter. The two nuclei contain a karyosome surrounded by a clear halo and are located one on each side of the axonemes. The karyosomes may appear distinct or may be fragmented. The axonemes extend into eight flagella, four located laterally, two in the ventral region, and two in the caudal section of the microorganism. The dark staining median bodies located below the two nuclei give these microorganisms the appearance of a smiling face. On a lateral view of the microorganisms it is possible to see the sucking disk, used for attachment to the mucosa, that gives the trophozoite the appearance of a "flying saucer" (not shown).




**15-17 *Giardia lamblia*. Cyst. Feces. Trichrome stain ( $\times 2000$ ).** The cysts have an oval shape measuring 10 to 15  $\mu\text{m}$  in diameter. The immature cysts have two nuclei that develop into four nuclei in the mature microorganism. These nuclei are usually margined toward the broader section of the cyst and have prominent karyosomes. Axonemes and fibrils are found in the cytoplasm. As in this figure, retraction of the wall in fixed specimens may give the appearance of a clear halo around the microorganism.




**15-18 *Chilomastix mesnili*. Trophozoite. Feces. Trichrome stain ( $\times 2500$ ).** This microorganism is pear-shaped with a pointed end and measures 10 to 25  $\mu\text{m}$  in length by 10 to 15  $\mu\text{m}$  in width. The nucleus in this figure is anteriorly located and has a poorly defined karyosome. The karyosome can be centrally or peripherally placed, and the chromatin is uniformly distributed. A cytosome, a mouth-like structure, may occupy up to one half of the microorganism.




**15-19 *Chilomastix mesnili*. Cyst. Feces. Trichrome stain. ( $\times 1500$ ).** The cyst of this organism has a lemon shape with a hyaline knob at one end. Fibrils present along the side of the cytostome may give the cyst the appearance of an open safety pin. These cysts measure approximately 7 to 9  $\mu\text{m}$  in diameter and are sometimes misidentified as amebas.



**15-20 *Trichomonas hominis*. Trophozoite. Feces. Trichrome stain ( $\times 2000$ ).** The trophozoites of *T. hominis* have a pyriform shape, measuring 15 to 20  $\mu\text{m}$  in length by 10 to 15  $\mu\text{m}$  in width. The cytoplasm, which may contain granules and a large nucleus, with evenly distributed chromatin, is located at the broad end of the microorganism. An undulating membrane and three to five flagella, four anterior and one posterior, give this microorganism an erratic, brisk movement in wet mount preparations.



**15-21 *Balantidium coli*. Trophozoites. Colon biopsy. Hematoxylin and eosin ( $\times 225$ ).** Specimen from a human colon containing trophozoites of *B. coli*. The trophozoites are oval with a tapering end, and measure 40 to 50  $\mu\text{m}$  in diameter by 50 to 100  $\mu\text{m}$  in length. The macronucleus is clearly visible in several microorganisms, and the cytostome can also be observed in some of the trophozoites. In viable specimens, the cilia are in constant movement.



**15-22 *Balantidium coli*. Cyst. Feces. Carmine stain ( $\times 1200$ ).** Cysts are oval to spherical, measuring 50 to 70  $\mu\text{m}$  in diameter. A thick refractive wall can be easily observed. A kidney bean-shaped macronucleus is very distinct in the middle of a fairly uniformly stained cytoplasm in this preparation. An indentation corresponding to the cytostome can be observed on the same side where the macronucleus is located.


## Coccidia and Microsporidia





**15-23 *Isospora belli*. Cyst. Feces. Iodine stain ( $\times 1200$ ).** The oocysts are the form more frequently found in the stool; trophozoites are rarely seen. The cyst has the shape of a football measuring 25 to 35  $\mu\text{m}$  in length by 10 to 20  $\mu\text{m}$  in width. The hyaline wall is double-layered, refractive, and is clearly visible in this preparation. In stool specimens, however, the wall may be difficult to observe. A large, centrally located spherical sporoblast is fairly prominent. Mature oocysts contain two sporocysts with four sporozoites each (not shown).




**15-24 *Isospora belli*. Cyst. Feces. Acid-fast stain ( $\times 600$ ).** In acid-fast stained preparations, the body of *I. belli* may appear with an ellipsoid halo around a prominent red-stained spherical sporoblast, as shown here. In other instances, the red stain precipitates along the hyaline wall.



**15-25 *Isospora belli*. Cyst. Feces. Wet preparation, Autofluorescence ( $\times 1250$ ).** One of the properties that can be used to identify *I. belli* in stool preparations is the fact that this microorganism autofluoresces under UV light. The hyaline wall is particularly prominent, and the sporoblasts can also be identified.




**15-26, 27 *Cryptosporidium parvum*. Oocysts. Feces. Acid-fast (15-26;  $\times 1250$ ) and fluorescence stain (15-27;  $\times 1250$ ).** In acid-fast stained preparations, the oocysts of *C. parvum* appear round to oval, measuring approximately 4 to 6  $\mu\text{m}$  in diameter. Several black granules can be observed, usually in the periphery of the cyst. Under fluorescence, the wall of the cysts may appear smooth or wrinkled.



**15-28, 15-29 *Microsporidium* spp. Feces. Modified trichrome (15-28;  $\times 1250$ ). Calcofluor white stain (15-29;  $\times 1250$ ).** The spores of *Microsporidium* spp. are round or oval and measure 1 to 3  $\mu\text{m}$ . Diagonal bands crossing the cell and corresponding to the polar tubule can sometimes be observed.

### Tissue Protozoa




**15-30 *Pneumocystis carinii*. Sputum. Gram stain ( $\times 1250$ ).** Both the trophozoite and the cyst form can be observed in respiratory specimens. The microorganisms multiply to form clusters that fill alveolar spaces and block air exchange. The trophozoites are round to oval, measure approximately 5  $\mu\text{m}$  in diameter, and have a nucleus that can be observed with certain stains such as Giemsa.





**15-31 *Pneumocystis carinii*, Bronchoalveolar lavage. Monoclonal antibody stain ( $\times 2000$ ).** Often called the "honeycomb" structure, this cluster consists of cysts that measure 5 to 8  $\mu\text{m}$  in size, and when mature, can contain up to eight trophozoites.




**15-32 *Pneumocystis carinii*. Lung biopsy. Gomori's methenamine silver stain ( $\times 1250$ ).** The cup-shape structure of collapsed cysts is clearly visible in this preparation.

**15-33 *Toxoplasma gondii*. Trophozoites. Culture. Giemsa stain ( $\times 1250$ ).** Trophozoites of *T. gondii* have a crescent shape, measuring 5 to 8  $\mu\text{m}$  in length by 2 to 3  $\mu\text{m}$  in width. The nucleus is relatively large and appears to be centrally located.




**15-34 *Toxoplasma gondii*. Cyst. Brain. Hematoxylin and eosin stain ( $\times 1250$ ).** Infections with *T. gondii* in humans can result in the formation of cysts in different tissues, including the brain. The cysts can range in size from 5 to 50  $\mu\text{m}$  and may contain up to several hundred microorganisms. The tissue stage of the microorganism is called cystozoite or bradyzoite.




## Malaria and *Babesia* spp.


**15-35-15-51 *Plasmodium falciparum*, *Plasmodium malariae*, *Plasmodium ovale*, and *Plasmodium vivax*. Blood films, thick and thin preparations. Giemsa stain ( $\times 1250$  except as noted).** Malaria is transmitted to humans by female *Anopheles* spp. mosquitoes. The inoculated sporozoites travel by the bloodstream to the liver where they infect the hepatocytes. After a period of 2 to 3 weeks, they mature into schizonts that release merozoites into the blood and infect the erythrocytes. Inside the red blood cells, the parasite may follow asexual (producing merozoites) or sexual (producing gametocytes) development. The asexual replication cycle repeats itself approximately every 48 hours in the case of *P. falciparum*, *P. vivax*, and *P. ovale*, and every 72 hours in the case of *P. malariae*. Once inside the red blood cells, merozoites form trophozoites. The trophozoites mature into schizonts that rupture and release merozoites that infect more red blood cells. Gametocytes, on the other hand, represent the sexual stage, which is infectious to mosquitoes. In the mosquitoes, the gametocytes develop into male and female gametes, which undergo fertilization and mature into sporozoites. All stages of the growth cycle of *P. malariae*, *P. vivax*, and *P. ovale* can be found in the circulating blood. In infections due to *P. falciparum* in the peripheral blood, only ring forms and the gametocytes with the pathognomonic banana, crescent, or half-moon shape are usually observed.




**15-35 *Plasmodium falciparum*, ring forms.** This species, which causes the most serious disease, often produces heavy parasitemia with double rings, and "headphone" forms. The microorganism does not cause enlargement of the cytoplasm of the parasitized erythrocyte, but it may induce formation of large, purplish Maurer's dots, as seen in the central red blood cell in this preparation.




**15-36 *Plasmodium falciparum*, ring form.** Visible is a "headphone" structure created by two dots of chromatin on the same ring.




**15-37 *Plasmodium falciparum*.** A later trophozoite stage is on the right. This stage is not usually observed on specimens from patients with a *P. falciparum* infection. However, it is frequently found in cases of infection with the other three species of *Plasmodium*.




**15-38 *Plasmodium falciparum*. Thick blood smear.** Multiple ring forms are visible.




**15-39 *Plasmodium vivax*.** The infected red blood cells can be enlarged up to twice their normal size when infected by *P. vivax* or *P. ovale*. Eosinophilic cytoplasmic stippling, called Schüffner's dots, are also present in erythrocytes infected with these two species. This infected red blood cell is enlarged and displays Schüffner's dots.




**15-40 *Plasmodium vivax*.** Infected cells show fine stippling of Schüffner's dots around the edges and the typical heavy chromatin dot.




**15-41 *Plasmodium vivax*.** Although double ring forms are suggestive of *P. falciparum*, the large cell size and obvious Schüffner's dots help to identify this infection as *P. vivax*.




**15-42 *Plasmodium ovale*.** The enlarged size and irregular edge of this infected erythrocyte are typical of *P. ovale* infection. Schüffner's dots, often present with *P. ovale* infections, are absent here.




**15-43 *Plasmodium ovale*.** Infected erythrocytes may develop a fimbriated edge, shown here.




**15-44 *Plasmodium malariae*.** Band-form trophozoite, seen almost exclusively in *P. malariae*. Schüffner's dots are absent in infections with this species.



**15-45 *Plasmodium malariae*.** Schizonts in a typical rosette formation. The schizonts of this species typically contain six to 10 merozoites. Hemazoin pigment, brown in color, is also visible. The pigment is present in schizonts of all four species but is most prominent in *P. malariae*.




**15-46 *Plasmodium vivax*.** Thick blood smear ( $\times 3500$ ). Mature schizonts of *P. vivax* have on the average 16 merozoites, although the number can range from 12 to 24. The pigment is usually golden-brown and not very prominent.





**15-47 *Plasmodium ovale*.** Thick blood smear ( $\times 3500$ ). Although schizonts of this species develop up to 16 merozoites, they can reveal fewer numbers during early stages of development. Hemazoin pigment is most difficult to see in schizonts of this species.





**15-48 *Plasmodium falciparum*.** Banana-shaped gametocyte seen only in infections due to this species. Peripheral blood reveals only gametocytes and ring forms in *P. falciparum* infections. At the gametocyte stage, the red cell membrane may be invisible.



**15-49 *Plasmodium falciparum*.** Thick blood smear. Crescent-shaped gametocyte displays prominent brown pigment.




**15-50, 51 *Plasmodium vivax*, macrogametocytes.** The micro- and macrogametocytes of *P. vivax*, *P. malariae*, and *P. ovale* are large, oval to round bodies. The chromatin of the macrogametocyte is usually more basophilic than that of the microgametocyte.




**15-52, 53 *Babesia microti*. Blood. Giemsa's stain ( $\times 1250$ ).** *B. microti* is transmitted in the United States by the tick *Ixodes scapularis*, the same vector that transmits *Borrelia burgdorferi*. The intraerythrocytic form of *Babesia* spp. are round to oval and measure 1 to 5  $\mu\text{m}$  in length. This microorganism reproduces by asexual budding into two to four daughter cells. As shown in these two figures, *B. microti* infects RBCs and produces small, ringlike forms, with a scant amount of cytoplasm and a minute chromatin dot. Some of the RBCs show the characteristic tetrads. In some preparations, these tetrads may appear with a "Maltese cross" configuration. This microorganism should be differentiated from the agent of malaria—in particular, *P. falciparum*. In addition to the presence of the tetrads, in *B. microti* infections, extracellular merozoites can be found, while the brown pigment deposits of hemozoin observed in *Plasmodium* spp. infections are absent.

## Leishmaniasis



**15-54 *Leishmania* spp. Promastigotes. Culture. Giemsa stain ( $\times 1250$ ).** Promastigotes of *Leishmania* spp. are found in the gut of the sandfly where they replicate and subsequently migrate to the proboscis. The vertebrate host is infected at the time of feeding. As shown in this figure, the promastigotes are cigar-shaped, measuring 10 to 12  $\mu\text{m}$  in length, and have a nucleus in the center of the body. The rodlike kinetoplast is located in the anterior part. The flagellum extends from the anterior end, and an undulating membrane can be observed at this stage of development.




**15-55 *Leishmania* spp. Amastigotes. Blood preparation. Giemsa stain ( $\times 1250$ ).** The amastigotes of *Leishmania* spp. are oval and measure 4 to 5  $\mu\text{m}$  in length by 2 to 3  $\mu\text{m}$  in diameter. A dark-staining kinetoplast can be observed close to the nuclei in some of the microorganisms. This is the only stage found in humans. This microorganism can be confused with *Histoplasma capsulatum*, *Toxoplasma gondii*, and *Trypanosoma cruzi*.

## Trypanosomiasis




**15-56, 57 *Trypanosoma cruzi*. Blood (15-56) and culture (15-57). Giemsa stain ( $\times 1250$ ).** Typical C-shape trypanomastigote of *T. cruzi* (15-56). This microorganism is the cause of the American trypanosomiasis, or Chagas' disease. It measures approximately 15 to 20  $\mu\text{m}$  in length, has a central nucleus, and a conspicuous kinetoplast. The free flagellum measures 5 to 10  $\mu\text{m}$ , and the undulating membrane is not as prominent as the one in the African trypanosomes. Figure 15-57 demonstrates trypanomastigotes of *T. cruzi* growing in NNN medium.





**15-58 *Trypanosoma brucei gambiense*.** Mouse blood. Giemsa stain ( $\times 1250$ ). *T. brucei gambiense* and *T. brucei rhodesiense* are the causative agents of African trypanosomiasis, or sleeping sickness. The trypomastigotes of *T. brucei gambiense* have a nucleus, a kinetoplast located at the blunt posterior end, and an undulating membrane with a flagellum. In general, the trypomastigotes are spindle-shaped and measure 20 to 30  $\mu\text{m}$  in length.

**Other Protozoa**




**15-59 *Acanthamoeba* sp. Cyst from agar culture. Wet mount ( $\times 500$ ).** The cysts of *Acanthamoeba* sp. are spherical and measure approximately 15 to 18  $\mu\text{m}$  in diameter. Typically the wall has two layers. The outer wall is wrinkled, while the internal wall may be smooth or may appear polygonal or spherical. The trophozoites of *Acanthamoeba* are large, measuring 20 to 40  $\mu\text{m}$  in diameter, and have thin extensions called acanthopodia (not shown). The cytoplasm appears irregular and contains different types of vacuoles. The karyosome is large and centrally located in the nucleus. No peripheral chromatin can be observed in this microorganism.



**15-60, 61 *Acanthamoeba* sp. Cysts from agar culture. Trichrome stain (15-60;  $\times 1250$ ). Calcofluor white (15-61;  $\times 500$ ).** With trichrome stain the membrane of *Acanthamoeba* sp. appears green-cyan with the typical wrinkled appearance. The cytoplasm stains red, and the nucleus has a prominent karyosome surrounded by a clear halo. With the calcofluor white stain, the typical wrinkled appearance of the cyst wall is very distinct.

## Nematodes

### Intestinal nematodes



**15-62 *Ascaris lumbricoides*. Fertile egg. Feces, iodine stain ( $\times 900$ ).** The typical fertile *A. lumbricoides* egg is yellow-brown, has a thick, mammillated shell, is oval in shape, and measures approximately  $60 \times 40 \mu\text{m}$ . These eggs may lack the external mammillated layer and can be differentiated from hookworm and pinworm eggs by their size and shape.




**15-63 *Ascaris lumbricoides*. Infertile egg. Feces, iodine stain ( $\times 600$ ).** The infertile eggs of *A. lumbricoides* are oval, large (approximately  $90 \times 45 \mu\text{m}$ ), lack one or more of the shell layers, and the contents have a heterogeneous appearance due to the presence of fat globules and refractive granules.




**15-64, 65 *Enterobius vermicularis*. Female (15-64) and male (15-65) adult worm.** The female adult worm of *E. vermicularis* (pinworm) is white and measures from 8 to 15 mm in length by approximately 0.4 mm in width. In contrast, the male adult worm is much smaller, measuring 2 to 3 mm in length by 0.1 to 0.2 mm in width. The dilated cephalic region is similar in both sexes, while the tail is pointed in the females and blunted in the males. Females are more frequently found in cellophane tape preparations than the males. The cephalic inflation of the cuticle and the muscular and bulbous portions of the esophagus separated by a narrow region can be observed in Figure 15-65.





**15-66 *Enterobius vermicularis*. Eggs. Feces. Iodine stain ( $\times 600$ ).** Occasionally the eggs of *E. vermicularis* can be found in the feces. The eggs are oval in shape, with one side flattened, and measure in the range of  $50 \times 25 \mu\text{m}$ . When the eggs are laid, they are partially embryonated and, as shown in this figure, no larvae can be observed. Eggs reach the infectious stage 5 to 10 hours after they are laid.





**15-67 *Enterobius vermicularis*. Eggs. Cellophane tape preparation ( $\times 225$ ).** The females lay eggs at night in the perianal region. Thus, the best time to collect the specimen with cellulose tape is early in the morning. The shell of the eggs is thick and hyaline. In some of the eggs, it is possible to observe the larvae (not shown).

**15-68 *Trichuris trichiura*. Male adult worm ( $\times 7.5$ ).** The female worms measure 40 to 50 mm, while the male worms usually range from 35 to 45 mm. In the male, the posterior end is coiled, while in the female it is straight. Both sexes have a whiplike overall shape with a slender anterior end and a thicker, short posterior region.




**15-69 *Trichuris trichiura*. Egg. Feces. Iodine stain ( $\times 800$ ).** The morphology of the *T. trichiura* eggs is fairly characteristic. The eggs are oval in shape, measuring  $50 \times 25 \mu\text{m}$ , with a thick hyaline wall that has two distinct mucoid condensations, or "plugs," at each end. The eggs are not embryonated in the stools. Embryonated eggs may be found in the soil 2 to 3 weeks after they are passed.






**15-70, 71 *Necator americanus*. First-stage rhabditiform larvae. Iodine stain (15-70;  $\times 100$ , and 15-71  $\times 500$ ).** The first-stage rhabditiform larvae of *Necator americanus* measure approximately 220 to 280  $\mu\text{m}$  in length by 16 to 18  $\mu\text{m}$  in width. The buccal canal (Figure 15-71) is long, while the genital primordium is small and difficult to see.




**15-72 Hookworm egg. Feces. Iodine stain ( $\times 900$ ).** The eggs produced by *Necator americanus* and by *Ancylostoma duodenale* are indistinguishable. The eggs are oval, large (70  $\times$  40  $\mu\text{m}$ ), and have a thin shell. These embryonated eggs have from four to eight cells at the time they are passed in the stool, and they can develop first-stage larvae in a day or two when left at room temperature. This is in contrast to the first-stage larvae of *Strongyloides stercoralis* that are found directly in the feces.



**15-73 *Strongyloides stercoralis*. First stage rhabditiform larvae. Feces. Iodine stain ( $\times 75$ ).** The first-stage larvae are found directly in the feces and are the diagnostic stage of this microorganism. They measure 200 to 400  $\mu\text{m}$  in length by 15 to 20  $\mu\text{m}$  in width. The parasitic adult female measures approximately 2 to 3 mm in length by 30 to 40  $\mu\text{m}$  in diameter. The pointed tail in the female is straight, and eggs can usually be observed in the genital tract. Females are parthenogenetic, and parasitic males do not exist. The male is smaller than the female worm and can be identified by its curved, pointed tail. The parasitic female produces embryonated eggs in the mucosal epithelium of the small intestine. The embryonated egg hatches in the mucosa of the epithelium, and the first-stage larvae migrate to the lumen and are passed in the feces. Subsequently, these larvae mature into the third-stage larvae that are infectious.




**15-74, 75 *Strongyloides stercoralis*. First-stage rhabditiform larvae. Feces. Iodine stain (15-74,  $\times 1250$ ; 15-75,  $\times 500$ ).** The first stage larvae of *S. stercoralis* have a short buccal canal (as shown here) in contrast to the hookworm larvae that have a long buccal canal. In addition, in the *S. stercoralis* first-stage larvae, the genital primordium is quite apparent (Figure 15-75), while in the hookworm larvae, the genital primordium cannot be seen.



**15-76, 77 *Strongyloides stercoralis*. Filariform larvae. Wet mount (15-76,  $\times 75$ ). Rhabditiform larvae. Sputum. Gram stain (15-77,  $\times 500$ ).** In immunocompromised patients, hyperinfection may result in multiplication of *S. stercoralis* in the intestinal tract, invasion of the wall, and migration of the third stage larvae to all the tissues. In this case the patient had invasion of the lung, and the larvae were observed on a wet mount preparation and on a Gram stain of the sputum. The filariform larva measures 400 to 500  $\mu\text{m}$  in length, the tail is notched (not shown), and the esophagus occupies close to half of the body length.

## Tissue nematodes

**15-78 *Trichinella spiralis*. Larvae, Muscle. Hematoxylin and eosin stain ( $\times 225$ ).** The adult female produces larvae in the intestinal mucosa that migrate to the muscle, where they become encapsulated by the host tissues. These larvae measure up to 1 mm by the time they mature, but initially they are approximately 100  $\mu\text{m}$  in length by 5  $\mu\text{m}$  in diameter.






**15-79, 80 *Toxocara canis*. Embryonated egg (15-79,  $\times 560$ ) and egg with larvae (15-80,  $\times 700$ ). Dog feces. Iodine stain.** The eggs of *T. canis* are spherical or oval, measuring  $85 \times 75 \mu\text{m}$ , with a thick and pitted shell. These eggs are passed in the feces by dogs. Persons coming into contact with these feces, or contaminated soil, can acquire the infection by ingesting the infective eggs. The infection in humans is diagnosed by serology or by identifying the larvae in histological sections.


### Filaria



15-81




15-82



15-83

**15-81-15-83 *Loa loa*. Microfilariae. Blood, Giemsa stain (15-81,  $\times 600$ ; 15-82 and 15-83,  $\times 1250$ ).** The microfilariae of *L. loa* measure approximately  $240 \times 7 \mu\text{m}$ . Although the sheath is not visible on Giemsa stain, a clear halo can sometimes be observed surrounding the microorganism. The cephalic region does not have nuclei (15-82). Typically four to six nuclei at the tip of the tail are spaced evenly and extend to the end (15-83). Deerflies of the genus *Chrysops* transmit the microfilariae, which appear in the blood during the day. The adult microorganism moves through the subcutaneous tissues, producing an inflammatory reaction termed *Calabar*, and can enter the eye and migrate through the conjunctiva. The female adult worms measure 50 to 70 mm by 0.5 mm, while the males are usually half that size.

15-84




15-85



15-86



**15-84-15-86 *Wuchereria bancrofti*. Microfilariae. Blood, Giemsa stain (15-84,  $\times 500$ ; 15-85 and 15-86,  $\times 1250$ ).** The adult worm measures 4 to 10 cm in length and lives in the lymphatics, producing fibrosis and obstruction. The microfilariae are transmitted by several kinds of mosquitoes including *Aedes* spp., *Culex* spp., and *Anopheles* spp. The microfilariae are sheathed and measure approximately  $250 \mu\text{m} \times 8 \mu\text{m}$  (Figure 15-84). The cephalic region is round (Figure 15-85), while the tail is pointed (Figure 15-86). Both areas have an end space devoid of nuclei. The microfilariae are found in the lungs during the day while they circulate in the peripheral blood at night.




**15-87 *Onchocerca volvulus*. Tissue section. Hematoxylin and eosin stain ( $\times 30$ ).** The microfilariae of *O. volvulus* are produced by the adult female worms and become distributed in the skin. The female blackflies, *Simulium* spp., ingest these larvae when they bite. Following a developmental period of 1 to 2 weeks in the fly, the larvae are again infective to humans when the fly next bites. The adult worms produce nodules in the subcutaneous tissues or in fascial planes. Onchocercal nodules are usually well demarcated by a thick band of connective tissue. Worms from both sexes are embedded in a chronic inflammatory infiltrate containing numerous blood vessels and giant cells. Microfilariae in the surrounding connective tissues produce an inflammatory infiltrate with plasma cells, eosinophils, and lymphocytes.


## Trematodes



**15-88 *Clonorchis sinensis*. Adult worm (flukes). Carmine stain ( $\times 18$ ).** The adult form of these trematodes lives in the biliary tract of humans. The adult worms measure approximately 10 to 30 mm in length by 2 to 5 mm in width. The microorganism has a ventral sucker and the coiled uterus, ovary, and two branched testes occupy most of the body.




**15-89, 90 *Clonorchis sinensis* eggs. Iodine (15-89,  $\times 1300$ ) and trichrome stains (15-90,  $\times 800$ ).** The eggs of *C. sinensis* are ovoid and measure approximately  $30 \times 15 \mu\text{m}$ . The shell is relatively thick, and there is a well-defined operculum at the narrow end, and a small knob at the opposite side. The eggs release a free-swimming miracidium when hatched.




**15-91, 92 *Fasciola hepatica*. Egg. Feces. Iodine stain ( $\times 500$ ).** The unembryonated eggs of *F. hepatica* are large, elongated, and measure  $150 \times 80 \mu\text{m}$ . The shell is thin with a small operculum that is difficult to see. The operculum can easily be opened by applying pressure to the coverslip (15-92). The adult worms are large, measuring  $30 \times 15 \text{ mm}$ , and live in the liver and bile duct where they produce eggs that are discharged into the feces.



**15-93 *Paragonimus westermani*. Egg. Feces. Iodine stain ( $\times 600$ ).** These eggs are large, usually  $100 \times 50 \mu\text{m}$ , oval, with a thick shell and a well-defined operculum at the broad end. The opposite end is thickened but lacks a knob such as the one present in *Diphyllobothrium latum* eggs. The adult fluke lives in the lungs of humans and measures  $10 \times 5 \text{mm}$  (not shown).



**15-94, 95 *Schistosoma japonicum*. Carmine stain. Adult male (15-94,  $\times 10$ ) and cercaria (15-95,  $\times 225$ ).** The adult schistosomes reside in blood vessels (Figure 15-94). The eggs produced by the adult female after reaching the water release miracidia that infect specific snails, their intermediate hosts. Following development in the snail, the cercariae (Figure 15-95) emerge and penetrate the skins of humans who are in direct contact with snail-infested water.




**15-96 *Schistosoma haematobium*. Egg. Feces. Iodine stain ( $\times 600$ ).** These eggs are large, measure  $150 \times 50 \mu\text{m}$ , and have a transparent shell with a conspicuous terminal spine.

**15-97 *Schistosoma japonicum*. Egg. Feces. Iodine stain ( $\times 700$ ).** The embryonated eggs of *S. japonicum* are ovoid, measure  $90 \times 50 \mu\text{m}$ , and have a thin shell and a small spine.



**15-98, 99** *Schistosoma mansoni*. Egg. Feces. Iodine stain ( $\times 500$ ). The eggs of *S. mansoni* are large,  $150 \times 60 \mu\text{m}$ , elongated with a thin shell and a distinctive lateral spine. At the time that the eggs are passed in the feces, they contain a miracidium that in fresh preparations can be seen moving. In iodine preparations, the miracidium may stain dark (Figure 15-99).

## Cestodes



**15-100, 101** *Diphylobothrium latum*. Egg. Feces. Iodine stain ( $\times 800$ ). The eggs of *D. latum* are ovoid, measure  $60 \times 50 \mu\text{m}$ , and have a relatively thick shell. The operculum is not very distinctive, and the small knob located at the opposite end in many instances cannot be observed. The operculum can rupture as shown in Figure 15-101. At the time they are passed in the feces, the eggs are unembryonated. It takes approximately 2 weeks before a ciliated embryo develops.







**15-102, 103 *Taenia saginata*. Proglottids. Carmine stain (15-102,  $\times 5$ , 15-103,  $\times 16$ ).** The gravid proglottids measure approximately 18 mm in length  $\times$  5 to 7 mm in width with the genital pore located at the lateral margin. The proglottids passed in the feces can be identified by injecting India ink through the lateral genital pore and counting the number of primary lateral branches of the uterus. Proglottids with 13 or fewer branches belong to *T. solium*, while those with more than 15 lateral branches are *T. saginata*.




**15-104 *Taenia solium*. Scolex. Carmine stain ( $\times 10$ ).** The scolex of *T. solium* measures 1 mm in diameter, has four suckers, and a rostellum with two rows of hooks. The adult worm can measure up to 5 meters. The life cycle is similar to that of *T. saginata*, although the eggs are directly infectious to humans.




**15-105 *Taenia solium*. Cysticercus. Carmine stain ( $\times 15$ ).** Infections of humans with *T. solium* occur in areas where pigs are bred. Ingestion of the eggs of *T. solium* may result in cysticercosis. The eggs hatch in the small intestine and release oncospheres, which are carried in the bloodstream to distant tissues, where they mature into cysticerci. Cysts containing cysticerci may be formed in the central nervous system, skin, muscle, and skin. These cysts are usually 1 to 4 cm in diameter and are filled with a clear fluid and a single scolex that is invaginated.




**15-106, 107 *Taenia solium* proglottids. Carmine stain. (15-106,  $\times 5$ ) (15-107,  $\times 15$ ).**  
The gravid proglottids of *T. solium* have fewer than 13 lateral branches on each side of the uterine central core. These proglottids measure around 10 mm in length by 5 mm in width.



**15-108 *Taenia* spp. Eggs. Iodine stain ( $\times 500$ ).**  
The eggs of *T. saginata* and *T. solium* are spherical, measuring approximately  $40\ \mu\text{m}$  in diameter. The shell is thick with a radial striation and can stain quite dark. In fresh eggs, a thin outer membrane can occasionally be observed. Inside the egg there is a six-hooked embryo. It is important to be able to observe these six hooks for final identification. These eggs have the same morphology as those of *Echinococcus* spp. The eggs of *T. saginata* are not infectious to humans, in contrast to those of *T. solium*. The eggs of *T. saginata* are ingested by cattle, and after several months they develop in the muscles into cysticerci, which are infectious for humans.




**15-109 *Hymenolepis nana*. Egg. Feces. Iodine stain ( $\times 1000$ ).**  
The eggs of *H. nana* are spherical to oval and measure  $30$  to  $45\ \mu\text{m}$  in diameter. The shell is thin, and the inner envelope has two polar thickenings that give rise to four to eight filaments located between the embryo and the shell. The embryo, or oncosphere, has six hooks.



**15-110 *Hymenolepis diminuta*. Egg. Feces. Iodine ( $\times 800$ ).** These eggs are oval, large, measuring  $80 \times 70 \mu\text{m}$ , with a thick shell and an inner membrane that surrounds the oncosphere. The oncosphere has six hooks and is surrounded by a membrane that does not have polar filaments.



**15-111 *Dipylidium caninum*. Scolex. Carmine stain ( $\times 80$ ).** The scolex of this microorganism contains four large suckers and a conical, retractile rostellum with several rows of small spines.





**15-112 *Dipylidium caninum*. Proglottids. Carmine stain ( $\times 15$ ).** Dogs are the usual host of this parasite, but humans can acquire the disease by ingestion of infected fleas. The adult tapeworms measure up to 70 cm. The gravid proglottids measure  $25 \times 8 \text{ mm}$ , and are divided into sections, each containing approximately 10 eggs with six hooks. Typically these proglottids have two genital pores, one on each lateral margin. The eggs can be found occasionally in the feces.




**15-113 *Echinococcus granulosus*. Adult worm. Carmine stain ( $\times 50$ ).** The adult worm of the canine tapeworm *E. granulosus* measures 4 to 6 mm in length, has a scolex with four suckers and a rostellum with hooklets, and three to five proglottids, one of them immature, one mature, and one or two gravid. Dogs release eggs in the feces, which are then ingested by cattle, sheep, and other animals. Humans acquire the disease following accidental ingestion of *E. granulosus* eggs.

15-114



15-115


15-116



15-117





**15-114-117 *Echinococcus granulosus*. Liver. Hematoxylin and eosin stain ( $\times 300$ ).** Hydatid cysts are found in humans mainly in the liver and lung. They consist of an outer wall with a laminated, non-nucleated layer, and an inner layer called *germinative membrane* that is nucleated (Figure 15-114). The germinal layer gives rise to brood capsules with scoleces. Degeneration of these structures results in a fluid material known as hydatid sand (Figure 15-115). In this figure, it is possible to see the wall of the hydatid cyst and four brood capsules with protoscolex. The scoleces are invaginated in their own bodies. At higher magnification, it is possible to observe the rostellum with the hooklets (Figure 15-116). Placing the hydatid sand in saline may cause the scoleces to evaginate (Figure 15-117).





**15-118 *Echinococcus granulosus*. Eggs. Iodine stain ( $\times 800$ ).** The eggs of *E. granulosus* are identical to the eggs of the *Taenia* spp. They are spherical with a thick, radially striated shell and measure 30 to 40  $\mu\text{m}$  in diameter.


## Structures Frequently Found in the Feces that Require Differentiation from Parasites




**15-119, 120 Yeast cells. Feces. Trichrome stain ( $\times 1250$ ).** Yeast cells are usually round to oval and may stain different colors ranging from green to blue to red, depending on status of the wall of the microorganism.




**15-121, 122 Plant cells (15-121: wet mount) and pollen grains (15-122: Trichrome stain). Feces ( $\times 1250$ ).** Plant cells and pollen grains can frequently be seen in fecal specimens. The plant cells can widely vary in size and morphology (Figure 15-121). Pollen usually ranges in size from 15 to 60  $\mu\text{m}$  and has a natural brown to green color since it does not pick up the stain. The wall may be smooth or, as in this case, may have a radiated structure (Figure 15-122).




**15-123 Plantar spiral. Feces. Trichrome stain ( $\times 1250$ ).** These structures are part of the spirovascular bundle of vegetables and are frequently found in the feces. They can have different sizes and degrees of coiling, and should not be confused with parasites.



**15-124 Ciliated columnar epithelial cell. Bronchoalveolar lavage. Trichrome stain ( $\times 1250$ ).** Epithelial cells from the respiratory mucosa can be seen relatively frequently in pulmonary specimens. The morphology may vary depending on how well the specimen was preserved.



**15-125 White blood cells. Feces. Trichrome stain ( $\times 300$ ).** Polymorphonuclear leukocytes and mononuclear cells can be observed in stool specimens. These cells can be confused with members of the *Entamoeba* spp. It is important to consider the size, the ratio of the nucleus to the cytoplasm (usually close to 1:1 in the white cells), and the internal structure of the nucleus. Degenerated macrophages that have ingested debris, including red blood cells, may be particularly difficult to differentiate from *E. histolytica*.




**15-126 Charcot-Leyden crystals. Feces. Trichrome stain ( $\times 1250$ ).** These crystals are by-products of eosinophils and are commonly found in patients with parasitic and allergic diseases. With the trichrome stain they usually stain a red color.

## CHAPTER 16 *Virology*

**V**iruses are the most common cause of human disease. Because of the difficulty in growing them in culture, they have only recently begun to be understood. In contrast to other microbes, mature viruses contain only one type of nucleic acid. DNA viruses that infect humans include pox viruses, some hepatitis viruses, wart viruses, and herpes viruses. RNA viruses include poliovirus, rabies virus, influenza virus, mumps virus, and many others. Viruses also have a protein coat, called a capsid, and some viruses may acquire a lipid envelope during egress from an infected cell or from the nucleus to the cytoplasm of a cell.

Viruses are either detected directly in clinical material using sensitive detection tests including immunological assays and nucleic acid detection methods, or they are cultured in cells and recognized by their damage to the cells, called cytopathic effect (CPE). Analogous to infecting viruses are also often used to diagnose viral infections.


### Instrumentation and Equipment



**16-1** Microscopes (Olympus Optical Co., Ltd., Japan). Compound microscopes consist of an eyepiece (ocular) shown on the right and an objective microscope as detailed. The fluorescence microscope is frequently used by the storage laboratory for the examination of clinical specimens and viruses in culture using transmitted-light optics. The inverted microscope illuminates the sample from the top, with the objectives under a slide to observe. This configuration allows the examination of the cell surface from the bottom of the chamber while the cells are covered by the nutrient medium. Inverted microscopes are particularly useful for the examination of tissue culture cells in bottles. Both electron microscopes are used.




**16-2** Tubes with various culture tubes. Different types of cell monolayers already noted in tubes can be prepared from commercial companies or can be prepared in the laboratory. The monolayers grow in the sterile substrate in contact with tissue culture medium when the tube is set in a shaker during incubation. The cells plus medium (10% to 20% acetone) after clinical specimens have been inoculated into the tube. The monolayers are observed under an inverted microscope for the presence of CPE.




**16-3 Roller drum (Bellco Glass, Inc., Vineland, N.J.).** This instrument holds, and continuously rotates, tubes containing cell lines. The drum is placed inside an incubator at 33°–37°C. The holder containing the tubes can be removed from the incubator so that the tubes can be checked for CPE. This apparatus is most frequently used for culturing respiratory specimens for viruses such as influenza and parainfluenza.



**16-4 Shell vial.** A shell vial is a 15 × 45 mm glass vial containing culture medium and a round 12 to 13 mm glass coverslip on which a monolayer of cells has been cultivated. To inoculate a clinical specimen, the medium is discarded, the specimen is added to the vial on top of the monolayer, and the vial is centrifuged. At the completion of the centrifugation, fresh culture medium is added and the vial is capped and incubated. Specimens suspected of containing certain viruses, including cytomegalovirus, herpes simplex virus, varicella zoster virus, and the bacterium *Chlamydia* spp., among others, are processed using this technique since it shortens the time for identification of the pathogen and may increase the recovery.



**16-5 Multiple well plates.** Plastic plates with 24, 48, or 96 wells are used instead of shell vials by some laboratories for the isolation of viruses and *Chlamydia* from clinical specimens. The main advantage of these plates is that they allow for the simultaneous processing of multiple specimens. Problems include the possibility of cross-contamination between specimens and a lower sensitivity of the assay due to, at least in part, the decrease in sample volume.



**16-6, 7 Tissue culture flasks.** Plastic or glass flasks of different sizes can be utilized to grow cell monolayers or cells in suspension. Once the flasks are seeded, they are placed inside incubators.





**16-8 Spinner bottles.** This type of bottle is used for producing large quantities of viruses in cells that are grown in suspension. In the center of the bottle there is a magnetic stir bar. When the bottle is placed on a stirring apparatus, the magnetic bar spins around and maintains the cells in suspension. The CPE can be visualized by taking an aliquot of the suspended cells and observing them under the microscope.




**16-9, 10 Roller bottles.** Roller bottles are also utilized for growing large quantities of viruses in cell monolayers (Figure 16-9). The cells are seeded into the bottle with medium and the bottles are placed in a roller bottle apparatus (Figure 16-10) located inside an incubator that continuously rotates the bottles. Once the cells are attached in a monolayer around the interior walls of the bottle, the viral inoculum is added and the bottle placed back on the rotating module (New Brunswick Scientific, Edison, N. J.). The cells can be harvested when the desired CPE is observed in the cell monolayer.




**16-11 Liquid nitrogen tanks (Union Carbide Cryogenic Equipment, Danbury, Ct.).** Liquid nitrogen containers are used for the long-term preservation of tissue culture cells and viruses. Specimens and isolates stored under liquid nitrogen, with a temperature of  $-195^{\circ}\text{C}$ , can be maintained in a viable state for indefinite periods of time. Although liquid nitrogen tanks are self-contained, (i.e., they do not need electricity to maintain the temperature) they need to be replenished with liquid nitrogen on a regular basis because the nitrogen is continuously being lost due to evaporation.

### Detection and Identification of Viruses





**16-12 Adenoviridae, Adenovirus. Direct fluorescent assay, Nasopharyngeal swab ( $\times 300$ ).** For the direct detection of adenovirus-infected cells, a mouse monoclonal anti-adenovirus antibody is layered onto and incubated with the fresh specimen on a slide. A goat antimouse fluorescein-labeled conjugate with Evans blue counterstain is then added. As shown here, positive specimens have an apple-green stippled staining of the nucleus.



**16-13, 14 Adenoviridae, Adenovirus. MRC-5 cells. Cytopathic effect (CPE). Phase contrast ( $\times 225$ ).** Adenoviruses grow well in a variety of human cell lines such as HeLa, HFp-2, KB, A549, and can also be isolated in primary human kidney cells and in fetal diploid fibroblast, such as MRC-5, cell cultures. Typically the initial CPE appears at 3 to 5 days and consists of rounding of individual cells (Figure 16-13) that subsequently become significantly enlarged and refractile, giving the appearance of balloons. As the CPE continues and involves the surrounding cells, it forms grape-like clusters (Figure 16-14). The CPE usually spreads to the rest of the monolayer over a period of days.




**16-15 Adenoviridae. Adenovirus. MRC-5 cells. Culture identification. Fluorescent stain ( $\times 500$ ).** To confirm the identification of an adenovirus, fluorescein-labeled monoclonal antibodies can be used. Here fluorescence is visible both in the cytoplasm and the nucleus of the infected cells. Adenoviruses replicate in the nucleus where bright fluorescent dots can be observed.




**18-16, 17 Adenoviridae. Adenovirus. Electron microscopy. Transmission (18-16,  $\times 27,500$ ) and negative staining (16-17,  $\times 100,000$ ).** Tissues or cells suspected of containing adenoviruses can be fixed, cut, and stained for observation by transmission electron microscopy. Adenoviruses can form crystalline arrays in the nucleus of the cell (Figure 16-16). Urine and feces specimens can be clarified and observed using the negative-staining technique, in which the background is stained dark to reveal the unstained virus particles (Figure 16-17). Adenoviruses measure 60 to 90 nm in diameter, do not have an envelope, are of cubic symmetry, and possess a capsid composed of 252 capsomers, each approximately 7 to 9 nm in diameter.




**16-18 Adenoviridae. Adenovirus. Antibody detection. Fluorescent assay ( $\times 800$ ).** Several techniques can be used for the detection of antibodies to adenoviruses including complement fixation tests, hemagglutination inhibition, enzyme immunoassays, and fluorescent methods. This figure shows the detection of antibodies to adenoviruses using an indirect immunofluorescent assay. Human epithelial cells are infected with adenovirus, mixed with non-infected cells, and then the cells are fixed to a glass slide. The serum from the patient is incubated with the specimen and washed; an antihuman fluorescein-labeled globulin is then added, and the specimen is observed using a fluorescence microscope. If antibody is present, as in this case, infected cells fluoresce apple-green.



**16-19 Calciviridae. Calcivirus. Feces. Electron microscopy. Negative staining ( $\times 400,000$ ).** Calciviruses have been identified as a frequent cause of gastroenteritis in young children. Diagnosis requires electron microscopy or immunoserological testing. These virions are spherical, 35 to 40 nm in diameter, with 32 cup-shaped depressions and icosahedral symmetry.



**16-21 Herpesviridae. Herpes simplex virus (HSV). Skin lesions.** HSV-1, in general, produces infections in the upper part of the body, while HSV-2 is frequently isolated from sites below the waist, including the genital area. Skin and mucocutaneous lesions, in the form of fluid-filled vesicles approximately 2 to 5 mm in diameter, evolve over a period of 7 to 14 days. At the time that vesicles are present, direct antigen detection from scrapings of tissue at the vesicle base using monoclonal antibodies is a highly specific and sensitive method.




**16-20 Coronaviridae. Coronavirus. Feces. Electron microscopy. Negative staining ( $\times 100,000$ ).** Coronaviruses have been associated with acute upper respiratory tract infections in children and diarrheal illnesses in children and adults. Identification of these viruses requires electron microscopy studies or culture in human embryonic trachea tissue cultures. The coronaviruses are enveloped, pleomorphic particles that measure 80 to 120 nm in diameter. This family of viruses is named for the club-shaped peptomers (virally encoded proteins) that project from the envelope and that, as shown here, form a thin border resembling a "solar corona" under electron microscopic examination of negatively stained preparations.




**16-22 Herpesviridae. Herpes simplex virus (HSV). Direct antigen test kit (Syva, Co., San Jose, Calif.).** A system frequently used for the detection of certain pathogens, such as HSV and *Chlamydia trachomatis*, is shown in this slide. Swabs are provided for collecting specimens from the eye, the female cervix and the male urethra, or other body sites. Once the specimen is collected, the swab is rolled onto the glass slide and fixative is added. The slide, enclosed in the provided container, is then submitted to the laboratory for staining with a fluorescein-labeled monoclonal antibody.


**16-23 Herpesviridae. Herpes simplex virus (HSV). Direct immunofluorescent assay (DFA) ( $\times 500$ ).** Direct staining with a specific monoclonal antibody of infected cells collected from the base of HSV vesicle allows for the rapid detection and identification of this virus. Although this technique may be overall less sensitive than culture, it has the advantage that a viable virus is not required for obtaining a positive result. An additional advantage is that the test results are available shortly after collecting the specimen. As shown in this DFA-stained slide, infected cells show an apple-green fluorescence in both cytoplasm and nucleus.





**16-24 Herpesviridae. Herpes simplex virus. Cervical swab. Horseradish peroxidase stain ( $\times 300$ ).** Cervical specimens, collected as for Pap smears, can be stained with specific antibodies to HSV. As shown here, infected cells stained with horseradish peroxidase have a brown color.



**16-25, 26 Herpesviridae. Herpes simplex virus (HSV). MRC-5 cells. CPE seen under phase contrast ( $\times 225$ ).** HSV 1 and 2 grow rapidly *in vitro*, producing CPE in 24 to 48 hours. Several cell lines can be used to isolate and identify these viruses, including human diploid fibroblasts such as MRC-5 or WI-38, and primary rabbit kidney cells. As shown in Figure 16-25, early CPE is characterized by large, round cells appearing as localized foci in several arcs of the monolayer. The CPE progresses rapidly and by 3 to 5 days usually involves the entire monolayer. Occasionally, formation of multinucleated syncytial giant cells can be observed (Figure 16-26).



**16-27, 28 Herpesviridae. Herpes simplex virus (HSV). MRC-5 cells. Shell vial culture. Fluorescein (16-27,  $\times 300$ ) and horseradish peroxidase stains (16-28,  $\times 225$ ).** The monolayers were stained after 48 hours of incubation. After fixation, a fluorescein-labeled monoclonal antibody (Figure 16-27) or a horseradish peroxidase-labeled antibody (Figure 16-28) was used to detect the infected cells. Specific monoclonal antibodies to HSV 1 and 2 can be used separately or in combination, thus allowing for typing of the isolate.




**16-29, 30 Herpesviridae. Herpes simplex virus (HSV). Transmission electron microscopy (16-29,  $\times 100,000$ ; 16-30,  $\times 35,500$ ).** The family of Herpesviridae has several members that are frequent human pathogens, including HSV 1 and 2, cytomegalovirus, varicella zoster virus, and Epstein Barr virus, in addition to the recently discovered HSV 6, HSV 7, and HSV 8. All have a similar morphological structure consisting of a cylindrical core containing the viral DNA, an icosahedral capsid that measures 90 to 110 nm in diameter, a granular zone that surrounds the capsid, and an external envelope. The complete viral particle measures 180 to 200 nm in diameter. As shown in Figure 16-29, the viral particles form in the host cell nucleus, but the envelope is acquired at the time of budding through the cellular membrane (Figure 16-30).





**16-31 Herpesviridae. Herpes simplex virus (HSV). Antibody detection. Fluorescent stain ( $\times 225$ ).** Human diploid cells infected with HSV 1 or 2 can be used as substrate to detect antibodies to these viruses in human serum. There is a significant amount of cross reactivity between antibodies to HSV 1 and HSV 2 and thus, these tests should not be used to identify the specific virus causing the infection. Western blot analysis and detection of antibodies to specific HSV 1 and HSV 2 antigens are now available for that purpose.




**16-32 Herpesviridae. Cytomegalovirus (CMV). MRC-5 cells. CPE seen under phase contrast ( $\times 225$ ).** CMV has species-specific growth requirements; thus, human diploid fibroblast cells, either MRC-5 or WI-38, are used for culture. The foci of CPE are slow to appear but usually are visible by 5 to 10 days, so cultures should be maintained for up to 21 days. CPE is characterized by the presence of round, large, refractile cells in elongated foci parallel to the long axis of the monolayer. The CPE spreads slowly and, as a result, the surrounding cells maintain a normal morphology for extended periods of time. In most instances, unless the initial inoculum is large, the entire monolayer is not involved.





**16-33, 34 Herpesviridae. Cytomegalovirus (CMV). MRC-5 cells. Shell vial. Immunofluorescent stain ( $\times 500$ ).** The use of the shell vial centrifugation method has allowed early detection of CMV in clinical specimens. Following centrifugation of the specimen onto the monolayer, the culture is incubated for 24 to 48 hours at 37°C and subsequently stained with a monoclonal antibody to one of the early CMV antigens. As shown in Figure 16-33, staining of the nuclei occurs early, before cytopathic effect can be detected. Monolayers stained at a later time after infection show both nuclear and cytoplasmic fluorescence. Early CPE can be observed in Figure 16-34.



**16-35, 36 Herpesviridae. Cytomegalovirus (CMV).** 16-35; Kidney. Hematoxylin and eosin (H&E) stain ( $\times 500$ ). 16-36; Lung. Horseradish peroxidase (HRP) stain ( $\times 500$ ). Tissues infected with CMV may have nuclear and cytoplasmic inclusions. The most characteristic nuclear inclusions have the appearance of an "owl's eye" as a result of the retraction of the tissue, producing the clear halo. Granular eosinophilic inclusions in the cytoplasm are not so distinctive. These inclusions are visible in H&E stained preparations (Figure 16-35) but the use of specific antibodies labeled with HRP allow for a more specific and sensitive identification (Figure 16-36).




**16-37 Herpesviridae. Cytomegalovirus (CMV). Antibody detection. Immunofluorescent assay ( $\times 600$ ).** Monolayers of human fibroblasts are infected with CMV for use as the substrate for this assay. It is important to have a mixture of infected and noninfected cells in the monolayer in order to be able to discriminate between specific and nonspecific staining. The specific CMV staining should be mainly nuclear, as shown here. Seroconversion from a negative to a positive status is a good indication of a primary infection. However, changes in antibody titers should only be interpreted in conjunction with other clinical and laboratory parameters.



**16-38, 39 Herpesviridae. Varicella zoster virus (VZV). Skin lesions. H&E stain** 16-38,  $\times 30$  and 16-39,  $\times 500$ . Intraepidermal vesicle resulting from a VZV infection (Figure 16-38). Several multinucleated giant cells with eosinophilic Cowdry type A (typical intranuclear) inclusions are visible (Figure 16-39).






**16-40 Herpesviridae. Varicella zoster virus (VZV). Skin lesion. Direct fluorescent assay ( $\times 300$ ).** Specific monoclonal antibodies can be used for the identification of VZV infected cells in clinical samples. As shown here, fluorescein-tagged antibodies produce apple-green fluorescence in the nucleus and cytoplasm of the infected cells.



**16-42 Herpesviridae. Varicella zoster virus (VZV). Shell vial. Fluorescent stain ( $\times 500$ ).** Human diploid lung cells are used to prepare shell vials for the detection of VZV. Commercially available fluorescein-conjugated monoclonal antibodies are used to stain the monolayer 48 to 72 hours after inoculation. This method greatly facilitates the identification of VZV and speeds up the process by several days compared with waiting for CPE. Nuclear and cytoplasmic fluorescence can be observed in this slide.

**16-44 Herpesviridae. Epstein-Barr virus (EBV). Commercial particle agglutination test for serological diagnosis.** EBV does not readily grow in tissue culture; therefore, the diagnosis of an EBV infection is frequently made using serological tests. A screening test for heterophile antibodies. IgM antibodies that react with antigens that are not related to the organism producing the antibody response, is positive in most patients with infectious mononucleosis. A latex agglutination test, "Monospot-type" test (Figure 16-44, Biokit USA, Inc., Lexington, Mass.), is often used to detect heterophile antibodies.




**16-41 Herpesviridae. Varicella zoster virus (VZV). A549 cells. CPE seen under phase contrast ( $\times 225$ ).** Human diploid fibroblast cell cultures are used in most laboratories for the isolation and identification of VZV. The CPE produced is slow to appear and does not spread readily to the rest of the monolayer. The cultures should be maintained for 21 days before they are discarded. Figure 16-41 shows VZV CPE in A549 cells, in which the initial CPE consists of swollen, refractile cells. As the CPE progresses, it acquires a doughnutlike shape containing a center of necrotic cells surrounded by large, refractile, giant cells. Spread of the foci occurs by infection of adjacent cells, since the virus is highly cell-associated.





**16-43 Herpesviridae. Varicella zoster virus. Antibody detection. Immunofluorescent assay ( $\times 300$ ).** Cells infected with VZV can be used as the antigen in an indirect immunofluorescent test for the presence of antibodies to this virus.








**16-45, 46 Immunofluorescent tests for antibodies to Epstein-Barr virus structural antigens.** Detection of antibodies to specific antigens of EBV can be accomplished using cell lines infected with EBV as substrates for indirect immunofluorescent methods. P3HR-1 cells expressing the EBV viral capsid antigen (VCA) are used for detecting IgM (Figure 16-45) and IgG (Figure 16-46) antibodies to VCA.



**16-47, 48 Immunofluorescent tests for antibodies to Epstein-Barr virus early antigens.** Antibodies to two types of early antigens (EA) are seen in infected Raji cells: diffuse (D), in which the antigen is distributed in the nucleus and cytoplasm (Figure 16-47); and restricted (R), in which the antigen is only in the cytoplasm. Figure 16-48 shows antibodies against both D and R early antigens.



**16-49 Immunofluorescent test for antibodies against nuclear antigen of Epstein-Barr virus (EBNA).** The Raji cell line is also used for the detection of antibodies to EBNA using an anti-complement indirect immunofluorescent staining technique.



**16-50 Orthomyxoviridae. Influenza A. Direct fluorescent assay ( $\times 1350$ ).** Nasopharyngeal specimens can be directly stained using fluorescein-labeled monoclonal antibodies to influenza A, B, and C viruses. An apple-green fluorescent granular pattern can be observed in the nucleus of this infected cell.





**16-51, 52 Orthomyxoviridae. Influenza A. African green monkey kidney cells. Hemadsorption. Phase contrast ( $\times 225$ ).** Although in the past, embryonated hens' eggs were used for the isolation of influenza viruses, continuous cells lines are currently utilized in most laboratories. Certain viruses, including influenza, parainfluenza, and mumps viruses, produce glycoproteins that are incorporated into the host cell membrane. Although these viruses cause minimal or no cytopathic effect, guinea pig red blood cells adhere to the host cell membranes when they are added to tissue culture cells infected with these viruses. This phenomenon, termed *hemadsorption*, is demonstrated in Figure 16-51. Figure 16-52 shows a control uninfected monolayer of African green monkey kidney cells.

**16-53 Orthomyxoviridae. Influenza. Shell vial. Fluorescent stain ( $\times 300$ ).** Inoculation of the specimen into a shell vial followed by staining at 48 to 72 hours with a specific fluorescein-labeled monoclonal antibody allows for the rapid detection and identification of influenza viruses from clinical specimens. Fluorescent staining appears both in the nucleus and cytoplasm of the infected cells.




**16-54 Orthomyxoviridae. Influenza. Antibody detection. Immunofluorescent assay, ( $\times 500$ ).** Complement fixation, enzyme immunoassays, hemagglutination inhibition, and fluorescent methods are used for the detection of antibodies to influenza viruses. In this slide, a monolayer of cells infected with influenza has been used as a substrate to test for the presence of anti-influenza antibodies using an indirect immunofluorescent method.







**16-55, 56 Orthomyxoviridae. Parainfluenza 3 (16-55) and 4 (16-56). Nasal washings. Immunofluorescent staining ( $\times 500$ ).** Fluorescein-labeled monoclonal antibodies are used for direct staining of clinical specimens for the detection of parainfluenza types 1, 2, 3, and 4 viruses. The fluorescent staining appears predominantly in the cytoplasm of the cells.




**16-57 Papovaviridae. JC virus. Progressive multifocal leukoencephalopathy (PML). Brain tissue section. Immunoperoxidase stain ( $\times 500$ ).** Fluorescent and enzymatic immunostaining methods and nucleic acid probes can be used for the detection of the JC virus in biopsies from patients suspected of having PML. In this case, the infected glial cells stain dark brown.



**16-58, 59 Paramyxoviridae. Measles. Phytohemagglutinin stimulated cord blood mononuclear cell culture. CPE; 16-58 ( $\times 500$ ), Fluorescent staining; 16-59 ( $\times 200$ ).** Primary monkey kidney and Vero cells are frequently used for the isolation of the measles virus. The resulting CPE may consist of spindle-shaped single cells or syncytial multinucleated giant cells. Peripheral blood mononuclear cells have also been used for detection of the measles virus where, as shown in Figure 16-58, formation of giant cells can be observed. The virus responsible for the CPE can be confirmed by staining the culture with fluorescein-labeled antibodies (Figure 16-59).



**16-60 Paramyxoviridae. Measles. Vero cell culture. Fluorescent staining ( $\times 300$ ).** The measles virus-infected cells are stained with a fluorescent-labeled specific antibody.



**16-61 Paramyxoviridae. Measles. Lung tissue section. H&E ( $\times 500$ ).** Measles virus infections in the lungs can produce a severe pneumonia. The histological section shown here depicts the formation of multinucleated giant cells (Warthin-Finkeldey cells).

**16-62 Paramyxoviridae. Measles. Antibody detection. Indirect immunofluorescent assay ( $\times 500$ ).** Antibodies to the measles virus can be detected using different techniques, including complement fixation, indirect fluorescent antibody stains, hemagglutinin inhibition, and enzyme immunoassays. In this slide, a measles-infected cell culture has been incubated with patient serum and stained with a fluorescein-conjugated antihuman immunoglobulin. Either IgG or IgM can be detected in this way.



**16-63 Paramyxoviridae. Respiratory syncytial virus (RSV). Direct fluorescent assay ( $\times 500$ ).** Fluorescein-labeled monoclonal antibodies can be used for the direct detection of RSV in clinical specimens. Typically an apple-green stippled fluorescence can be observed in the nucleus and cytoplasm of the infected cells, as shown in this slide. RSV is sensitive to temperature and dry conditions and as a result, it quickly loses its viability unless it is collected and transported to the laboratory under optimal conditions. Thus, in certain clinical situations, the direct assays may be more sensitive than culture.






**16-64, 65 Paramyxoviridae. Respiratory syncytial virus. CPE in MRC-5 cells. Phase contrast ( $\times 225$ ).** HEp-2, HeLa, and MRC-5 cell lines are frequently used for the isolation of RSV. CPE usually appears by 5 to 10 days and is characterized by the formation of multinucleated giant syncytial cells (Figure 16-64). Occasionally, on heavily positive specimens, CPE may appear by 2 to 4 days; in this case the monolayer may show only rounded up cells that can be confused with a toxic or degenerative effect (Figure 16-65).




**16-66 Paramyxoviridae. Respiratory syncytial virus. MRC-5 cells in shell vial. Culture identification. Fluorescent stain ( $\times 300$ ).** RSV can be identified in cell culture using fluorescent-labeled monoclonal antibodies. As shown in this slide, an apple-green bright speckled stain can be observed in the cytoplasm of the infected cells.




**16-67, 68 Picornaviridae. Echovirus 11. CPE in MRC-5 cells. Phase contrast ( $\times 225$ ).** The family Picornaviridae includes four genera: enterovirus, rhinovirus, cardiovirus, and aphthovirus. The first two genera are significant human pathogens. The enteroviruses comprise the coxsackieviruses, echoviruses, and polioviruses. Most of the enteroviruses replicate well in primary monkey kidney cells and produce CPE in 2 to 5 days. As shown in Figure 16-67, the infected cells round up, shrink, and become refractile. The cells quickly degenerate, with marked pyknosis (shrinking) of the nuclei, and detach from the surface of the container (Figure 16-68). Identification of the specific enterovirus isolate is performed by blocking the infectivity with pools of neutralizing sera.




**16-69 Picornaviridae. Rhinovirus. CPE in MARC-5 cells. Phase contrast ( $\times 225$ ).** Human diploid fibroblasts are often utilized for the isolation of rhinoviruses from clinical specimens. It is recommended that these cultures be incubated in roller drums inside incubators at 33 to 34°C. The CPE is characterized by the formation of rounded cells of uneven size and a refractile or ground glass appearance. These foci of CPE usually become evident during the first week and subsequently spread throughout the rest of the monolayer.



**16-70 Poxviridae. Vaccinia virus. Tissue culture preparation. Transmission electron microscopy ( $\times 60,000$ ).** Poxviruses are large and very complex viruses with a typical brick-shape morphology measuring 250  $\times$  200 nm. The internal biconcave core contains the DNA genome. Two lateral bodies are arranged along the concavities of the core, and the virion is surrounded by a viral membrane. Inside the host cell, the virion may have a double membrane that is lost when the virus is extruded from the cell.




**16-71 Reoviridae. Rotavirus. Feces. Electron microscopy. Negative staining ( $\times 100,000$ ).** Rotaviruses are one of the most common agents causing acute gastroenteritis in children. Visualization by electron microscopy was the standard method for detecting these viruses in fecal specimens until sensitive enzyme immunoassays were developed; now EIA is the most frequently used method for diagnosis. Rotavirus particles measure approximately 70 nm in diameter and have the appearance of a wheel, thus the name *Rota* from the Latin word for *wheel*. The core contains 11 segments of double-stranded RNA and is surrounded by a capsid that has 92 capsomers.



**16-72 Retroviridae. Human immunodeficiency virus type 1 (HIV-1). Tissue culture. Transmission electron microscopy ( $\times 100,000$ ).** HIV has a cylindrical core surrounded by the viral envelope. In this slide, transverse sections of the core appear round and electron dense. The virions measure 80 to 100 nm in diameter and have glycoprotein projections, as shown here, that measure 8 to 10 nm in diameter. The virions bud from the cytoplasmic membrane of the host cell. (Special thanks to Dr. Edward Robinson for growing this preparation of HIV-1).



**16-73 Retroviridae. HIV-1. Antibody detection. Enzyme immunoassay.** Screening for the presence of antibodies to HIV-1 is often performed by EIA. The antigen used in these tests can be purified HIV-1 particles or HIV-1 recombinant proteins that have been expressed in *E. coli* or in another type of vector. The antigen is attached to a solid phase (such as the plastic bead shown here) and the patient samples are added. Following incubation and washing, an enzyme-labeled antihuman antibody is added. The sample is again incubated and washed. Addition of the substrate results in the formation of a color reaction (first tube on the left) that can be measured with a spectrophotometer. A positive EIA test should be confirmed by Western blot or an immunofluorescent assay.




**16-74, 75 Retroviridae. HIV-1. Antibody detection. Western blot test.** Immunoblotting is similar to enzyme immunoassays (EIA). The microorganism of interest, in this case HIV 1, is disrupted using a detergent and heat, and the various viral components (antigens) are separated in a gel by electrophoresis. The antigens on the gel are transferred (blotted) to a nitrocellulose membrane using a special instrument (Figure 16-74, BioRad, Hercules, Calif). Serum from the patient is incubated with strips cut from the membrane. After binding has occurred, enzyme-labeled antihuman globulin and the substrate for the enzyme are added. A colored enzymatic end product is produced in the regions (bands) of the strip where the antibodies from the sample were bound to the antigenic components of the pathogen (Figure 16-75).




**16-76 Retroviridae. HIV-1. Antibody detection. Immunofluorescent assay ( $\times 300$ ).** The substrate for this test consists of a mixture of HIV-1 infected and uninfected H9 cells (a T cell line), which are used to prepare a smear on a microscope slide. The serum to be tested is added and incubated. This is followed by a second incubation with an antihuman immunoglobulin conjugated to fluorescein isothiocyanate. Positive IFA reactions are defined as diffuse cytoplasmic fluorescence, but in certain instances, positive staining patterns may be focal or limited to membrane staining. Nonspecific staining can be adsorbed out using noninfected H9 cells.









**17-3 Direct EIA.** The basic principle of EIA tests uses an antigen or an antibody bound to a solid phase (often a plastic well or a plastic bead) to which the patient sample is added. If the appropriate homologous component is present, antigen-antibody binding occurs. In the direct EIA, a labeled antibody (called a conjugate) against the bound substance is subsequently added to the reaction. The label is usually an enzyme that catalyzes a reaction yielding a colored end product. Addition of the enzyme's substrate results in color development.



**17-4 Indirect EIA.** For this test, an unlabeled specific antibody is added to the specimen before the second labeled conjugate is added. Thus a single conjugate can be used to detect antibodies of the same Fc type directed against numerous different specific antigens. Figure 17-4 shows a positive reaction in wells on the top row (yellow end product of horseradish peroxidase enzyme and orthophenylenediamine substrate) and negative reactions in the wells on the bottom row. Other conjugates utilizing chemiluminescent, fluorogenic, or radioactive substrates are also in use.



**17-5 Complement fixation method.** In the past, antibodies to certain infectious disease agents, primarily viruses, could be detected only by complement fixation. In this system, an antibody directed against red blood cells (RBC) is able to lyse the RBC only in the presence of complement. The antigen against which the antibodies to be detected have been produced is introduced into the system, along with serum being tested. If the serum contains antibodies, they will bind to the specific antigen, and the subsequent antigen-antibody complexes will bind, or "fix" complement. Without free complement, the anti-RBC antibody cannot lyse the RBCs, and they form a button at the bottom of a well in a microtiter plate or a tube.




**17-6, 7 Direct and indirect immunofluorescence (IF) methods.** IF methods can be utilized for the detection of specific microbial antigens and antibodies. Direct (Figure 17-6) or indirect (Figure 17-7) approaches can be used.



**17-8 Indirect fluorescent antibody method for *Toxoplasma gondii* antibodies.** *T. gondii* cysts on the surface of slides act as the antigen. The serum samples to be tested are layered over the cysts, and the slides are incubated to allow binding of specific anti-*Toxoplasma* antibody present in the serum. Subsequently antihuman immunoglobulin conjugated with a fluorescent dye, fluorescein isothiocyanate, is added, and unbound conjugate is washed off. In a positive assay result (Figure 17-8), the microorganisms fluoresce.



**17-9 Microimmunofluorescence (MIF) test for rickettsial antibodies.** Suspensions of six rickettsial antigens: *Rickettsia rickettsii*, *R. akari*, *R. typhi*, *R. prowazekii*, and *Coxiella burnetii* phase I and II microorganisms are placed as dots on a multiwell glass slide. The top row of wells has been marked with a black pen to demonstrate how several antigen dots can be placed into one well for the MIF test.




**17-10 Positive MIF test for rickettsial antibodies.** After patient serum and the fluorescein labeled immunoglobulin are added, a positive reaction reveals strongly fluorescing microorganisms of the specific species against which the patient has produced antibody, and may show weakly fluorescent results for the other species.

**17-11 Indirect immunofluorescent test for antibodies to *Borrelia burgdorferi*.** *B. burgdorferi*, the spirochetal agent of Lyme disease, is used as substrate. Both IgG and IgM antibodies can be detected.


**17-12 *Coccidioides immitis*. Antibody detection. Passive immunoprecipitation test.** The *C. immitis* antigen is placed in the central well of the gel, and patient and control sera are placed in the six peripheral wells. Immunoprecipitation of the antigen-antibody complex occurs following passive diffusion of the antigens and antibodies toward each other within the gel. To ensure specificity, the positive sample must produce a line of identity with the positive control (e.g., the precipitation line formed by the patient's sample merges with that of the positive control).




**17-13 Indirect hemagglutination (IHA) test for detection of antibodies to *Entamoeba histolytica*.** In this test, *E. histolytica* antigen is adsorbed onto the surface of tanned human erythrocytes. Sera having specific IgG or IgM antibodies to this microorganism will agglutinate the RBC, which will result in the appearance of a mat in the bottom of the U-shaped wells in a microtiter plate. If the sample does not have antibodies to *E. histolytica*, the RBCs will form a compact cell button (negative control, row 2). In general, patients with extra-intestinal amebiasis develop antibodies 1 to 2 weeks into their symptomatic phase. Positive and negative serum controls and nonsensitized cells should be used in the test to ensure the specificity and sensitivity of the reaction.




**17-14 Microhemagglutination assay for the detection of antibodies to *Treponema pallidum* (MHA-TP).** This test is used to detect the presence of specific antibodies to *T. pallidum*, present in patients with syphilis. The antigen for the test consists of formalinized tanned sheep RBC sensitized with *T. pallidum* (Nichols strain). The sample is first adsorbed with nonpathogenic Reiter treponemes to remove cross-reactive antibodies. Sensitized sheep RBC are then added. Samples containing antibodies to *T. pallidum* will react with the RBCs, which will form a smooth mat of agglutinated cells at the bottom of the plate (wells 1 through 5, row 1). On the other hand, if there are no antibodies to *T. pallidum*, the unagglutinated RBCs will form a button (row 3).






### 17-15 Detection of antibodies to *Coccidioides immitis* by the complement fixation (CF) test.

The CF test is based on the ability of complement to lyse RBC in the presence of RBC-specific antibody. Dilutions of the patient sample, *C. immitis* antigen, and guinea pig complement are added to a microtiter plate and allowed to react before a standard amount of sensitized sheep RBC and goat anti-sheep RBC antibody, or hemolysin, is added to the wells. Hemolysis of the sheep RBCs occurs only in those wells that contain free complement, indicating the absence of a specific antigen-antibody reaction (row 3, for example). No hemolysis occurs in the wells where complement fixing antigen-antibody complexes utilized the complement. Thus, samples with antibody will not lyse the RBCs, indicating a positive test (wells 1 through 3, row 1).



### 17-16 Rapid plasma reagin test (RPR).

This test detects the presence of "reagin" antibodies in sera of patients infected with *Treponema pallidum*. It is thought that "reagin" is an antibody developed against tissue lipids. *T. pallidum* damages the tissue of the host causing lipoidal fractions to combine with spirochetal proteins, which subsequently stimulate the production of antibodies. Reagin antibodies are not specific for *T. pallidum*, but they can be used in a sensitive screening test, such as RPR. When the "reagin" binds to the RPR antigen, macroscopic flocculation occurs. The antigen consists of cardiolipin, lecithin, and cholesterol bound to charcoal particles. Choline chloride is added to block complement so that it is not necessary to heat-inactivate the sample. Reactive sera show a clumping of the charcoal particles (top row).




### 17-17, 18 Latex particle agglutination test for the detection of *Haemophilus influenzae*, *Streptococcus pneumoniae*, *Neisseria meningitidis*, and group B streptococci (Becton Dickinson Microbiology Systems, Cockeysville, Md.).

Agglutination tests using polystyrene beads coated with specific antibody or antigen are easy and rapid to perform. Antigen-antibody complexes present in a positive test show clumping.

## CHAPTER 18 Molecular Techniques

**B**ecause clinical specimens often contain very low numbers of an infecting pathogenic microbe, amplification methods can be used to increase the sensitivity of the test and the likelihood of detection. Such methods remove the need to grow the organism *in vitro*, a slow and variable process. New methods to generate multiple copies of a specific section of nucleic acid of the microorganism being sought allow not only detection of the organism, but detection of specific genetic characteristics, such as a resistance gene. Problems with these techniques include contamination from previously amplified product, inhibition of amplification by specimen components, inability to differentiate viable from nonviable microbes, and interpretation of results. Although the final utility of these methods in clinical microbiology is not yet well-defined, there is no doubt that such techniques will be important elements of the future microbiology laboratory.

**18-1 DNA hybridization.** The DNA hybridization technique is a useful molecular method for the detection and identification of pathogens. The DNA in the specimen is extracted (e.g., the two strands of the DNA molecule are separated using heat or alkali treatment). A labeled single-stranded piece of DNA complementary to a DNA sequence unique to the pathogen being sought (the “probe”) is added and hybridization (i.e., one of the DNA strands of the specimen if it has a complementary sequence to the probe) is present. The hybridization is detected by various methods, including enzymatic and fluorescent techniques.




**18-2 Radioisotoped probe detection of herpes virus.** An epidemiologist needs to determine how long herpes simplex virus (HSV) or VZV. It also has been hybridized with  $^{32}\text{P}$  labeled specific HSV-1 or HSV-2 probes (left and right panels) or without any probes simultaneously (right panel).






**18-3 DNA hybridization of human papilloma virus probes.** Specimen material treated with radioactive isotope labeled probes for the detection of human papillomaviruses types 6, 11, 16, 31, 33, and 35.




**18-4 Gel apparatus (International Biotechnology, Inc., New Haven, Conn.).** Restriction fragment length polymorphism (RFLP) analysis allows for the detection of genetic differences between organisms. Restriction endonucleases are enzymes that recognize, bind and cleave specific DNA sequences. When the DNA of different organisms is digested with these restriction enzymes, fragments of different sizes may result. After enzyme treatment, DNA fragments are separated on a gel apparatus.




**18-5 Herpes simplex virus DNA restriction patterns.** DNA from five different herpes simplex virus isolates was cleaved with three different restriction endonucleases (five rows per endonuclease), electrophoresed on an agarose gel, and stained with ethidium bromide.





**18-6 Southern blot.** This technique allows for the detection of specific DNA fragments in complex DNA samples. To perform this test, the DNA is first cut with restriction endonucleases. The different size DNA fragments are then separated in a gel based on their electrophoretic mobility (size and charge). The DNA fragments are subsequently transferred, or blotted, onto a nylon membrane or to a similar flat support matrix. The membrane is then hybridized with a specific labeled probe containing the DNA sequence homologous to the one sought. If that sequence is present in the sample DNA, the labeled probe will hybridize to it, and the specific labeled band will be detected in the membrane.



**18-7 DNA sequencing apparatus (International Biotechnology, Inc., New Haven, Conn.).** Four DNA oligonucleotides (short sequences) complementary to one of the strands are synthesized using specific DNA primers. In each of the reactions, a different dideoxynucleoside triphosphate (ddNTP) is randomly incorporated. At the point where a ddNTP, rather than a dNTP, is added, DNA polymerase cannot continue the synthesis, and the strand is truncated. By labeling and separating these fragments based on their length, the DNA sequence of a 300 to 500 base fragment can be determined.



**18-8 Labeled DNA fragments from a sequencing reaction.** This technique allows for the identification of specific microbial isolates and for the determination of the genetic relatedness between different microorganisms.



**18-9 Nucleic acid amplification thermocycling instrument (Coy Laboratory Products, Inc., Ann Arbor, Mich.).** The new molecular approaches use biochemical techniques to amplify *in vitro* the nucleic acids of microorganisms of interest in patient specimens. *In vitro* amplification of DNA is mediated by DNA polymerase and amplification of RNA is accomplished after conversion by a reverse transcriptase enzyme to DNA. Several methods for the amplification of nucleic acids require repetitive changes in the temperature of the sample, held in a programmable thermocycling instrument.




**18-10 Positive displacement pipet used in nucleic acid amplification tests (Gilson Medical Electronics, Villiers-le-Bel, France).** Amplification methods can result in the production of billions of copies of the same molecule, all concentrated in a very small sample volume. As a result, it is very easy to cross-contaminate specimens during handling of the samples. In order to avoid cross-contamination during pipeting, positive displacement pipets or pipet tips with special filters are used.




**18-12 PCR reaction products on agarose gel.** The amplified DNA strands in the PCR reaction mixture are separated electrophoretically and stained with ethidium bromide for visualization.




**18-13 Automated PCR instrument (Perkin Elmer, Norwalk, Conn.).** The system includes a thermocycler and an EIA reader to perform and read the result of a PCR that uses an enzyme-labeled detection system.




**18-11 Polymerase chain reaction (PCR) basic principle.** This amplification reaction requires two oligonucleotides, 20 to 30 bases in length, complementary to each one of the two strands of a DNA molecule, and separated by approximately 200 to 400 bases in the target molecule. The DNA in the specimen is first denatured by heat, and the oligonucleotides are added to hybridize to their complementary strands. Addition of a DNA polymerase (Taq) that copies the two original strands using the oligonucleotides as primers completes the cycle. The cycle is repeated 30 to 40 times until enough copies of the target DNA can be detected.



**18-14 Ligase chain reaction (LCR) or oligonucleotide ligation amplification (OLA).** This reaction relies on the enzymatic activity of a DNA ligase to amplify the target nucleic acid. In the LCR, two oligonucleotides homologous to adjacent sequences on the target DNA are joined together by the ligase only when their ends are brought into close proximity by hybridization to the template DNA. Once the ligase has connected the two oligonucleotides, the product of the ligation is denatured from the target DNA so that both the target DNA and the ligated primers can serve again as templates for amplification.



**18-15 Principle of Q $\beta$  replicase amplification system.** In the Q $\beta$  system, the enzyme Q $\beta$  replicase amplifies the signal of the probe and not the target nucleic acid itself. The probe is an RNA molecule known as the midvariant (MDV-1) and a specific RNA fragment corresponding to a complementary piece of RNA or DNA in the organism to be detected. The probe RNA binds to the complementary target, and after removing the unbound MDV-1 molecules, the Q $\beta$  replicase is added to replicate the RNA probe.



**18-18 Q $\beta$  amplification for HIV-1 infected cells detected by fluorescent label.** Very low numbers of HIV-1 infected cells can be detected with this method. (From Pritchard CG, Stefano JE: *Detection of viral nucleic acids by Q $\beta$  amplification*. In Medical Virology 10: de la Maza LM and Peterson EM, editors. New York, 1991. Plenum Press.)



**18-17 Self-sustained sequence replication (3SR).** The avian myoblastosis virus (AMV) reverse transcriptase, the T7 bacteriophage RNA polymerase, and RNase H are utilized in the isothermic amplification of target nucleic acid. The promoter sequence for the T7 RNA polymerase and a region complementary to the target sequence are included in the reaction to allow the polymerase to amplify the target and to provide specificity. The RNase denatures the RNA complement of the RNA-DNA intermediate formed during each amplification phase, allowing sequential amplification to proceed without further high temperature denaturation.



**18-18 Nucleic acid sequence-based amplification (NASBA) method.** Using components similar to those described for 3SR, this system also amplifies nucleic acids at room temperature.

# Index

- A**
- A 129 + 66, 129**
- Ablate varietal of 129**
- Ablate varietal**  
histological stain blue preparation of 129  
on *Sarcosina dentata* eggs 129
- Acanthamoeba cysts of 162**
- Acanthamoeba Kerker's 53**
- Acanthamoeba 74**
- Acanthamoeba kerkeri on MacOcker eggs 64**
- Acanthamoeba submicroscopica var. kerkeri on MacOcker eggs 26**
- Acanthamoeba Owen's stain of 29**
- Acanthamoeba**  
varieties of 26  
1 agent *microscopica* stain stain of 15
- Acanthamoeba vesalii** histological identification of 15
- Adenovirus 129 (67)**
- Adenovirus**  
cellular detection of 129  
culture identification of with *Sarcosina* stain 149  
cytopathic effect of 129  
direct fluorescent assay of 129  
electron microscopy histology stain 149  
immunofluorescence microscopy of 149
- Adenovirus type 129 (serogroup B) 144**
- Adenovirus 12**  
Gross stain of 26
- Aerobic cultures method on fish cultures method**
- Aerobiosis**  
microaerophilic 1  
incubation of 29-32
- Aeromonas 37**
- Aeromonas histological stain blue preparation of 179**
- Aerobius 194-195**
- Aerobic acid hydrolysis reactions**  
of *Sarcosina arbutus* 94  
of *Sarcosina brucei* 94  
of *Sarcosina* 94
- Aeromonas 89-94**  
morphology 93  
culturing for 92  
identification by microscopical data of 92  
on system 92  
growth 93
- Aeromonas hydrophila** egg of 92
- Aerobic detection 194-195**
- Aeromonas susceptibility testing 192-193**  
summary 193  
with *microscopica* 192  
with *hydrophila*, 193  
1149 197 193
- Aeromonas hydrophila**  
histological stain blue preparation of 69
- Aeromonas hydrophila var. aeromonas 147**
- Aeromonas hydrophila - for identification features 1**
- Aeromonas 45**
- Aeromonas hydrophila var. aeromonas varieties of 37**  
Gross stain stain  
histology of gross stain 161  
histology of gross stain 161
- Aeromonas Gross stain of 31**  
histology Gross stain 174
- Aeromonas Gross**  
histological stain blue preparation of 129  
on Gross blood eggs 129
- Aeromonas hydrophila**  
Gross stain of 129  
histological stain blue or colorless white preparation of 131  
morphology of 131  
on Gross blood eggs 129
- Aeromonas Gross**  
histological stain blue preparation of 129  
morphology of 129  
on *Sarcosina dentata* eggs 129
- Aeromonas Gross**  
histological stain blue preparation of 129  
on *Sarcosina dentata* eggs 129
- Aeromonas hydrophila**  
histological stain blue preparation of 129  
morphology of 129  
on *Sarcosina dentata* eggs 129
- Aeromonas hydrophila water 94**
- Aeromonas 9**
- B**
- Bacterial blood stain of 159**
- Bacterial culture media using stain preparation testing of 157**
- Bacteria**  
varieties of 66-67  
egg, pink eggs 67  
growth hydrolysis test for 67  
Gross stain 27-68
- Bacteria culture** colonies of 66
- Bacterial susceptibility test 19 29**
- BATCC blood culture bottles 19 29**
- Bacterial culture media eggs**  
Gross stain eggs on 66  
194-195 *microscopica* on 66
- Bacterial culture**  
on *Sarcosina* 194-195 eggs (BAC) 19, 21  
on *Sarcosina* blood eggs 66-68  
Gross stain 66  
194-195 193

- Bacteroides urealyticus*, identification of, 90
- Balantidium coli*  
cyst of, in feces, 151  
trophozoites of, in feces, 151
- Bay agar, 53
- Beta-lactamase  
cefinaise for detecting production of, 112  
positive test for, 112
- Bile esculin slant, 44  
of *Enterococcus*, 45
- Bile solubility test for *Streptococcus pneumoniae*, 42
- Bio-Bag Environmental System, 11
- Biological safety cabinet, 6, 7
- BioMerieux Vitek colorimeter, 111
- Bipolaris*  
germ tube preparation of, 127  
lactophenol cotton blue preparation of, 127
- Blastocystis hominis* in feces, 149
- Blastomyces dermatitidis*  
colonies of, 113  
Gomori methanamine silver stain, 114  
lactophenol cotton-blue, preparation of, 113  
mycelial phase of, 113  
yeast form of, in lung tissue, 114
- Blood culture bottles, 12, 13, 23, 24
- Blood specimens  
collection of, 12, 13
- Bordetella bronchiseptica*  
on sheep blood agar, 76  
urea identification of, 77
- Bordetella pertussis*  
recovery of, 23  
specimen collection of, 10
- Bordetella versutus Pasteurella*, 76
- Borrelia* sp., Gram stain, 30
- Borrelia burgdorferi*, indirect immunofluorescent test for antibodies to, 197  
Gram stain, 30
- Botryomyces caespitosus* in chromoblastomycosis, 129
- Bottles  
roller, 178  
spinner, 178
- Branhamella* (see *Moraxella*)
- Broth microdilution method, 110
- Brucella*, 72  
on chocolate agar, 77  
Gram stain of, 29, 76  
on sheep blood agar, 77
- Burkholderia (Pseudomonas) cepacia* on MacConkey agar, 87
- C**
- Caliciviridae, 181
- Calicivirus, electron microscopy, negative stain of, in feces, 181
- CAMP test  
for group B streptococci, 40  
for *Clostridium perfringens*, 93
- Campylobacter*, 72  
Gram stain of, 30, 80  
nalidixic acid susceptibility, 80  
recovery of, 23  
on TSI slant, 80
- Campylobacter jejuni*  
nalidixic acid susceptibility testing for, 80  
on TSI slant with lead acetate strip, 80
- Candida*, Gram stain of, 31
- Candida albicans*  
on corn meal agar with trypan blue, 123  
germ tube test for, with calcofluor white, 123  
Gram stain of blood culture of, 123  
on sheep blood agar, 122
- Candida glabrata* (see *Torulopsis glabrata*)
- Candida tropicalis*, Gram stain of blood culture of, 123
- Capnocytophaga*, 72  
on chocolate agar, 79  
on sheep blood agar, 79
- Cardiobacterium*, Gram stain of, 29
- Catalase test, 33
- Cefinase, detection of beta-lactamase production by, 112
- Centrifuge tubes, safety of, 7, 8
- Cestodes, 146, 169–173
- Charcot-Leyden crystals versus parasites in stool, 175
- Chilomastix mesnili*  
cyst of, in feces, 151  
trophozoite of, in feces, 150
- Chlamydia*, specimen collection of, 10  
transport medium, 11
- Chlamydia pneumoniae*  
shell vial culture of, 103  
transmission electron microscopy of, 104
- Chlamydia psittaci*, shell vial culture of, 103
- Chlamydia trachomatis*  
collection kit, 10  
direct antigen detection of, 103  
direct antigen test of, 181  
direct fluorescent assay of, 102  
shell vial culture of, 103  
transmission electron microscopy of, 104
- Chloramphenicol, susceptibility testing of, 109, 110
- Chromoblastomycosis  
fungi associated with, 129  
Gomori's methanamine silver stain of, 129
- Ciliates, 150–151
- Citrate, utilization of, 60, 64
- CIT/RIIAM medium, reactions to, 64
- Cladosporium*, 128
- Cladosporium carrionii*  
in chromoblastomycosis, 129  
lactophenol cotton blue preparation of, 127  
on Sabouraud's dextrose agar, 127
- Clindamycin, susceptibility testing of, 109
- Clonorchis sinensis*  
adult, 167  
eggs of, 167

- Clostridium**  
 on egg yolk agar, 89  
 Gram stain of, 93  
 identification of, 90
- Clostridium difficile***  
 on Brucella blood agar, 94  
 cytotoxin cell culture assay for, with MRC-5 cells, 102  
 on egg yolk agar, 94
- Clostridium parvum***, Gram stain of, 94
- Clostridium perfringens***  
 on Brucella blood agar, 93  
 on egg yolk agar, 93, 94  
 reverse CAMP test for, 93
- Clostridium tertium*** on Brucella blood agar, 94
- Clostridium tetani***  
 on Brucella blood agar, 94  
 Gram stain of, 27
- Clothing, protective, 1, 2, 3, 4
- Coagulase test, *Staphylococcus* differentiation with, 34  
 slide test, 34  
 tube test, 34
- Coats, laboratory, 3
- Coccidia, 152–153
- Coccidioides immitis***  
 arthroconidia, 155  
 complement fixation test for detecting antibodies to, 199  
 lactophenol cotton-blue preparation of, 115  
 on Gomori's methenamine silver preparation, 115  
 passive immunoprecipitation test for detecting antibodies to, 198  
 on Sabouraud's dextrose agar, 114  
 wet mount and calcofluor white preparations of, 114
- Complement fixation test, 196  
 for detecting antibodies to *Coccidioides immitis*, 199
- Coronaviridae, 181
- Container  
 puncture-proof, 5  
 safety bucket, 7  
 mailing, 8
- Coronavirus, electron microscopy, negative stain of, in feces, 181
- Corynebacterium***, 46  
 Gram stain of, 27, 48
- Corynebacterium diphtheriae***, colonies of, 48, 49  
 Loeffler's methylene blue stain, 49  
 Tinsdale agar, 49
- Corynebacterium jeikeium***, colonies of, 48
- Cryptococcus neoformans***  
 on chocolate agar, 115  
 differential identification of, 115  
 India ink and Gram stain of, 31, 116  
 mixed culture of, with *Histoplasma capsulatum*, 118  
 urease test for, 116
- Cryptosporidium parvum***, oocysts of, in feces, 153
- CTA (see cystine-trypticase agar)
- Culture(s)  
 aerobic  
 collection and transport of, 9  
 incubation of, 21  
 anaerobic  
 collection and transport of, 11  
 incubation of, 21  
 media, 16, 23
- Cysticercus**, *Taenia solium* in, 170
- Cystine-trypticase agar (CTA) for identifying *Neisseria* species, 73
- Cytomegalovirus**, 183  
 antibody detection of, 185  
 cytopathic effect of, 184  
 hematoxylin and eosin, 185  
 horseradish peroxidase, 185  
 in kidney and lung, 185  
 shell vial of, with immunofluorescent stain, 184
- Cytotoxin cell culture assay with MRC-5 cells for *Clostridium difficile*, 102
- D**
- Decarboxylase-dihydrolase test, 62
- Diarrhea, antibiotic-associated, 102
- Dientamoeba fragilis***, trophozoite of, in feces, 150
- Diphyllobothrium latum***, egg of, in feces, 169
- Diplococci, gram-positive, 26
- Dipylidium caninum***  
 proglottids of, 172  
 scolex of, 172
- Disinfection, 4
- Disk diffusion method, 107, 108  
 with *Escherichia coli* ATCC 25922, 108  
 with *Pseudomonas aeruginosa* ATCC 27853, 109  
 with resistant strain of *Pseudomonas aeruginosa*, 109  
 with *Staphylococcus aureus* ATCC 25923, 108  
 with *Streptococcus pneumoniae*, 109
- Disposal, waste, 4, 5
- DNA hybridization, 200  
 of herpes simplex virus, 200  
 of human papilloma virus, 201
- DNA probes, fungi, 145
- DNA restriction patterns in herpes simplex virus, 201
- DNA sequencing, 202
- DNase agar, 86
- E**
- E test method for *Streptococcus pneumoniae*, 109, 110
- Echinococcus granulosus***  
 adult worm, 172  
 eggs of, 173  
 hydatid sand, 173  
 hydatid cyst in liver, 173
- Echovirus 11, cytopathic effect of, in MRC-5 cells, 191
- EIA (see enzyme immunoassay)
- Eikenella***, 72  
 Gram stain of, 29  
*Eikenella corrodens* on sheep blood agar, 78
- Electron microscopy  
 of *Adenovirus*, 180  
 of *Calicivirus*, 181  
 of *Chlamydia pneumoniae*, 104

## Electron microscopy (continued)

- of *Chlamydia trachomatis*, 104
  - of *Cornavirus*, 181
  - of Herpes virus, 183
  - of HIV-1, 192
  - of Rotavirus, 192
  - of *Ureaplasma urealyticum*, 106
  - of *Vaccinia virus*, 192
- ELISA (see enzyme immunoassays)
- Endocarditis, agents of, 29
- Endolimax nana*  
cyst of, in feces, 149  
trophozoite of, in feces, 148
- Entamoeba coli*  
cyst of, in feces, 148  
trophozoite of, in feces, 148
- Entamoeba hartmanni*  
cyst of, in feces, 148  
trophozoites of, in feces, 147
- Entamoeba histolytica*  
cysts of, in feces, 147  
indirect hemagglutination test for detecting antibodies to, 198  
trophozoites of  
in liver, 147  
in feces, 146
- Enterobacter aerogenes*, characteristic reactions of, in API 20 E strips, 68, 69
- Enterobacter cloacae*, characteristic reactions of, in API 20 E strips, 68
- Enterobacteriaceae, 56–71  
characteristics of, 56  
citrate utilization test for, 60  
decarboxylase-dihydrolase test for, 62  
Gram stain, 30  
on Hektoen enteric agar, 58  
indole spot test for, 59  
on MacConkey agar, 58  
methyl red-Voges-Proskauer test for, 59  
ONPG test for, 59  
oxidase spot test for, 59  
PAD test for, 61  
phenylalanine deaminase test for, 60  
rapid spot test identification of, 58  
on sheep blood agar, 56  
triple sugar iron test for, 62  
urea test for, 61  
urease test for, 61
- Enterobius vermicularis*  
adult female and male, 161  
eggs of  
cellophane tape preparation of, 162  
in feces, 162
- Enterococcus*, 37, 44  
bile esculin agar, 45  
colonies of, 45  
NaCl broth, 6.5%, 45
- Enzyme immunoassays (EIA), 194  
direct, 195  
indirect, 195
- Epidermophyton floccosum*  
lactophenol cotton blue preparation of, 132  
on Sabouraud's dextrose agar, 132

## Epithelial cells, ciliated columnar, versus parasites, in feces, 175

- Epstein-Barr virus, 183  
commercial particle agglutination test for serological diagnosis of, 186  
heterophile antibody, 186  
immunofluorescent tests for antibodies to nuclear antigens of, 187  
immunofluorescent tests for antibodies to early and structural antigens of, 187
- Erysipelothrix*, 46
- Erysipelothrix rhusiopathiae* on TSI agar, 50  
colonies of, 50
- Erythromycin, susceptibility testing of, 109
- Escherichia coli*  
characteristic reactions of, 67  
in API 20 E strip, 68  
Gram stain of, 30  
on MacConkey agar, 57  
in Micro-ID, 67  
on sheep blood agar, 57
- Esculin hydrolysis test, 48
- Exophiala (Wangiella) dermatitidis*, lactophenol cotton blue preparation of, 128
- Exophiala jeikei*  
lactophenol cotton blue preparation of, 128  
on Sabouraud's dextrose agar, 128
- Exserohilum*, lactophenol cotton blue preparation of, 128
- Eye wear, 2

## F

- Face mask  
conditions requiring, 2  
improperly fitted, 3
- Face shield, 3
- Fasciola hepatica*, egg of, in feces, 167
- Filaria, 165–166
- Flagellates, 150–151
- Flasks, tissue culture, 177
- Flatworms; see Trematodes
- Flavobacterium*, 72
- Flavobacterium meningosepticum* on sheep blood agar, 86
- Fonsecaea compacta* in chromoblastomycosis, 129
- Fonsecaea pedrosoi*  
in chromoblastomycosis, 129  
lactophenol cotton blue preparation of, 129
- Fungi; see also Mycology; specific fungi  
classification and characteristics of, 113  
identification of, 113
- Fusarium*  
lactophenol cotton blue preparation of, 140  
on potato dextrose agar, 140
- Fusobacterium*  
egg yolk agar for identifying, 89  
Gram stain of, 30  
identification of, 90  
spot indole disk test of, 89
- Fusobacterium nucleatum*  
on Brucella blood agar, 92  
Gram stain of, 30, 92

## G

- Gamma hemolysis, 38
- Gardnerella*, 46
- Gardnerella vaginalis*  
colonies of, 50  
Gram stain of, 50  
hippurate hydrolysis test for, 50
- GasPak  
anaerobic jar system, 22  
anaerobic pouch, 23
- Gelatin hydrolysis test, 47
- Gemella*, 37
- Geotrichum* on corn meal agar with trypan blue, 124
- Giardia lamblia*  
cyst of, in feces, 150  
trophozoite of, in feces, 150
- Globicatella, 37
- Gloves, removal and disposal of, 1, 2
- Goggles, conditions requiring, 2
- Gram stain, 25–31  
preparation of, 24
- Graphium conidia*, lactophenol cotton blue stain of, 130

## H

- HACEK, Gram stain of, 29
- Haemophilus*, 72  
Gram stain of, 29, 30
- Haemophilus aphrophilus*, Minitex sugar identification of, 76
- Haemophilus influenzae*  
on chocolate agar, 74  
Latex particle agglutination test for detecting, 199  
porphyrin production test for, 75  
on X and V strips, 75  
satelliting, 75
- Haemophilus parainfluenzae*  
on chocolate agar, 75  
Minitex sugar identification of, 76
- Hektoen enteric agar, Enterobacteriaceae on, 58
- Helminths, classification of, 146
- Hemolysis  
alpha, 37  
gamma, 38
- Hepatitis C, EIA detection of, 194
- Herpes simplex virus  
antibody detection of, 184  
cervical swab of, 182  
cytopathic effect of, 182  
direct antigen test of, 181  
direct immunofluorescent assay of, 182  
DNA probe, 200  
DNA restriction patterns in, 201  
horseradish peroxidase stain, 182  
shell vial culture of, 183  
skin lesions with, 181  
transmission electron microscopy of, 183
- Herpes virus  
collection kit, 10  
radiolabeled probe detection of, 200  
specimens of, 10

- Herpesviridae, 181–186
- Heterophile antibody, 186
- Hippurate hydrolysis test, 40  
for *Gardnerella vaginalis*, 50  
for group B streptococci, 40  
for *Listeria monocytogenes*, 48
- Histoplasma capsulatum*  
on buffy coat preparation (Wright's stain), 117  
colonies of, 116  
differential identification of, 115  
on direct bone marrow preparation stained with calcofluor white, 117  
lactophenol cotton blue stained colonies of, 117  
lung section showing, with Gomori's methenamine silver, 117  
mixed culture of, with *Cryptococcus neoformans*, 118  
yeast form of, 116
- Homogenization of specimens, 17
- Hookworm, egg of, 163
- Human immunodeficiency virus type 1  
antibody detection of  
with enzyme immunoassay, 193  
with immunofluorescent assay, 193  
with Western blot testing, 193  
positive Q $\beta$  amplification for, 204  
transmission electron microscopy, 192
- Human papilloma virus, DNA hybridization of, 201
- Hymenolepis diminuta*, egg of, in feces, 172
- Hymenolepis nana*, egg of, in feces, 171

## I

- Immunofluorescence, direct and indirect, 196
- Immunoserology, 194–199
- Indole spot test, 59  
for anaerobes, 89
- Influenza A virus  
African green monkey kidney cell isolation of, 188  
direct fluorescent assay of, 187  
hemadsorption, 188
- Influenza virus  
antibody detection of, 188  
shell vial culture of, 188
- Inoculum  
preparation of, 107  
turbidity meter for adjusting, 107
- Interstate Shipment of Etiologic Agents, 8
- Iodamoeba bütschlii*  
cyst of, in feces, 149  
trophozoite of, in feces, 149
- Iron uptake test, 101
- Isolatorysis-centrifugation blood culture system, 12
- Isospora belli*, cyst of, in feces, 152
- J
- JC virus in progressive multifocal leukoencephalopathy, brain tissue section of, 189
- Jemec plate, 11



- K**
- Kingella*, 72  
Gram stain of, 29
- Kingella kingae* on sheep blood agar, 78, 79
- Kinyoun acid-fast stain, 95  
modified, for *Nocardia*, 53  
modified, for *Rhodococcus equi*, 52
- Klebsiella oxytoca*, characteristic reactions of, to API 20 E strips, 69
- Klebsiella pneumoniae*  
on MacConkey agar, 57  
urease test for, 61
- Klebsiella versus Rhodococcus equi*, 51
- L**
- Laboratory coats, 3
- Laboratory safety, 1–8  
and biological safety cabinets, 6, 7  
and disinfection, waste disposal, and sterilization, 4, 5, 6  
mailing containers and, 8  
miscellaneous practices and, 7, 8  
and protective clothing and equipment, 1, 2, 3, 4
- Lactobacilli, Gram stain of, 27
- Lactobacillus*, 46  
colonies of, 49  
Gram stain of, in vaginal secretions, 49
- Lactococcus*, 37, 44
- Latex particle agglutination test for detecting *H. influenzae*, *S. pneumoniae*, *N. meningitidis*, and group B streptococci, 199
- LCR (ligase chain reaction)
- Lead acetate strip, *Campylobacter jejuni* on, 80
- Legionella pneumophila* on buffered charcoal-yeast extract agar, 87
- Leishmania*  
amastigotes of, 159  
promastigotes of, 159
- Leuconostoc*, 37, 44  
colonies of, 45
- Leukoencephalopathy, progressive multifocal, JC virus in, 189
- Ligase chain reaction (LCR), 203
- Liquid nitrogen tanks, 179
- Listeria*, 46
- Listeria monocytogenes*  
colonies of, 47  
esculin hydrolysis test for, 48  
Gram stain of, 27  
hippurate hydrolysis, 48
- Liver, *Echinococcus granulosus* in, 173
- Loa loa*, microfilariae of, 165
- Loeffler methylene blue stain of *C. diphtheriae*, 49
- M**
- MRVP broth test (see Methyl red-Voges-Proskauer)
- Mailing containers, safety requirements for, 8
- Malaria, 155–158  
development of, 155
- Malassezia furfur*  
Gram stain and calcofluor white preparations of, 133  
on Sabouraud's dextrose agar, 132
- Mannitol salt agar, *Staphylococcus aureus* differentiation with, 35
- Measles virus  
antibody detection of, with indirect immunofluorescent assay, 190  
in lung tissue, 190  
phytohemagglutinin stimulated cord blood mononuclear cell culture of, 189  
Vero cell culture of, 190
- Methyl red-Voges-Proskauer (MRVP) broth test, 59
- MIC, defined, 107
- Micrococcaceae, 32–36  
morphology of, 32
- Micrococcus*  
Gram stain of, 33  
infections due to, 32  
versus *Staphylococcus*, 32
- Microhemagglutination assay for detecting antibodies to *Treponema pallidum*, 198
- MICRO-ID strip, *E. coli* reactions to, 67
- Microimmunofluorescence test  
for *Mycoplasma* and *L. repleans*, 106  
for rickettsial antibodies, 197
- Microscopes  
brightfield, with epiluminescence, 176  
inverted, 176
- Microsporidium*, 152–153  
in feces, 153
- Microsporium audouinii*  
lactophenol cotton blue preparation of, 133  
on Sabouraud's dextrose agar, 133
- Microsporium canis* var. *canis*  
lactophenol cotton blue preparation of, 134  
on Sabouraud's dextrose agar, 134
- Microsporium gypseum*  
lactophenol cotton blue preparation of, 135  
on Sabouraud's dextrose agar, 134
- Minitek sugars for identifying *Haemophilus*, 76
- Molds, specimen preparation and identification systems for, 144–145
- Molecular techniques, 200–204
- Monospot test, 188
- Moraxella*, 72
- Moraxella (Branhamella) catarrhalis* on sheep blood agar, 73
- Morganella morganii*  
characteristic reactions of, in API 20 E strips, 70  
PAD and urea tests for, 61
- MRC-5 cells, 102
- Mucor*  
Gomori's methenamine silver stain of, 126  
lactophenol cotton blue preparation of, 126  
structure of, 124
- Multichannel pipettor, 194
- Mycobacterium*, 95–101

- Mycobacterium avium*  
on Löwenstein-Jensen agar, 97  
on Middlebrook 7H11 agar, 97  
rough colony of, 97  
smooth colony of, 97
- Mycobacterium chelonae*  
on chocolate agar, 99  
colony of, 99  
on Middlebrook 7H11 agar, 99
- Mycobacterium fortuitum*  
colony of, 100  
iron uptake test for differentiating, 101  
on Löwenstein-Jensen agar, 100  
on Middlebrook 7H11 agar, 99
- Mycobacterium goodii*  
on Löwenstein-Jensen agar, 100  
on Middlebrook 7H11 agar, 100, 101
- Mycobacterium kansasii*  
colony of, 98  
on Löwenstein-Jensen agar, 98  
on Middlebrook 7H11 agar, 98
- Mycobacterium simiae*, niacin accumulation testing of, 101
- Mycobacterium tuberculosis*  
colonies of, 96  
on Löwenstein-Jensen agar slant, 96  
on Middlebrook 7H11, 96  
niacin accumulation testing of, 101
- Mycology, 113–145
- Mycoplasma*  
antibody detection of, by microimmunofluorescence test, 106  
characteristics of, 105
- Mycoplasma hominis* colonies, 105
- Mycoses  
deep-seated, 113–118  
miscellaneous, 140–143  
opportunistic, 113, 119–126  
specimen preparation and identification systems for, 144–145  
subcutaneous, 127–131  
superficial, 132–139  
systemic, 113
- N**
- Nalidixic acid, *Campylobacter jejuni* susceptibility to, 80
- Nasopharyngeal  
collection system, 10  
-urogenital swab, 9
- National Committee for Clinical Laboratory Standards, 107
- Necator americanus*  
egg of, 163  
first-stage rhabditiform larvae of, 163
- Needles, disposal of, 5
- Neisseria*, 72  
API Quad-FERM+ identification of, 74  
carbohydrate utilization identification of, 73  
Gram stain of, 29, 72  
*Neisseria/Haemophilus* identification card for identifying, 74
- Neisseria gonorrhoeae*  
Gram stain of, 29  
on modified Thayer-Martin agar, 73  
recovery of, 23  
specimen collection of, 9

- Neisseria meningitidis*  
on chocolate agar, 72  
Gram stain of, 29  
latex particle agglutination test for detecting, 199  
on sheep blood agar, 73
- Neisseria/Haemophilus* identification (NHI) card for identifying *Neisseria*, 74
- Nematodes, 146, 161–166  
intestinal, 161–164  
tissue, 164–166
- NHI card (see *Neisseria/Haemophilus* identification card)
- Niacin accumulation test, 101
- Nitrate test, 77  
disk test, 89
- Nocardia*  
biochemical reactions for identifying, 54  
colonies of, 53  
Gram stain of, 53  
on Bay agar, 53
- Nocardia asteroides*  
amino acid hydrolysis reactions of, 54  
colonies of, 53
- Nocardia brasiliensis*, amino acid hydrolysis reactions of, 54
- Novobiocin susceptibility test, 36
- Nucleic acid amplification thermocycling instrument, 202
- Nucleic acid sequence-based amplification method, 204
- Nucleic acidification tests, positive displacement pipette for, 203

**O**

- ONPG test (see O-nitrophenyl- $\beta$ -D-galactosidase test)
- Occupational Safety and Health Act, safety provisions of, 1
- Orskovia*, 46  
colonies of, 52  
Gram stain of, 28
- Oligonucleotide ligation amplification, 203
- Ouchocerca volvulus*, microfilaria, 166
- O-nitrophenyl- $\beta$ -D-galactosidase (ONPG) test, 59
- Optochin susceptibility test  
for *Streptococcus pneumoniae*, 42  
of viridans streptococci, 43
- Orthomyxoviridae, 187–189
- OSHA; see Occupational Safety and Health Act
- Oxacillin  
screen plate, 35  
*Staphylococcus aureus* resistance to, detection of, 35  
susceptibility testing of, 109
- Oxidase spot test, 59
- Oxidative-fermentative medium, 33, 82  
glucose oxidation in, 33  
glucose fermentation in, 33

**P**

- PAF test (see phenylalanine deaminase test)
- Paezilomyces variotti*  
lactophenol cotton blue preparation of, 141  
on Sabouraud's dextrose agar, 141

- Papovaviridae, 189
- Paracoccidioides brasiliensis*  
in adrenal tissue, stained with Gomori's methenamine silver, 118  
differential identification of, 115  
stained with lactophenol cotton blue, 118
- Paragonimus westermani*, egg of, in feces, 168
- Parainfluenza 3 and 4 viruses, immunofluorescent staining of, 189
- Paramyxoviridae, 189–191
- Para-Pak Parasitology Collection Kit, 12
- Parasites, specimen collection of, 12
- Parasitology, 146–175
- Pasteurella*, 72
- Pasteurella multocida* on sheep blood agar, 78
- PCR (polymerase chain reaction)
- Pediococcus*, 37  
colonies of, 44  
vancomycin susceptibility testing of, 44
- Penicillin  
*Streptococcus pneumoniae* resistance to, 43
- Penicillium*  
*disiagram* of, 141  
on Sabouraud's dextrose agar, 142  
stained with calcofluor white, 142
- Penicillium marneffei*  
lactophenol cotton blue and calcofluor white preparations of, 143  
on Sabouraud's dextrose agar, 142
- Peptostreptococcus*, spot indole disk test of, 89
- Peptostreptococcus anaerobius* on Brucella blood agar, 92  
identification of, 90
- Peptostreptococcus magnus*  
on Brucella blood agar, 92  
SPS resistance of, 92
- Peptostreptococcus anaerobius*, sodium polyanethol sulfonate disk test of, 89
- Personal Protective Devices, 1, 2, 3, 4
- Phenylalanine deaminase (PAD) test, 60, 61
- Phialophora*, 128
- Phialophora verrucosa*  
in chromoblastomycosis, 129  
lactophenol cotton blue preparation of, 129
- Picornaviridae, 191–192
- Pinworm, 161
- Pipet, positive displacement, 203
- Pipettor, multichannel, 194
- Plant cells versus parasites, in feces, 174
- Plantar spiral versus parasites in feces, 175
- Plasmodium falciparum*  
gametocyte appearance and, 157  
ring forms of, 155  
thick blood smear of, 155, 157  
trophozoites of, 155
- Plasmodium malariae*  
schizonts and, 157  
trophozoites of, 156
- Plasmodium ovale*  
erythrocyte appearance and, 156  
Schüffner's dots and, 156  
thick blood smear of, 157
- Plasmodium vivax*  
macrogametocytes, 158  
microgametocytes, 158  
versus *P. falciparum*, 156  
Schüffner's dots and, 156  
thick blood smear of, 157
- Plates, multiple well, 177
- Pneumocystis carinii*  
from bronchoalveolar lavage, 154  
from lung biopsy, 154  
in sputum, 153
- Pollen versus parasites in feces, 174
- Polymerase chain reaction (PCR)  
automated instrument for, 203  
basic principle of, 203  
products of, on agarose gel, 203
- Polymorphonuclear leukocytes, 26, 27  
Gram stain of, 28, 53  
with gram-positive bacilli, 75  
poorly differentiated, 30
- Porphyrin production test for *Haemophilus influenzae*, 75
- Porphyromonas*, spot indole disk test of, 89
- Positive displacement pipette for nucleic amplification tests, 203
- Poxviridae, 192
- Prevotella*, spot indole disk test of, 89  
identification of, 90
- Prevotella melaninogenica*  
on Bacteroides bile esculin agar, 91  
on Brucella blood agar, 91  
Gram stain of, 91
- Progressive multifocal leukoencephalopathy, JC virus in, 189
- Propionibacterium*, spot indole disk test of, 89
- Propionibacterium acnes* on Brucella blood agar, 93
- Proteus*  
Gram stain of, 30  
on sheep blood agar, 56
- Proteus mirabilis*  
characteristic reactions of  
in API 20 E strips, 69, 70  
PAD and urea tests for, 61  
urease test for, 61
- Proteus vulgaris*  
characteristic reactions of, 66  
to API 20 E strips, 70  
with API 20 E strips, 69  
PAD and urea tests for, 61
- Protozoans, 146–160  
coccidia and microsporidia, 152–153  
flagellates and ciliates, 150–151  
intestinal, 146–149  
in tissue, 153–154
- Providencia rettgeri*, PAD and urea tests for, 61
- Pseudallescheria boydii*  
on Sabouraud's dextrose agar, 130  
sexual state of, lactophenol cotton blue stain of, 130

Pseudomonadaceae, 82–86

*Pseudomonas*, 72

Gram stain of, 30

*Pseudomonas aeruginosa*

GNF, urea, and Pseudosel agar for identifying, 84

GNF tube identification of, 83

on MacConkey agar, 83

pigment production by, 83

resistant strain of, disk diffusion method for, 109

on sheep blood agar, 82

on TSI slant, 83

Uni-N/F-Tek plate for identifying, 84

*Pseudomonas aeruginosa* ATCC 27853, disk diffusion method for, 109

*Pseudomonas cepacia*, (see *Burkholderia cepacia*)

*Pseudomonas paucimobils*, (see *Sphingomonas paucimobils*)

*Pseudomonas putrefaciens*, (see *Shewanella putrefaciens*)

PVR test, 39

## Q

Q $\beta$  replicase amplification

for HIV-1 infected cells detected by fluorescent label, 204

principle of, 204

## R

Radiolabeled probe for detecting herpesvirus, 200

Rapid plasma reagin (RPR) test, 199

Rapid spot test for Enterobacteriaceae, 58

r/b<sub>2</sub> medium, reactions to, 64

r/b tube

for Enterobacteriaceae, 62

reactions to, 63

Reoviridae, 192

Respiratory syncytial virus

culture identification of, in MRC-5 cells in shell vial, 191

cytopathic effect of, in MRC-5 cells, 191

direct fluorescent assay of, 190

Restriction fragment length polymorphism analysis, 201

Retroviridae, 192–193

Rhamnose, fermentation of, 64

*Rhinocladella*, 128

*Rhinocladella aquaspersa* in chromoblastomycosis, 129

Rhinovirus, cytopathic effect of, in MRC-5 cells, 192

Rhizopus

lactophenol cotton blue preparation of, 126

structure of, 124

*Rhodococcus*, 46

colonies of, 51

microscopic appearance of, 52

*Rhodococcus equi*, classification of, 51

acid fast stain, 51

*Rickettsia*, microimmunofluorescence test for detecting antibodies to, 197

Rifampin, susceptibility testing of, 109

Roller bottles, 178

Roller drum, 177

Rotavirus, electron microscopy of, in stool, 192

Roundworms: see Nematodes

RPR (see rapid plasma reagin test)

## S

Safety: see Laboratory safety

Safety cabinet, biological, 6, 7

Saliva versus sputum, 18

*Salmonella*

characteristic reactions of, 65

on Hektoen enteric agar, 58, 71

triple sugar iron agar test for, 62

*Salmonella enteritidis*, colonial morphology versus biochemical identification of, 71

*Salmonella typhi*, characteristic reactions of, 65

Salt broth, 6.5%, 45

Samples: see Specimens

*Scedosporium apiospermum*, lactophenol cotton blue stain of, 130

*Scedosporium prolificans*, lactophenol cotton blue stain of, 131

Schöffner's dots, 156

*Schistosoma haematobium*, egg of, 168

*Schistosoma japonicum*

adult male of, 168

cercaria, 168

egg of, 168

*Schistosoma mansoni*, egg of, in stool, 169

Schlichter test, 112

*Scopulariopsis*

lactophenol cotton blue preparation of, 143

on Sabouraud's dextrose agar, 143

Self-sustained sequence replication, principle (SSSR), 204

Septi-Chek biphasic blood culture bottle, 12

*Serratia marcescens*, characteristic reactions of, with API 20 E strips, 69

*Serratia* on MacConkey agar, 57

Serum inhibitory titer, 112

Shell vial, 177

*Shewanella*, 72

*Shewanella (Pseudomonas) putrefaciens*

on sheep blood agar, 85

TSI reaction of, 85

*Shigella*, characteristic reactions of, 67

*Shigella sonnei*, characteristic reactions of, to API 20 E strips, 70

Skin, protection of, 4

Sleeping sickness, 159–160

Slide culture method for fungi, 144

Slides

preparation of, 19, 20

wet, 20

Smears

fixing and staining of, 24

preparation of, 19

Sodium chloride test

6.5%, 45

tolerance test, 101

- Sodium polyanethol sulfonate disk test  
for *Peptostreptococcus anaerobius*, 89  
for *Peptostreptococcus magnus*, 92
- Southern blot test, 202
- Specimens  
anaerobic  
media for, 16  
prereduced, 22
- Specimens (continued)  
blood, 12, 13  
of *Bordetella pertussis*, 10  
chlamydial, 10  
collection containers for, 9–13  
containers for, 15  
culture of, media for, 16  
discrepancy between microscopic and gross appearances of, 19  
disposal of, 5  
entering into computer system, 14, 15  
of expectorated sputa, 15  
homogenizing of, 17  
inspection of, before inoculation of media, 18  
labeling of, 14  
mailing containers for, 8  
of *Neisseria gonorrhoeae*, 9, 10  
processing of, 14–24  
atmospheric conditions for, 21, 22  
equipment for, 17  
procedure manual for, 16  
rapid direct tests for, 24  
quality of, 18, 19  
in syringes, 21  
transport of, 9  
urethral, collection of, 9  
viral, 10, 11
- Sphingomonas*, 72  
on sheep blood agar, 84
- Sphingomonas (Pseudomonas) paucimobilis* on sheep blood agar, 84
- Spinner bottles, 178
- Sporothrix schenckii*  
lactophenol cotton blue preparation of, 131  
on potato dextrose agar, 131
- Spot indole disk test, uses of, 89
- Sputum  
versus saliva, 18  
specimens of, 15, 19
- Squamous epithelial cells, vaginal, 27
- SSSR (see self-sustained sequence replication)
- Staphylococcus saprophyticus*, colonial morphology of, 36
- Staphylococcus*  
coagulase-negative, 32  
infections due to, 32  
versus *Micrococcus*, 33  
versus *Streptococcus*, 75  
with *Streptococcus*, 27  
susceptibility testing for, 107  
in upper respiratory tract, 38
- Staphylococcus aureus*, 32  
antimicrobial susceptibility testing of, 35  
colonies of, 32  
differentiation of  
by coagulase test, 34  
with mannitol salt agar, 35  
Gram stain of, 26  
methicillin (and oxacillin)-resistant, detection of, 35
- Staphylococcus saprophyticus*, novobiocin susceptibility differential testing of, 36
- Stenotrophomonas*, 72
- Stenotrophomonas (Xanthomonas) maltophilia*  
on sheep blood agar, 85  
Uni-Ni-Tek for identifying, 85
- Sterilization, 6
- Stomatococcus*  
Gram stain of, 26  
infections due to, 32
- Streptococcaceae, 37–45
- Streptococcus*, 37  
alpha-hemolytic, 37  
beta-hemolytic, 37  
gamma hemolysis, 37  
Gram stain, 27  
group A, 38, 39  
PYR test of, 39  
group A beta-hemolytic, bacitracin susceptibility test of, 39  
group B, 40  
CAMP test for, 40  
hippurate hydrolysis test of, 40  
latex particle agglutination test for detecting, 199  
group D, bile esculin slant and NaCl broth testing of, 45  
nonhemolytic, 38  
versus *Staphylococcus*, 75  
with *Staphylococcus*, 27  
in vaginal secretions, 40  
viridans  
conventional biochemical reaction testing of, 43  
Gram stain of, 43  
optochin susceptibility test of, 43
- Streptococcus bovis*, bile esculin slant testing of, 44
- Streptococcus pneumoniae*  
antimicrobial susceptibility testing of, 43  
bile solubility test for, 42  
colonies of, 41, 42  
E test method for, 109, 110  
Gram stain of, 41  
latex particle agglutination test for detecting, 199  
optochin susceptibility test for, 42  
secondary test battery for, 109
- Streptococcus pyogenes*, colonies of, 51
- Streptomyces*  
amino acid hydrolysis reactions of, 54  
Gram stain of, 55
- Strongyloides stercoralis*  
filariform larvae of, 164  
first-stage rhabditiform larvae of, 163, 164
- Swab, nasopharyngeal-urogenital, 9
- Syringes, specimens in, 21

## T

- Taenia*, eggs of, 171
- Taenia saginata*, proglottids of, 170  
eggs of, 171

- Taenia solium*  
 cysticercus, 170  
 eggs of, 171  
 proglottids of, 171  
 scolex of, 170
- Tapeworms; see Cestodes
- Tease mount technique, 144
- Tinsdale agar, 49
- Tissue culture flasks, 177
- Tissue culture tubes for viral detection, 176
- Tissue processing, 17
- Torulopsis (Candida) glabrata*, differential identification of, 115
- Toxocara canis*, embryonated egg and egg with larvae of, 165
- Toxoplasma gondii*  
 cyst of, in brain tissue, 154  
 indirect fluorescent method for detecting antibodies to, 197  
 trophozoites of, 154
- Transgrow bottle, 11
- Trematodes, 146, 167–169
- Treponema pallidum*, microhemagglutination assay for detecting antibodies to, 198  
 RPR test, 199
- Trichinella spiralis*, larvae of, in muscle, 164
- Trichomonas hominis*, trophozoite of, in feces, 151
- Trichophyton mentagrophytes*  
 Christiansen urea agar slant for demonstrating urease production by, 136  
 lactophenol cotton blue preparation of, 135  
 on Sabouraud's dextrose agar, 135
- Trichophyton rubrum*  
 Christiansen urea agar slant for demonstrating urease production by, 136  
 lactophenol cotton blue preparation of, 136  
 on Sabouraud's dextrose agar, 135
- Trichophyton schoenleinii*  
 lactophenol cotton blue preparation of, 136  
 on Sabouraud's dextrose agar, 136
- Trichophyton tonsurans*  
 lactophenol cotton blue preparation of, 137  
 nutritional test for, 138  
 on Sabouraud's dextrose agar, 137
- Trichophyton verrucosum*  
 lactophenol cotton blue preparations of, 138  
 on Sabouraud's dextrose agar, 138
- Trichophyton violaceum*  
 lactophenol cotton blue and calcofluor white preparations of, 139  
 on Sabouraud's dextrose agar, 139
- Trichuris trichiura*  
 adult male, 162  
 egg of, in feces, 162
- Triple sugar iron agar  
 for Enterobacteriaceae, 62  
 reactions to, 63
- Trypanosoma brucei gambiense*, blood of, 160
- Trypanosoma brucei rhodesiense*, blood of, 160
- Trypanosoma cruzi*, trypanomastigote of, 159
- Trypanosomiasis, 159–160
- TSI slant  
*Campylobacter jejuni* on, 80  
 Enterobacteriaceae on, 62, 63  
*Erysipelothrix rhusiopathiae* on, 50  
*Pseudomonas aeruginosa*, 83  
*Shewanella putrefaciens*, 85
- Turbidity meter, 107
- ## U
- Uni-N/F-Tek for identifying *Sienotrophomonas (Xanthomonas) maltophilia*, 85
- Ureaplasma*, antibody detection of, by microimmunofluorescence test, 106
- Ureaplasma urealyticum*  
 colonies of, 105  
 transmission electron microscopy of, 106
- Urease test, 61  
 for *Bordetella*, 77  
 for *Brucella*, 77  
 for *Cryptococcus neoformans*, 116  
 for Enterobacteriaceae, 61  
 for *Trichophyton mentagrophytes*, 136  
 for *Ureaplasma urealyticum*, 105
- ## V
- Vaccinia virus, transmission electromicroscopy, 192
- Vancomycin  
 anaerobe resistance to, 90  
*Penicillium* resistance to, 44  
*Peptostreptococcus magnus* resistance to, 92
- Vancomycin susceptibility test of *Penicillium*, 44
- Varicella zoster virus, 183  
 antibody detection of, 186  
 cytopathic effect, 186  
 shell vial of, fluorescent stain of, 186  
 in skin lesions  
 direct fluorescent assay of, 186  
 hematoxylin and eosin stain of, 185
- Veillonella, identification of, 90
- Vials, shell, 177
- Vibrio*, 72
- Vibrio cholerae*  
 on sheep blood agar, 81  
 on TCBS, 81  
 on thiosulfate-citrate-bile salts-sucrose (TCBS), 81
- Vibrio parahaemolyticus*  
 on sheep blood agar, 81  
 on TCBS, 81
- Virology, 176–193
- Viruses; see also specific viruses  
 collection and transport of, 10, 11  
 detection and identification of, 176–193
- VITEK system, 111
- VITEK yeast biochemical card, 145
- ## W
- Waste disposal, 5
- Wuchereria bancrofti*, microfilariae of, 166

## X

X and V strips, *Haemophilus influenzae* on, 75

*Xanthomonas* (see *Stenotrophomonas*)

Xylose-lysine-desoxycolate, Enterobacteriaceae on, 58

## Y

Yeast

form of, 25

versus parasites, in feces, 174

Yeast cells, Gram stain of, 31

*Yersinia enterocolitica*, characteristic reactions of, 66  
in API 20 E strips, 70

## Z

*Zygomycetes*, diagram of, 124